



00-276-009-007
Centrifuge with
Fixed Angle Rotor



Drucker Diagnostics 00-276-009-007 Centrifuge with Fixed Angle Rotor User Guide

[Home](#) » [Drucker Diagnostics](#) » Drucker Diagnostics 00-276-009-007 Centrifuge with Fixed Angle Rotor User Guide 

Contents

- [1 Drucker Diagnostics 00-276-009-007 Centrifuge with Fixed Angle Rotor](#)
- [2 Specifications](#)
- [3 TrueBond 00-376-009-005 Centrifuge](#)
- [4 TrueBond Duet](#)
- [5 FAQs](#)
- [6 Documents / Resources](#)
 - [6.1 References](#)



Drucker Diagnostics 00-276-009-007 Centrifuge with Fixed Angle Rotor



A centrifuge for every animal health need.

Specifications

- **Brand:** TrueBond
- **Models:**
 - TrueBond Set 3 00-276-009-007
 - TrueBond 00-376-009-005
 - TrueBond Duet 00-056-009-000
 - TrueBond TriFLEX Triac Replacement
- **Rotor Types:** Fixed Angle, Dual Rotor, Microhematocrit
- **Interface:** 3-setting, 10-setting digital, 5-setting
- **Max Speed RPM:** Varies by model
- **Max RCF:** Varies by model
- **Dimensions:** Varies by model

TrueBond Centrifuges Multipurpose Veterinary Processing

Model	Rotor Type	Interface	Preset Cycles Time / Speed / RCF	Max Speed RPM / RCF	Up to 16 x 100 mm or 10 mL	Up to 16 x 125 mm or 15 mL	Microhematocrit	Dimensions (W x D x H)
TrueBond Set 3 00-276-009-007	Fixed Angle	3-Setting	10 min / 3900 RPM / 1650 xg 5 min / 2000 RPM / 425 xg 6 min / 3800 RPM / 1550 xg	3900 / 1850	6	6	8	12 x 14 x 9 in 30 x 36 x 23 cm
TrueBond 00-376-009-005	Fixed Angle	10-Setting, Digital	10 min / 3900 RPM / 1650 xg 5 min / 2000 RPM / 500 xg 5 min / 1600 RPM / 300 xg 6 min / 3800 RPM / 1550 xg	3900 / 1850	6	6	8	12 x 14 x 9 in 30 x 36 x 23 cm

(adapters available for a variety of tube sizes)

This centrifuge is suitable for processing up to 16 x 100mm or 10 mL tubes. It has preset cycles for different speeds and relative centrifugal forces (RCF).

TrueBond 00-376-009-005 Centrifuge



Model	Rotor Type	Interface	Preset Cycles Time / Speed / RCF	Max Speed RPM / RCF	1.5 - 2 mL Fixed Angle Rotor	Microhematocrit MHCT Rotor	Dimensions (W x D x H)
TrueBond Duet 00-056-009-000	Dual Rotor	5 Setting	30 sec / 15800 RPM / 12000 xg 1.5 min / 15800 RPM / 12000 xg 2.5 min / 15800 RPM / 12000 xg 45 sec / 9800 RPM / 3900 xg 2 min / 15800 RPM / 13700 xg	15800 / 13700 (microhematocrit) 15800 / 12000 (1.5 mL - 2 mL)	2	12	7 x 7 x 7 in 18 x 18 x 18 cm

This centrifuge is designed for up to 16 x 125mm or 15 mL tubes. It features a fixed-angle rotor and a 10-setting digital interface for precise control.

TrueBond Duet

Compact Dual Rotor Centrifuge

This centrifuge is ideal for smaller volumes, with a focus on 1.5 – 2mL tubes. It offers dual rotors and preset cycles for efficient processing.

TrueBond TriFLEX



Triac Replacement Centrifuge

Model	Rotor Type	Interface	Preset Cycles Time / Speed / RCF	Max Speed RPM / RCF	Up to 16 x 100 mm or 10 mL	Up to 16 x 125 mm or 15 mL	Microhematocrit	Dimensions (W x D x H)
TriFLEX In Development	Combination	10 Setting, Digital	5 min / 3500 RPM / 1500 xg 6 min / 2300 RPM / 600 xg 5 min / 1571 RPM / 300 xg 3 min / 10400 RPM / 12600 xg	10400 / 12600	8 Fixed Angle	8 Fixed Angle	12 Horizontal	15 x 18 x 11 in 38 x 46 x 28 cm

This centrifuge is versatile, accommodating various tube sizes with its combination of rotor types. It features a 10-setting digital interface for customizable cycles.

(adapters available for a variety of tube sizes)

All centrifuges are available for a free 30-day evaluation.

Learn more at druckerdiagnostics.com/evaluation-program



+1-[866-265-1486](tel:866-265-1486) (U.S. only)
+1-[814-692-7661](tel:814-692-7661)

druckerdiagnostics.com
email to: sales@druckerdiagnostics.com

Designed, built, and supported in the USA

FAQs

Can I use different tube sizes in these centrifuges?

Yes, adapters are available for a variety of tube sizes to ensure compatibility with different sample containers.


How can I evaluate a centrifuge before purchasing?

All centrifuges are available for a free 30-day evaluation. Visit druckerdiagnostics.com/evaluation-program for more information.

How can I contact customer support?

You can reach customer support in the USA at +1-[866-265-1486](tel:866-265-1486) (U.S. only) or at +1-[814-692-7661](tel:814-692-7661). For email inquiries, contact sales@druckerdiagnostics.com.

Documents / Resources



[Drucker Diagnostics 00-276-009-007 Centrifuge with Fixed Angle Rotor](#) [pdf] User Guide 00-276-009-007, 00-376-009-005, 00-056-009-000, 00-276-009-007 Centrifuge with Fixed Angle Rotor, 00-276-009-007, Centrifuge with Fixed Angle Rotor, with Fixed Angle Rotor, Fixed Angle Rotor, Angle Rotor

References

- [Evaluation Program - Drucker Diagnostics](#)
- [User Manual](#)

Manuals+. [Privacy Policy](#)

This website is an independent publication and is neither affiliated with nor endorsed by any of the trademark owners. The "Bluetooth®" word mark and logos are registered trademarks owned by Bluetooth SIG, Inc. The "Wi-Fi®" word mark and logos are registered trademarks owned by the Wi-Fi Alliance. Any use of these marks on this website does not imply any affiliation with or endorsement.