



dynamic BIOSENSORS 10X BUFFER TE140 PH 7.4 Running Buffer User Manual

[Home](#) » [dynamic BIOSENSORS](#) » [dynamic BIOSENSORS 10X BUFFER TE140 PH 7.4 Running Buffer User Manual](#) 

Contents

- [1 dynamic BIOSENSORS 10X BUFFER TE140 PH 7.4 Running Buffer](#)
- [2 Product Information](#)
- [3 Product Usage Instructions](#)
- [4 Product Description](#)
- [5 Preparation](#)
- [6 Contact](#)
- [7 Documents / Resources](#)
 - [7.1 References](#)
- [8 Related Posts](#)

dynamic
BIOSENSORS

dynamic BIOSENSORS 10X BUFFER TE140 PH 7.4 Running Buffer



Product Information

Specifications

- **Product Name:** heliX+ 10X BUFFER TE140 PH 7.4
- **Manufacturer:** Dynamic Biosensors GmbH & Inc.
- **Order Number:** BU-TE-140-10
- **Composition:** Running buffer
- **Amount:** 50 mL
- **Storage:** For research use only. Limited shelf life. Check expiry date on label.

Product Usage Instructions

1. Ensure the product is labeled BU-TE-140-10 v2.1.
2. Store the product in the recommended conditions specified on the label.
3. Refer to the manufacturer's website for any updates or additional information.

Contact Information

Dynamic Biosensors GmbH: Perchtinger Str. 8/10,
81379 Munich, Germany

Dynamic Biosensors, Inc.: 300 Trade Center,
Suite 1400, Woburn, MA 01801, USA

Order Information:

order@dynamic-biosensors.com

Technical Support:

support@dynamic-biosensors.com

FAQ

What is the shelf life of heliX+ 10X BUFFER TE140 PH 7.4?

The product has a limited shelf life. Please check the expiry date on the label for specific information.

Can this product be used for purposes other than research?

No, this product is intended for research use only.

10X BUFFER TE140 PH 7.4

running buffer

Dynamic Biosensors GmbH & Inc.

BU-TE-140-10 v2.1

Product Description

Order Number: BU-TE-140-10

Table 1. Contents and Storage Information

Material	Composition	Amount	Storage
10x Buffer TE140 pH 7.4	100 mM Tris, 1.4 M NaCl, 500 µM EDTA, 500 µM EGTA and 0.5 % Tween20; 0.2 µm sterile filtered	50 mL	2-8°C

For research use only.

This product has a limited shelf life, please see expiry date on label.

Preparation

Dilute the complete solution 10x Buffer TE140 pH 7.4 (50 mL) by mixing with 450 mL ultrapure water.

After dilution TE140 Buffer is ready for use (10 mM Tris, 140 mM NaCl, 50 µM EDTA, 50 µM EGTA and 0.05 % Tween20).

The diluted buffer should be stored at 2-8°C.

Contact

Dynamic Biosensors GmbH

Perchtinger Str. 8/10

81379 Munich

Germany

Dynamic Biosensors, Inc.
300 Trade Center, Suite 1400
Woburn, MA 01801
USA

Order Information order@dynamic-biosensors.com
Technical Support support@dynamic-biosensors.com
www.dynamic-biosensors.com

Instruments and chips are engineered and manufactured in Germany.
©2024 Dynamic Biosensors GmbH | Dynamic Biosensors, Inc. All rights reserved.
www.dynamic-biosensors.com

Documents / Resources

	dynamic BIOSENSORS 10X BUFFER TE140 PH 7.4 Running Buffer [pdf] User Manual BU-TE-140-10, 10X BUFFER TE140 PH 7.4 Running Buffer, 10X BUFFER TE140 PH 7.4, Running Buffer, Buffer
---	--

References

- [📄 HomePage | Biosensors International Ltd](#)
- [🔍 Home - Dynamic Biosensors](#)
- [🔍 Home - Dynamic Biosensors](#)
- [📖 User Manual](#)

[Manuals+](#), [Privacy Policy](#)

This website is an independent publication and is neither affiliated with nor endorsed by any of the trademark owners. The "Bluetooth®" word mark and logos are registered trademarks owned by Bluetooth SIG, Inc. The "Wi-Fi®" word mark and logos are registered trademarks owned by the Wi-Fi Alliance. Any use of these marks on this website does not imply any affiliation with or endorsement.