



dynamic BIOSENSORS BU-HE-40-10 v2.1 Running Buffer User Manual

[Home](#) » [dynamic BIOSENSORS](#) » dynamic BIOSENSORS BU-HE-40-10 v2.1 Running Buffer User Manual 

dynamic BIOSENSORS BU-HE-40-10 v2.1 Running Buffer



Contents

- [1 Product Description](#)
- [2 Preparation](#)
- [3 Contact](#)
- [4 CUSTOMER SUPPORT](#)
- [5 Documents / Resources](#)
 - [5.1 References](#)
- [6 Related Posts](#)

Product Description

Order Number: **BU-HE-40-10**

Table 1. Contents and Storage Information

Material	Composition	Amount	Storage
10x Buffer HE40 pH 7.4	100 mM HEPES, 400 mM NaCl, 500 µM EDTA, 500 µM EGTA and 0.5 % Tween 20; 0.2 µm sterile filtered	50 mL	2-8°C

For research use only.

This product has a limited shelf life, please see expiry date on label.

Preparation

Dilute the complete solution 10x Buffer HE40 pH 7.4 (50 mL) by mixing with 450 mL ultrapure water.

After dilution HE40 Buffer is ready for use (10 mM HEPES, 40 mM NaCl, 50 µM EDTA, 50 µM EGTA and 0.05 % Tween 20).

The diluted buffer should be stored at 2-8°C.

Contact

Dynamic Biosensors GmbH	Dynamic Biosensors, Inc.
Perchtinger Str. 8/10	300 Trade Center, Suite 1400
81379 Munich	Woburn, MA 01801
Germany	USA

Order Information order@dynamic-biosensors.com

Technical Support support@dynamic-biosensors.com

CUSTOMER SUPPORT

www.dynamic-biosensors.com

Instruments and chips are engineered and manufactured in Germany.

©2024 Dynamic Biosensors GmbH | Dynamic Biosensors, Inc. All rights reserved.



Documents / Resources



[dynamic BIOSENSORS BU-HE-40-10 v2.1 Running Buffer](#) [pdf] User Manual
BU-HE-40-10, BU-HE-40-10 v2.1 Running Buffer, BU-HE-40-10 v2.1, Running Buffer, Buffer

References

- [HomePage | Biosensors International Ltd](#)
- [Home - Dynamic Biosensors](#)
- [Home - Dynamic Biosensors](#)
- [User Manual](#)

[Manuals+](#), [Privacy Policy](#)

This website is an independent publication and is neither affiliated with nor endorsed by any of the trademark owners. The "Bluetooth®" word mark and logos are registered trademarks owned by Bluetooth SIG, Inc. The "Wi-Fi®" word mark and logos are registered trademarks owned by the Wi-Fi Alliance. Any use of these marks on this website does not imply any affiliation with or endorsement.