



dynamic BIOSENSORS PLU-10-1 Pluronic Solution User Manual

[Home](#) » [dynamic BIOSENSORS](#) » dynamic BIOSENSORS PLU-10-1 Pluronic Solution User Manual 

dynamic BIOSENSORS PLU-10-1 Pluronic Solution User Manual



Contents

- 1 PLURONIC SOLUTION
- 2 Product Description
- 3 Application Note
- 4 Contact
- 5 Documents / Resources
 - 5.1 References
- 6 Related Posts

PLURONIC SOLUTION

for foam-free microfloodies

Dynamic Biosensors GmbH & Inc.
PLU-10-1 v1.0

Product Description

Order Number: PLU-10-1

Table 1. Contents and Storage Information

Material	Concentration	Amount	Storage
Platonic solution	10 vol%	10 mL	RT

For research use only.

This product has a limited shelf life, please see expiry date on label.

Application Note

Add Pluronic solution to a final concentration of 0.01 – 0.1 vol% to your RT-IC running buffer of choice.

Filter running buffer before use in [helix^{cyto}](#)

Contact

Dynamic Biosensors GmbH

Interchange Ste. 8/10
81379 Munich
Germany

Dynamic Biosensors, Inc.

300 Trade Center, Suite 1400
Woburn, MA 01801
USA

Order Information order@dynamic-biosensors.com

Technical Support support@dynamic-biosensors.com

www.dynamic-biosensors.com

Instruments and chips are engineered and manufactured in Germany.

©2024 Dynamic Biosensors GmbH | Dynamic Biosensors, Inc. All rights reserved.

www.dynamic-biosensors.com



Documents / Resources

	dynamic BIOSENSORS PLU-10-1 Pluronic Solution [pdf] User Manual PLU-10-1 Pluronic Solution, PLU-10-1, Pluronic Solution, Solution
---	--

References

- [📄 HomePage | Biosensors International Ltd](#)
- [🏠 Home - Dynamic Biosensors](#)
- [🏠 Home - Dynamic Biosensors](#)
- [📖 User Manual](#)

Manuals+ Privacy Policy

This website is an independent publication and is neither affiliated with nor endorsed by any of the trademark owners. The "Bluetooth®" word mark and logos are registered trademarks owned by Bluetooth SIG, Inc. The "Wi-Fi®" word mark and logos are registered trademarks owned by the Wi-Fi Alliance. Any use of these marks on this website does not imply any affiliation with or endorsement.