



Educational Insights EI-8801 GeoSafari Jr. Talking Microscope User Manual

[Home](#) » [Support](#) » Educational Insights EI-8801 GeoSafari Jr. Talking Microscope User Manual 

Contents

- [1 Educational Insights EI-8801 GeoSafari Jr. Talking Microscope](#)
- [2 Parts Description](#)
- [3 How to Use](#)
- [4 Facts](#)
- [5 Quiz](#)
- [6 Battery Instructions](#)
- [7 Cleaning Instructions](#)
- [8 Troubleshooting](#)
- [9 Changing Languages](#)
- [10 Frequently Asked Questions](#)
- [11 Video- Usage Instructions](#)
- [12 References](#)
- [13 Related Posts](#)



Educational Insights EI-8801 GeoSafari Jr. Talking Microscope



WARNING: CHOKING HAZARD – Small parts. Not for children under three (3) years.

Requires 3 AAA batteries (not included)



Parts Description



How to Use

Here's How to Use the Talking Microscope

1. Install 3 AAA batteries, following the instructions on page 3.


2. Press  to turn on the Talking Microscope and select your volume by pressing .

3. Pick a mode: Press  for Facts. Press  for Quiz.

4. Insert a slide into the microscope stage as shown below. If you place the slide correctly, you will see the specimen in the center of the viewing area when looking through the eyepiece. Your microscope will also begin reading facts or questions aloud. There are 20 slides in total; each slide contains three images.



1. Slide going in
2. Slide in place
3. Slide coming out

5. Turn the Talking Microscope off by pressing the  button.

Facts

- In fact mode, you will hear the Talking Microscope say fun facts about the slide you have chosen. Look into the microscope as you listen. Then gently pull the slide out as shown in the diagram above, shift it over to the next image, and click it into place. Each slide has a total of five facts. Each individual image may have one, two, or three facts.

Quiz

In quiz mode, the Talking Microscope will ask questions. You can enter your answer by pressing an answer button. The Talking Microscope will tell you if the answer was correct. When you adjust the slide to see the next image, you'll hear new questions. Each slide has a total of five questions. Each individual image may have one, two, or three questions.

- Some questions are True or False.
- Press 1 for True ✓. Press 2 for False ✗.
- Some questions are multiple-choice.
- Press 1, 2, or 3 to answer the question. Look into the microscope to see if the image will help you to answer the question.

Tips


- Hold your slides on the sides. Smudges will make the specimens harder to see.
- To clean the slides, wipe them gently with a damp cloth. Then store them in the drawer at the back of the microscope.
- Forget to turn it off? The Talking Microscope will power down on its own after 2 minutes.

Find more things to look at!

- The Talking Microscope may also be used to view other microscope slides. Simply place the slide on the viewing area and push for light.

Battery Instructions

1. Use a screwdriver to carefully remove the screw on the battery door.
2. Open the battery door.
3. Install 3 AAA batteries as shown in the diagram.
4. Close the battery door.
5. Use a screwdriver to secure the battery door.



Note:

- Do not mix old and new batteries.
- Do not mix different types of batteries: alkaline, standard (carbon zinc), or rechargeable (nickel-cadmium) batteries.
- Do not recharge non-rechargeable batteries.
- Remove rechargeable batteries from the toy before charging.
- Only charge rechargeable batteries under adult supervision.
- Only use batteries of the same or equivalent type as recommended.
- Insert batteries with the correct polarity.
- Remove exhausted batteries from the toy.
- Do not short-circuit the supply terminals.

Cleaning Instructions

Clean product with a damp or dry cloth. Do not immerse or spray water or other liquids on the product.

Troubleshooting

If the microscope should malfunction, try removing and reinserting the slide or replacing the batteries.

Changing Languages

1. Follow the Battery Instructions to carefully open the battery door.
2. Locate the language button inside the battery compartment as shown in the diagram on page 3.
3. Press the button once for English, twice for Spanish (Español), three times for French (Francais), or four times for German (Deutsche). The Talking Microscope will say the name of the selected language after each press to confirm your choice.
4. Follow the Battery Instructions to close the battery door.

Frequently Asked Questions

What age range is the Educational Insights EI-8801 GeoSafari Jr. Talking Microscope suitable for?

The Educational Insights EI-8801 GeoSafari Jr. Talking Microscope is suitable for children above the age of three (3) years.

What type of batteries does the Educational Insights EI-8801 GeoSafari Jr. Talking Microscope require?

The Educational Insights EI-8801 GeoSafari Jr. Talking Microscope requires 3 AAA batteries (not included).

How do I turn on the Educational Insights EI-8801 GeoSafari Jr. Talking Microscope?

To turn on the Talking Microscope, press the power button.

How many modes does the Educational Insights EI-8801 GeoSafari Jr. Talking Microscope have, and how do I switch between them?

The Educational Insights EI-8801 GeoSafari Jr. Talking Microscope has two modes: Fact mode and Quiz mode. You can switch between them by pressing the corresponding buttons on the microscope.

How many slides are included with the Educational Insights EI-8801 GeoSafari Jr. Talking Microscope, and how many images does each slide contain?

There are 20 slides included with the Educational Insights EI-8801 GeoSafari Jr. Talking Microscope. Each slide contains three images.

Can I clean the slides of the Educational Insights EI-8801 GeoSafari Jr. Talking Microscope, and if so, how?

Yes, you can clean the slides by wiping them gently with a damp cloth. Then, store them in the drawer at the

back of the microscope.

What should I do if I forget to turn off the Educational Insights EI-8801 GeoSafari Jr. Talking Microscope?

The Educational Insights EI-8801 GeoSafari Jr. Talking Microscope will power down on its own after 2 minutes.

Can I use other microscope slides with the Educational Insights EI-8801 GeoSafari Jr. Talking Microscope?

Yes, you can use other microscope slides with the Educational Insights EI-8801 GeoSafari Jr. Talking Microscope by placing them on the viewing area and pushing for light.

How do I change the language settings of the Educational Insights EI-8801 GeoSafari Jr. Talking Microscope?

To change the language settings, follow the instructions provided in the user manual.

What should I do if the Educational Insights EI-8801 GeoSafari Jr. Talking Microscope malfunctions?

If the microscope malfunctions, try removing and reinserting the slide or replacing the batteries.

Can the Educational Insights EI-8801 GeoSafari Jr. Talking Microscope be used without batteries?

No, the Educational Insights EI-8801 GeoSafari Jr. Talking Microscope requires 3 AAA batteries to operate.

Is the Educational Insights EI-8801 GeoSafari Jr. Talking Microscope suitable for outdoor use?

The Educational Insights EI-8801 GeoSafari Jr. Talking Microscope is designed for indoor use and may not be suitable for outdoor environments.

Video- Usage Instructions



[Download The PDF User Manual- Educational Insights EI-8801 GeoSafari Jr. Talking Microscope User Manual](#)
[GeoSafari-Jr.-Talking-Microscope-video.mp4](#)

References

- [User Manual](#)