



Home » Excelitas » Excelitas XYLISTM II Broad Spectrum LED Illumination System Owner's Manual 75

Contents [hide]

- 1 Excelitas XYLISTM II Broad Spectrum LED Illumination System
- **2 TECHNICAL SPECIFICATIONS**
- 3 Flexibility to Suit Application Needs
- 4 Exceptional Control with Built-in Options
- 5 Potential Cost and Energy Savings with LEDS
- 6 X-Cite XYLIS II Spectra
- 7 Documents / Resources
 - 7.1 References



Excelitas XYLISTM II Broad Spectrum LED Illumination System



TECHNICAL SPECIFICATIONS

	Model: XT730S	Model: XT730L	
Wavelength Range	360-770nm	380-770nm	
LED Peaks (nm)	365, 430, 475, 545, 635, 735	385, 430, 475, 545, 635, 7 35	
Input Power Supply	Universal input 100-240VAC, 50/60Hz		
Current	2.4 – 1.0a		
LED ON/OFF Response Times	100 μs TTL / 1 ms USB		
Control Options	speedDIAL ON/OFF – TTL compatible RS-232 commands (SDK available), USB Foot Pedal (optional)		
I/O Connections	BNC input USB (B-type) 3.5mm stereo plug		
Dimensions (W x H x D)	130mm x 250mm x 260mm (5.1" x 9.8" x 10.3")		
Weight	5.2 kg (11.6 lbs)		
Certifications	CE, Rohs, UKCA, KC, PSE compliant		
Warranty	LEDS: 25,000 hours or 3 years All other X-Cite XYLIS II components: 1 year, parts and lab our (excluding LLG)		

Patents	X-Cite XYLIS II incorporates technology protected by	
	patent US#9,239,133	

Better Performance, Less Noise

Introducing X-Cite XYLIS™ II. With improvements across its broad spectrum and sound levels 10 db quieter than the original XYLIS, the only thing you'll notice is great performance. Researchers can enjoy the benefits of LED technology without compromising on price, flexibility, or performance. No more hesitation, no more excuses.

Brightness of an Arc Lamp

Specially selected LEDS built into the X-Cite XYLIS II are powerful enough to replace arc lamps on both compound and stereomicroscopes. The impressive output and low maintenance of XYLIS II can help breathe new life into underused microscopes and make better use of laboratory resources.

Broad Spectral Coverage

Like the original XYLIS, XYLIS II covers an extended spectral range, providing excitation from DAPI to Cy7. Spectral highlights include:

- DAPI: Two X-Cite XYLIS II models are available to provide a choice of UV excitation.
 XT730S has a 365nm LED for a closer match to arc lamp output and compatibility with the narrow 365 DAPI filter sets which come standard in most microscopes. XT730L has a 385nm LED for use with sensitive specimens and 385 DAPI filter sets which are becoming increasingly common.
- TRITC/Tx Rd/mCherry: X-Cite XYLIS II incorporates Excelitas' patented and award-winning Laserled Hybrid Drive® technology, utilising high-efficiency lasers to excite a phosphor layer and generate light from 500nm to 600nm. The resulting intense, broad peak ensures plenty of power in this critical part of the spectrum.
- Cy7: X-Cite XYLIS II includes a 735nm peak for Cy7 excitation. Labs no longer have to choose between the benefits of LEDS and keeping their spectral options open.

Flexibility to Suit Application Needs

In addition to its powerful output and broad DAPI to Cy7 spectral range, X-Cite XYLIS

Il offers the ultimate in flexibility. Delivering light through a light guide, XYLIS II can be used with more than a dozen microscope adaptors or, where the microscope design allows, with the light guide alone.

- X-Cite XYLIS II can be installed on just about any new imaging system or used to
 retrofit the microscopes that labs have depended on for years. Offered in two models
 with a choice of UV wavelengths (365nm or 385nm), labs may choose the one that is
 suitable for their preferred or existing DAPI filter sets.
- "UV Off" mode disables the UV LED when it is not required at all, protecting sensitive samples and extending light guide lifetime. The whisper-quiet operation of the XYLIS II makes it ideal for use in multi-station FISH labs and any imaging facility wishing to maintain a peaceful environment for conducting research.

Exceptional Control with Built-in Options

When it comes to control options, all are built into every

- X-Cite XYLIS II each system includes manual fingertip control with speed dial, hands-free operation with a foot pedal, as well as USB and TTL inputs for automated applications.
- X-Cite XYLIS II's ergonomic speedDIAL is designed with a large speed-sensitive
 intensity dial that doubles as an ON/OFF button, allowing for quick and intuitive control
 of illumination. With a simple double tap on speed dial, users can quickly jump to a
 favourite intensity setting.
- Take full advantage of LED instant ON/OFF capability to limit photobleaching and phototoxicity with ultra-fast PC control or TTL triggering. X-Cite XYLIS II can be driven by commercial imaging software, and an SDK is available to develop customised control solutions.

X-Cite Costs & Energy Savings

X-Cite XYLIS II allows researchers to reduce the amount of hardware required by an imaging system – replace an arc lamp, separate shutter and neutral density filters – all with a single device. All systems include high speed shuttering, 1% power adjustment, and multiple manual/automated control options. Enjoy the long-term benefits of using LED technology, including long lifetimes, consistent output, lower maintenance, fewer consumables, reduced energy use, and zero mercury waste.

Potential Cost and Energy Savings with LEDS

Table 1: Cost of Ownership (per 20,000 hours of "ON time")

	НВО	X-Cite 120Q	X-Cite XYLIS II
Replacement Lamps	100	10	_
Mercury Content	1100 mg	200 mg	_
Lamp Costs	\$15,000	\$7,000	_
Replacement Light Guid es	_	5	2
Light Guide Costs	_	\$2,225	\$890
Bulb Disposal (\$5/bulb)1	\$500	\$50	_
Maintenance Costs (bulb, \$20/hr)2	\$1,000	\$17	_
TOTAL	\$16,500	\$9,292	\$890
Hourly Cost3	\$0.825	\$0.465	\$0.045
Annual Cost4	\$1,650	\$929	\$11

Table 2: Energy Consumption (per day)

НВО	X-Cite 120Q	X-Cite XYLIS II
-----	-------------	-----------------

Electricity Cost (per year) 4	\$46.50	\$72.00	\$9.00
Electricity Cost (per day) (\$ 0.15/kWh)6	\$0.19	\$0.29	\$0.04
Energy Used5	1.24 kWh	1.92 kWh	0.24 kWh
Lamp Lifetime Used	4%	0.4%	0.004%
ON Time3	8 Hours	8 Hours	1 Hour

Notes

- 1. Mercury-Free Microscopy white paper: www.mygreenlab.org.
- 2. Assumes 30 min to change/align HBO lamp, 5 min for X-Cite 120q.
- 3. Assumes 8 8-hour days, 4 x 15-minute imaging sessions. Arc lamps are left on for the day, and LEDS are on continuously during each session.
- 4. Assumes a 5-day week x 50 weeks.
- 5. Calculated based on published technical specifications.
- 6. Typical rate. Actual rates will vary by region and/or time of day.

With instant ON/OFF capabilities, intense output and whisper-quiet operation, X-Cite XYLIS II is the ultimate partner to advance your fluorescence imaging research.

X-Cite XYLIS II Spectra

X-Cite XYLIS II Spectra with Common Excitation Bands FITC TXRD Cv7 0.8 X-Cite XYLIS II XT730S (365) Relative Intensity 0.6 ... X-Cite XYLIS II XT730L (385) 0.2 0 300 350 400 450 550 600 650 700 750 800

Filter bands shown: 377/50, 470/40, 545/25, 560/40, 620/60, 710/70.

Wavelength (nm)

- www.excelitas.com
- x-cite@excelitas.com
- Telephone: +1 905 821-2600
- Toll Free (USA and CAN): +1 800 668-8752 Fax: +1 905 821-2055

For a complete listing of our global offices, visit www.excelitas.com/locations. © 2023 Excelitas Technologies Corp. X-Cite®LaserledD Hybrid Drive® is are registered trademarks, and Xylis™ is a trademark of Excelitas Technologies Corp. All rights reserved. The Excelitas logo and design are registered trademarks of Excelitas Technologies Corp. All other trademarks not owned by Excelitas Technologies or its subsidiaries that are depicted herein are the property of their respective owners. Excelitas reserves the right to change this document at any time without notice and disclaims liability for editorial, pictorial typographical errors.

Documents / Resources



Excelitas XYLISTM II Broad Spectrum LED Illumination System [pdf] Own er's Manual

XYLISTM II, 110LED, TETREMTM, NOVEMTM, XLED1, XYLISTM II Broad Spectrum LED IIIumination System, XYLISTM II, Broad Spectrum LED II Iumination System, LED IIIumination System, Illumination System, System

References

- ► Homepage | Excelitas
- User Manual
- Excelitas
- ♦ 110LED, Broad Spectrum LED Illumination System, Excelitas, Illumination System, LED Illumination System, NOVEMTM, System, TETREMTM, XLED1, XYLISTM II, XYLISTM II Broad Spectrum LED Illumination System
 - —Previous Post

excelitas TETREM X-Cite LED Light Source for Fluorescence Microscopy Owner's Manual

Leave a comment

Your email address will not be published. Required fields are marked *

Save my name, email, and website in this browser for the next time I comment.

Comment *	
Name	
Email	
Website	

Search:

Post Comment

e.g. whirlpool wrf535swhz

Search

Manuals+ | Upload | Deep Search | Privacy Policy | @manuals.plus | YouTube

This website is an independent publication and is neither affiliated with nor endorsed by any of the trademark owners. The "Bluetooth®" word mark and logos are registered trademarks owned by Bluetooth SIG, Inc. The "Wi-Fi®" word mark and logos are registered trademarks owned by the Wi-Fi Alliance. Any use of these marks on this website does not imply any affiliation with or endorsement.