

# FreeStyle Libre 3 Continuous Glucose System Instructions


Home » FreeStyle » FreeStyle Libre 3 Continuous Glucose System Instructions

# **Contents**

- 1 FreeStyle Libre 3 Continuous Glucose System
- 2 Specifications
- **3 Product Usage Instructions**
- **4 MOBILE DEVICE & OS COMPATIBILITY**
- 5 FAQs
- 6 Documents / Resources
  - **6.1 References**
- **7 Related Posts**



FreeStyle Libre 3 Continuous Glucose System



# **Specifications**

Product Name: FreeStyle Libre 3
Versions: 3.4, 3.5, 3.5.1, 3.6.1

# **Product Usage Instructions**

#### Compatibility

 The FreeStyle Libre 3 is compatible with select mobile devices and operating systems. Please refer to the list below for compatibility details:

#### iOS Devices

- iOS Versions: 15.6, 15.6.1, 15.7, 16, 17, 17.0.1, 17.0.2, 17.0.3, 17.1, 17.1.1, 17.1.2, 17.2
- Compatible Devices: iPhone 7 and newer models including SE (2020), SE (2022)

# Android Devices

- Android Versions: 10\*, 11, 12, 13, 14
- Compatible Devices: Various Samsung Galaxy models, Google Pixel devices, Motorola phones, Sony Xperia, Xiaomi, Huawei, and more

# • NFC Scan Location

• Refer to the diagram provided in the user manual to locate the NFC scan area on your mobile device.

#### App Installation


Before installing a new OS version on your phone or using the App with a new phone, it is recommended
to check the compatibility guide provided.

#### Important Note

• Beta OS versions are not supported. Jailbroken and rooted devices are also not compatible with the App.

# **MOBILE DEVICE & OS COMPATIBILITY**

Popular mobile devices and operating systems (OS) are regularly tested to evaluate NFC scan performance, Bluetooth connectivity, and app compatibility with Sensors. Refer to the diagram for the NFC scan location on your mobile device. We recommend checking this guide before installing a new OS version on your phone or before using the App with a new phone.



Compatible devices and operating systems

# FreeStyle Libre 3 (version 3.4, 3.5, 3.5.1, 3.6.1)

os	DEVICE	NFC LOCATION
	iPhone 7, 7 Plus, 8, 8 Plus, X, XR,	
	XS, XS Max, 11, 11 Pro, 11 Pro	
iOS: 15.6, 15.6.1, 15.7, 16, 17,	Max, 12, 12 Mini, 12 Pro, 12 Pro	
17.0.1, 17.0.2, 17.0.3, 17.1, 17.1.1,	Max, SE (2020), 13, 13 Mini, 13	1
17.1.2, 17.2	Pro, 13 Pro Max, 14, 14 Plus, 14	
	Pro, 14 Pro Max, 15, 15 Plus, 15	
	Pro, 15 Pro Max, SE (2022)	

FreeStyle Libre 3 (version 3.4.2, 3.5.1, 3.6.0)

os	DEVICE	NFC LOCATION
Android: 10*, 11, 12, 13, 14	Samsung Galaxy A13 5G, A15 5G, A22, A 22 5G, A23, A32, A32 5G, S10 Lite, Note1 0 Lite	1
	Dialn G65	
	Google Pixel 4XL	
	Motorola moto g60s	
	Sony Xperia 1 III	
	Xiaomi 11T, Redmi Note 10 5G	
	Samsung Galaxy S9, S10, S10+, S10e, S2 0 FE 5G, S21 5G, S21+ 5G, S21 FE 5G, S 22, S22+, S22 Ultra, S23, S23 Ultra, S24 U ltra, Note9, Note10, Note10+, Note20 5G, Note20 Ultra 5G, A20, A21, A42 5G, A52s 5G, Z Fold2 5G, Z Fold3 5G	2
	Google Pixel 4, 4a, 4a 5G, 5, 5a, 6, 6a, 6 P ro, 7, 7 Pro, 8 Pro	
	Huawei P30	
	Motorola moto g51 5G, one 5G ACE	
	Xiaomi Mi Note 10 Lite, 11 Lite 5G NE, 11T Pro, Poco F3	
	Samsung Galaxy S9+, S20+ 5G, S20 Ultra 5G, S21 Ultra 5G, Note10+ 5G, Z Flip3 5G, Z Flip4, Z Fold4	3

- This table will be updated as other devices and operating systems are evaluated.
- The contents of this guide are the most up to date and may be different from what is displayed within the FreeStyle Libre 3 app. Contact Customer Service for more information on devices and operating systems other than those listed.
- Beta OS versions are not evaluated or supported. The App cannot be used with jailbroken and rooted devices.
- Please check that your mobile device and OS continues to be supported by the manufacturer. Use of an unsupported mobile device or OS version may affect App security and functionality.
- You are responsible for the risks associated with the use of the App on an unsupported mobile device or OS version.

The sensor housing, FreeStyle, Libre, and related brand marks are marks of Abbott. Other trademarks are the property of their respective owners. ©2024 Abbott ART44628-004 Rev. S 09/2

#### **FAQs**

- Q: Can I use the App with any mobile device?
  - A: The App is only compatible with the listed devices and operating systems. Unsupported devices may compromise security and functionality.

- Q: How do I find the NFC scan location on my phone?
  - A: Refer to the diagram in the user manual or contact Customer Service for assistance.

# **Documents / Resources**



FreeStyle Libre 3 Continuous Glucose System [pdf] Instructions 3.4, 3.5, 3.5.1, 3.6.1, 3.4.2, 3.5.1, 3.6.0, Libre 3 Continuous Glucose System, Libre 3, Continuous Glucose System, Glucose System, System

#### References

User Manual

Manuals+, Privacy Policy

This website is an independent publication and is neither affiliated with nor endorsed by any of the trademark owners. The "Bluetooth®" word mark and logos are registered trademarks owned by Bluetooth SIG, Inc. The "Wi-Fi®" word mark and logos are registered trademarks owned by the Wi-Fi Alliance. Any use of these marks on this website does not imply any affiliation with or endorsement.