



## grafco 3618A Funnel Owner's Manual

[Home](#) » [grafco](#) » grafco 3618A Funnel Owner's Manual 

### Contents

- [1 grafco 3618A Funnel](#)
- [2 Features](#)
- [3 Sterilization Guidelines](#)
- [4 Overall Capacity](#)
- [5 Top Diameter](#)
- [6 ABOUT COMPANY](#)
- [7 Documents / Resources](#)
  - [7.1 References](#)
- [8 Related Posts](#)



**grafco 3618A Funnel**



## Features

- Made of linear polyethylene
- Stackable and dishwasher safe
- Proudly made in the USA

## Sterilization Guidelines

- **Boiling:** 212°F (100°C) 15 minutes maximum
  - **Gas:** Ethylene Oxide 86°F – 140°F (30°C – 60°C)
    - There should be no distortion with Ethylene Oxide Gas Sterilization temperature range.
    - Please allow ample time for the funnels to dissipate the absorbed gas after sterilization.
- Note:** Sterilization guidelines were provided by the raw material supplier, Dow Chemical Company.

## Overall Capacity

	oz	ml
3618A	2	59.1
3618C	8	236.6
3618D	16	473.2


## Top Diameter

inch      cm	UOM
2.75      7	1 ea
4.25      10.8	1 ea
5      12.7	1 ea

## ABOUT COMPANY

- GF Health Products, Inc. is an ISO 13485:2016 / MDSAP Certified Company
- +1 770.368.4700
- Information contained herein is subject to change without notice.
- The most current and complete product information can be found on our website.
- [www.grahamfield.com](http://www.grahamfield.com)
- © 2011, GF Health Products, Inc. All Rights Reserved.
- Grafcó is a registered trademark of GF Health Products, Inc.

## Documents / Resources

	<a href="#">grafco 3618A Funnel</a> [pdf] Owner's Manual 3618A Funnel, 3618A, Funnel
---	---

## References

- [Gf | Graham-Field | Home](#)