



# immedia IM426 SmartMove Support Pad Instruction Manual

[Home](#) » [immedia](#) » immedia IM426 SmartMove Support Pad Instruction Manual 

# immedia

by Etac

IM426 SmartMove Support Pad  
Instruction Manual



## Contents

- [1 IM426 SmartMove Support Pad](#)
- [2 Illustrations](#)
- [3 Documents / Resources](#)
  - [3.1 References](#)
- [4 Related Posts](#)

**IM426 SmartMove Support Pad**

Article no.	Description
IM426	SmartMove w/4 handles



#### Material:

**Outer surface:** Nylon, polyester, polyurethane

**Inner surface:** Polypropylene, polyethylene

**Handles:** Polyester.

Never try to repair a product yourself.

Never use a defective product.

#### Washing instructions:


This product can be wiped off with soapy water or cleaned with a solvent-free cleaning agent with a pH level between 5-9, or a 70% disinfection solution.

The product should be checked regularly, preferably each time it is used and especially after wash. Check that there is no damage to seams or fabric.

#### Illustrations



## Check



SmartMove  
<http://q.etac.com/003e>



Creating Possibilities

**Etac A/S**

Parallelsvej 3

DK-8751 Gedved, Denmark

[www.etac.com](http://www.etac.com)

IMM1027\_SI.VER9.2021-10-18

## Documents / Resources

	<p><a href="#">immedia IM426 SmartMove Support Pad</a> [pdf] Instruction Manual IM426 SmartMove Support Pad, IM426, SmartMove Support Pad, Support Pad, Pad</p>
---	---

## References

-  [Etac - A world-leading developer of assistive devices](#)