

# **LABCONCO 3970423 Purifier Vertical Clean Bench with Airflow Monitor Owner's Manual**

Home » LABCONCO » LABCONCO 3970423 Purifier Vertical Clean Bench with Airflow Monitor Owner's Manual



#### **Contents**

- 1 LABCONCO 3970423 Purifier Vertical Clean Bench with Airflow Monitor
- **2 Product Information**
- 3 Overview
- **4 Specifications**
- **5 Description**
- 6 FAQ
- 7 Documents / Resources
- 7.1 References
- **8 Related Posts**



LABCONCO 3970423 Purifier Vertical Clean Bench with Airflow Monitor



### **Product Information**

## **Specifications**

• Weight: 195.0 lbs (88.0 kg)

• **Dimensions**: 48" w x 28.6" d x 36.8" h

Electrical: 230V, 50/60 Hz, 3 A
Built-in Options: Airflow monitor

• Laminar Flow: Vertical

• Nominal Width:4'

Region: International

• Conformance: CE

• Power Cord & Plug: 230V, 10A, no plug

#### **Description**

The Purifier Vertical Clean Bench with Airflow Monitor is designed to provide a clean and controlled work environment for various laboratory applications. It comes in a 4' width size and features a built-in Guardian Airflow Monitor for continuous monitoring of downflow airflow.

#### **Features**

- · Guardian Airflow Monitor with LED lights
- 254 nm UV lamp with 3-way switch
- · Available in different widths
- International electrical configurations

# **Options**

- · Supporting work surface
- Adjustable Height Base Stands

#### Compliance

- UL Standard 3101-1/61010-1 (115 volt models)
- CAN/CSA Standard C22.2 No. 1010.1 conformity (115 volt models)
- CE Conformity marking (230-volt models)

#### **Contact Information**

If you need any product assistance, you can contact Labconco Corporation at the following address:

Labconco Corporation 8811 Prospect Avenue Kansas City, MO 64132-2696

• Phone: 816-333-8811

View online: <a href="https://www.labconco.com/product/4-purifier-vertical-clean-bench-with-airflow-monitor-2/1587">https://www.labconco.com/product/4-purifier-vertical-clean-bench-with-airflow-monitor-2/1587</a>

#### Overview

Purifier Vertical Clean Benches direct HEPA-filtered air downward over the work area to provide ISO Class 5 conditions and protect your work from particulate contamination. Applications for clean benches include pharmacy compounding, plant tissue culture, media plate preparation, electronics inspection, medical device assembly and more. Because they do not provide protection to the user, Vertical Clean Benches should not be used in conjunction with biohazardous material, toxins or radionuclides. Purifier Vertical Clean Benches are available in 2′, 3′, 4′ and 6′ widths. Labconco also offers Purifier Horizontal Clean Benches that direct HEPAfiltered air horizontally over the work area and Purifier Filtered PCR Enclosures, which are specialized vertical clean benches designed for polymerase chain reactions.

#### **Specifications**

• Weight: 195.0 lbs

• Weight metric: 88.0 kg

• **Dimensions**: 48" w x 28.6" d x 36.8" h

Electrical: 230V, 50/60 Hz, 3 A
Built-in Options: Airflow monitor

Laminar Flow: Vertical
Nominal Width: 4'
Region: International

• Conformance: CE

Power Cord & Plug: 230V, 10A, no plug

# **Description**

#### **Features**

- Nominal downflow velocity of 45-65 fpm
- Aerodynamic air foil
- Supply HEPA filter, 99.99% efficient on particles 0.3 micron
- ISO Class 5 air (formerly Class 100 per Federal Standard 209E)

- · Disposable prefilter
- Variable speed built-in blower with motorized impeller and solid state control rated for 10 amps
- Intrinsically-sound, negative pressure design
- · Minihelic pressure gauge
- Air diffuser
- UV-resistant, glacier white and gray, powder-coated aluminum and steel construction
- UV-opaque and resistant 1/4" thick tempered safety glass sash and sides
- Fluorescent lamp
- · Light and blower switches
- 9.4" high sash opening
- Ergonomic 10° angled and hinged sash that pivots upward, locks to a loading height of 20" and has wiping seals to promote ISO Class 5 conditions
- Two utility ports with iris openings and plugs
- ETL Listing (115 volt, 60 Hz models)
- One year warranty on parts and service

### **Options**

- Built-in Guardian Airflow Monitor that continuously monitor downflow airflow and has two LED lights to indicate when airflow is within or outside set point range.
- 254 nm UV lamp with 3-way switch (fluorescent, UV and off) and protection panel interlocked with the UV lamp that closes the sash opening
- 2', 3', 4' and 6' widths
- · International electrical configurations

#### **Required Accessories**

- Supporting work surface. See Solid Epoxy Dished Work Surfaces
- · Supporting base. See Adjustable Height Base Stands

#### **Optional Accessories**

- · General Accessories
- IV Bar Kit

# Compliance

- UL Standard 3101-1/61010-1 (115 volt models)
- CAN/CSA Standard C22.2 No. 1010.1 conformity (115 volt models)
- CE Conformity marking (230 volt models)

#### For Product Assistance



**International Customer Service** 



Eric Liu

Regional Sales Manager, Asia

# **SCAN THE QR CODE**



• Labconco Corporation 8811 Prospect Avenue Kansas City, MO 64132-2696

• Phone: 816-333-8811

### **FAQ**

### Q: What is the purpose of the Guardian Airflow Monitor?

A: The Guardian Airflow Monitor continuously monitors downflow airflow and indicates when the airflow is within or outside the set point range, ensuring a controlled environment.

### Q: Can the UV lamp be turned off?

A: Yes, the UV lamp has a 3-way switch that allows you to choose between fluorescent, UV, or turning it off completely.

#### Q: Are there different sizes available for the Purifier Vertical Clean Bench?

A: Yes, the Purifier Vertical Clean Bench is available in 2', 3', 4', and 6' widths to suit different laboratory needs.

# **Documents / Resources**



# LABCONCO 3970423 Purifier Vertical Clean Bench with Airflow Monitor [pdf] Owner's Man

3970423 Purifier Vertical Clean Bench with Airflow Monitor, 3970423, Purifier Vertical Clean Bench with Airflow Monitor, Vertical Clean Bench with Airflow Monitor, Airflow Monitor, Monitor

### References

- TCPDF
- 4' Purifier Vertical Clean Bench with Airflow Monitor Labconco
- International Customer Service Labconco
- L Eric Liu Labconco
- User Manual

#### Manuals+, Privacy Policy

This website is an independent publication and is neither affiliated with nor endorsed by any of the trademark owners. The "Bluetooth®" word mark and logos are registered trademarks owned by Bluetooth SIG, Inc. The "Wi-Fi®" word mark and logos are registered trademarks owned by the Wi-Fi Alliance. Any use of these marks on this website does not imply any affiliation with or endorsement.