



Home » Logic Mark » Logic Mark Alert 911 Plus Provides Reliable Owner's Manual 📆



Logic Mark Alert 911 Plus Provides Reliable Owner's Manual



Contents [hide]

- 1 Congratulations!
- 2 Important Product Information
- 3 Device Features
- 4 INCLUDED ACCESSORIES
- 5 Setting Up a New Device
- 6 Tips and Reminders
- 7 Customer Support
- **8 WARRANTY INFORMATION**
- 9 Documents / Resources
 - 9.1 References

Congratulations!

You're a few steps away from staying connected to emergency services

Important Product Information

§ 15.19 Labelling requirements. This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) This device may not cause harmful interference, and(2) this device must accept any interference received, including interference that may cause undesired operation.

§ 15.21 Information to user. Any Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

§ 15.105 Information to the user. Note: This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures: -Reorient or relocate the receiving antenna. -Increase the separation between the equipment and receiver. -Connect the equipment into an outlet on a circuit different from that to which the receiver is connected. -Consult the dealer or an experienced radio/TV technician for help. This equipment complies with FCC radiation exposure limits set forth for an uncontrolled environment. This equipment should be installed and operated with minimum distance 20cm between the radiator & your body.

IC Warning: This device complies with Industry Canada license-exempt RSS standard(s).

Operation is subject to the following two conditions: (1) this device may not cause interference, (2) this device must accept any interference, including interference that may cause undesired operation of the device.

Device Features



Two-way voice communication

Call emergency services at the push of a button

Powered by 4G LTE Technology



STATUS LIGHT CHART

SOLID RED: Fully Charged

SOLID BLUE: Placing a Call

FLASHING RED: Charging

FLASHING BLUE: Fall Detection is ON

INCLUDED ACCESSORIES

Charging Base

Lanyard

• Belt Clip

Setting Up a New Device

DEVICE ACTIVATION

If your device was purchased anywhere other than Logic Mark. com, please call our dedicated customer support team at 1 800-519-2419. They will assist you in activating your device and purchasing any necessary subscriptions.

SETTING UP A DEVICE

STEP 1: Plug the charging base into a power outlet.

STEP 2: Place your device into the charging base with the charging contacts touching

inside the base. Guardian Alert 911 Plus will say "Your Device is Powering On". It will

then say "Your Device is Ready".

STEP 3: Charge Guardian Alert 911 Plus fully.

STEP 4: Perform a test call by pressing and holding the emergency button until it

vibrates. After a short pause you will speak with an emergency services representative.

CONFIRMING BATTERY STATUS

STEP 1: Place your device in the charging base until it announces "Charging".

STEP 2: Remove it from the base and press the battery test button.

STEP 3: Your device will announce "Battery Low" or "Battery Okay".

Tips and Reminders

 Charge your Guardian Alert 911 Plus device nightly to ensure it always charged in an emergency.

 It is essential to test your Guardian Alert 911 Plus every month ensure that it is functioning properly.

 Your Guardian Alert 911 Plus has fall detection included as a feature. No fall detection device can detect 100% of falls. You should always push the emergency button to call for help whenever possible.

Customer Support

CONTACT US VIA PHONE, EMAIL OR ONLINE

To change device settings Customer Support must be contacted.

Call: 1 800-519-2419

Email: info@logicmark.com

Go to: https://www.logicmark.com/customer-support/

WARRANTY INFORMATION

To register your FREE 1-year limited warranty, go to:

https://www.logicmark.com/warranty-manuals/

DISCLAIMERS

Privacy Policy: https://www.logicmark.com/privacy-policy/

Terms and Conditions: https://www.logicmark.com/terms/

Documents / Resources



Logic Mark Alert 911 Plus Provides Reliable [pdf] Owner's Manual PT102A, UDV-201606, Alert 911 Plus Provides Reliable, Alert 911, Plus Provides Reliable, Provides Reliable, Reliable

References

- User Manual
- **■** Logic

Mark

◆ Alert 911, Alert 911 Plus Provides Reliable, Logic Mark, Plus Provides Reliable, Provides Reliable, PT102A, RELIABLE, UDV-201606

Leave a comment

Your email address will not be published. Required fields are marked *
Comment *
Name
Email
Walasita
Website

☐ Save my name, email, and website in this browser for the next time I comment.

Post Comment

Search:

e.g. whirlpool wrf535swhz Search

Manuals+ | Upload | Deep Search | Privacy Policy | @manuals.plus | YouTube

This website is an independent publication and is neither affiliated with nor endorsed by any of the trademark owners. The "Bluetooth®" word mark and logos are registered trademarks owned by Bluetooth SIG, Inc. The "Wi-Fi®" word mark and logos are registered trademarks owned by the Wi-Fi Alliance. Any use of these marks on this website does not imply any affiliation with or endorsement.