


LoliTrack 5 Video Tracking and Behavior Analysis Software



LoliTrack 5 Video Tracking and Behavior Analysis Software User Guide


[Home](#) » [LoliTrack](#) » LoliTrack 5 Video Tracking and Behavior Analysis Software User Guide 

Contents

- [1 LoliTrack 5 Video Tracking and Behavior Analysis Software](#)
- [2 Product Information](#)
- [3 Product Usage Instructions](#)
- [4 FAQs](#)
- [5 FEATURES](#)
- [6 PARAMETERS](#)
- [7 Documents / Resources](#)
 - [7.1 References](#)

LoliTrack 5

LoliTrack 5 Video Tracking and Behavior Analysis Software



Product Information

- **Specifications:**

- **Product Name:** LoliTrack 5
- **Type:** Video tracking and behaviour analysis software
- **Features:** Multiple Individuals Tracking, Kinematics Analysis, Heat Mapping, Heart Rate Monitoring
- **Demo:** Free online demo available
- **Website:** www.loligosystems.com
- **Email:** mail@loligosystems.com

Product Usage Instructions

- **Installation:**

- Download the LoliTrack 5 software from the official website. Follow the installation instructions to set up the software on your computer.

- **Setup:**

- Launch the LoliTrack 5 software and configure the settings according to your tracking and analysis requirements.

- **Video Tracking:**

- Import the video files you want to analyze into the software. Use the video tracking feature to track multiple individuals simultaneously.

- **Behavior Analysis:**

- Utilize the behaviour analysis tools provided in LoliTrack 5 to analyze the movements and interactions of the tracked individuals.

- **Kinematics Analysis:**

- Access the kinematics analysis feature to study the motion patterns and dynamics of the subjects in the videos.

- **Heat Mapping:**

- Generate heat maps based on the movement data collected during tracking to visualize hotspots and activity patterns.

- **Heart Rate Monitoring:**

- If applicable, integrate heart rate monitoring devices with the software to analyze heart rate data alongside behaviour and movement patterns.

FAQs

- **Q: Is technical support available for LoliTrack 5?**

- A: Yes, technical support is provided by contacting support at mail@lologosystems.com

- **Q: Can I track more than one individual in a single video?**

- A: Yes, LoliTrack 5 supports tracking multiple individuals simultaneously.

- **Q: Are there tutorials available for using the software?**

- A: Yes, you can find tutorials and guides on the LoliTrack website to help you get started with the software.

FEATURES

- No subscription fees. Free updates. Install on multiple computers.
- Track single and multiple free-ranging animals
- Contrast or background subtraction tracking (no markers needed)
- Any size or shape of the arena, e.g., maze, tank, pen, chamber, etc.
- Works with most video cameras and file formats
- Track tiny organisms in 24-, 48-, and 96-well microtiter plates
- Built-in video recorder and editor
- Setup multiple zones of interest
- Analyze entire videos or just user-defined intervals
- Data output for individuals, treatments, and zones
- Track specific body points for kinematic studies or biomechanics
- In-depth analysis tools with customizable graphs
- Export data (to Excel), images, and videos
- Frequency tool for heart rate, ventilation, and tail beat frequency
- 3D tracking with exportable coordinates and trajectories
- Batch analyzes multiple videos to save time
- Correct tracking errors with ease
- Fast performance on regular office PCs
- Genuine Windows 11 64-bit software

APPLICATIONS

- Microplate analysis (e.g., ZF embryos)
- Behavioural phenotyping (focal animals)
- Rodent assays (e.g., maze, open field)


- Group behaviour (e.g., insects, fish)
- Locomotory behaviour (e.g., kinematic)
- Social and sexual behaviours

PARAMETERS

- Active/Inactive
- Acceleration
- Bend
- Direction of movement
- Direction of Orientation
- Distance moved
- Distance to the centre of the zone
- Frequencies (tailbeat, heart rate etc.)
- Inter-Individual Distance
- Nearest Neighbor Distance
- Speed
- Turning rate movement
- Turning rate orientation
- Velocity
- X/Y position (raw and calibrated)



Documents / Resources



[LoliTrack LoliTrack 5 Video Tracking and Behavior Analysis Software](#) [pdf] User Guide
LoliTrack 5 Video Tracking and Behavior Analysis Software, LoliTrack 5, Video Tracking and Behavior Analysis Software, Tracking and Behavior Analysis Software, Behavior Analysis Software , Analysis Software

References

- [loligosystems.com](https://www.loligosystems.com)
- [User Manual](#)

Manuals+, Privacy Policy

This website is an independent publication and is neither affiliated with nor endorsed by any of the trademark owners. The "Bluetooth®" word mark and logos are registered trademarks owned by Bluetooth SIG, Inc. The "Wi-Fi®" word mark and logos are registered trademarks owned by the Wi-Fi Alliance. Any use of these marks on this website does not imply any affiliation with or endorsement.