



# LUMEX 7-2000 Convoluted Foam Overlay Unpacking Instructions

[Home](#) » [LUMEX](#) » LUMEX 7-2000 Convoluted Foam Overlay Unpacking Instructions 

## Contents

- [1 LUMEX 7-2000 Convoluted Foam Overlay Unpacking](#)
- [2 UNPACKING INSTRUCTIONS](#)
- [3 Documents / Resources](#)
- [4 Related Posts](#)



## LUMEX 7-2000 Convoluted Foam Overlay Unpacking



## UNPACKING INSTRUCTIONS


## FOAM MATTRESS UNPACKING INSTRUCTIONS

1. Check for any obvious damage to the packaging or its contents. If damage is evident, DO NOT USE. Contact carrier / distributor for further instruction.
2. Open overlay, unfold it gently, and spread it out on a flat surface.
3. Once unwrapped, the overlay's compressed foam may take 72 hours to fully recover. Environmental factors, such as cooler temperatures, may cause expansion to take longer. Mattress will expand to  $\pm 1"$  dimensions stated below.
4. Overlay dimensions will be smaller than mattress-top dimensions to allow a better fit for sheets and bedding.

- Twin.....34" x 72"
- Full.....52" x 72"
- Queen.....58" x 78"
- King.....76" x 78"

© 2011, GF Health Products, Inc. All Rights Reserved. Lumex is a trademark of GF Health Products, Inc., One Graham-Field Way, Atlanta GA 30340-3140 GF Health Products, Inc. is an ISO 13485 Certified Company.

## Documents / Resources

	<p><b>LUMEX 7-2000 Convolted Foam Overlay Unpacking</b> [pdf] Instructions</p> <p>7-2000 Convolted Foam Overlay Unpacking, 7-2000, Convolted Foam Overlay Unpacking, Foam Overlay Unpacking, Overlay Unpacking, Unpacking</p>
--	---