

REPLICATING NATURE.

Product Portfolio Fall 2018



Signia.

Taking the future of a sound legacy to the next level.


We invent the future of better hearing. Beginning with Siemens, this has been our consistent mission for 140 years and counting – the greatest legacy in the industry. As hearing aid pioneers, we have introduced countless world’s firsts and now we are continuing this success story by expanding the portfolio of our game-changing platform Signia Nx and enhancing the world of hearing aids.

Hearing solutions for everyone.

With Signia hearing aids, you can offer your customers a wide selection from many different models and performance levels for every requirement and all degrees of hearing loss. To facilitate orientation and help in the selection of the right hearing solution, our hearing aids are placed in three categories according to their construction models. This lets you find out quickly which model best suits a customer’s wearing preferences.

Contents

Hearing solutions for everyone	2	Run	30
Hearing aids overview	4	Fun	31
Signia Nx feature overview	5	Fast	32
primax feature overview	6	SecureEar	33
Essential segment feature overview	7	Apps feature overview	34
Lotus segment feature overview	8	App overview	35
Signia Styletto	9	Apps, accessories and hearing aid compatibility	36
Pure Charge&Go Nx	12	Accessory overview	37
Pure 312 Nx	13	Color overview	38
Pure 13 Nx	14	Hearing care professional support	40
Pure 10 Nx	15	Overview SLIM-RIC/RICs	42
Motion Charge&Go Nx	16	Overview BTEs	44
Motion 13 Nx	17	Overview Click ITEs/ITEs	48
Silk Nx	18		
Insio Nx	19		
Ace primax	20		
Pure primax	21		
Cellion primax	22		
Carat primax	23		
Motion primax	24		
Silk primax	25		
Insio primax	26		
Intuis 3 BTEs	27		
Intuis 3 Click ITEs and ITEs	28		
Prompt	29		



The degrees of hearing loss indicated above serve as a rough guide. For precise information, please refer to the individual product pages.

Hearing aids overview.

We offer a broad range of hearing aids in many different models and performance levels to help hearing care professionals meet their clients' needs, preferences and demands. Take advantage of this diversity to enhance your customers' satisfaction.

Models		Performance levels	RIC		BTE		ITE				Fitting
			Slim	Regular	LifeTube/ ThinTube	Earhook	IIC	CIC	ITC	ITE	Connexx
Styletto	Styletto New!	7 5 3 Powered by Signia Nx	●								8
Premium	Pure 10 Nx	7Nx 5Nx 3Nx		●							8
	Pure 312 Nx	7Nx 5Nx 3Nx		●							8
	Pure 13 Nx	7Nx 5Nx 3Nx		●							8
	Pure Charge&Go Nx	7Nx 5Nx 3Nx		●							8
	Motion 13 Nx	7Nx 5Nx 3Nx			●	●					8
	Motion 13P Nx	7Nx 5Nx 3Nx			●	●					8
	Motion Charge&Go Nx New!	7Nx 5Nx 3Nx			●	●					8
	Silk Nx	7Nx 5Nx 3Nx						●*			8
	Insio Nx New!	7Nx 5Nx 3Nx					●	●	●	●	8
	Ace primax	7px 5px 3px 2px 1px		●							8
	Pure primax	7px 5px 3px 2px 1px		●							8
	Cellion primax	7px 5px 3px 2px		●							8
	Carat primax	7px 5px 3px		●							8
	Motion primax	7px 5px 3px 2px 1px			●	●					8
Essential	Intuis 3				●	●		●	●	●	8
Lotus	Prompt				●	●		●	●	●	8
	Run				●	●		●*	●*		MobileFitting App
	Fun					●					Manual Fit
	Fast					●					Manual Fit

* Only available as Click ITE.

Signia Nx feature overview.

Our game-changing platform Signia Nx delivers the most natural own voice and the best speech understanding in noise, while combining the most advanced connectivity systems.

Features	Information	Performance levels		
Audiology		7Nx	5Nx	3Nx
OVP™ (Own Voice Processing)	Own Voice Processing (OVP™) utilizes real-time recognition of the wearer's voice to deliver a natural own-voice impression. Available for RIC/BTE devices with Ultra HD e2e.	●●●●●	●●●●●	●●●●●
3D Classifier	Unlike conventional classification systems which rely on acoustic data, 3D Classifier analyzes multi-dimensional information from the wearer's environment and communication behavior. Information on the wearer's voice activity, motion state, and acoustic data are used to tailor the hearing aid's response to the wearer's needs. Motion detection available for Bluetooth-equipped devices.	●●●●●	●●●●●	●●●●
Sound Clarity	Sound Clarity provides lifelike sound which covers the full spectrum of wearing situations. Exceptional clarity is achieved with features such as: HD Spatial, extended dynamic range, extended bandwidth, EchoShield, HD Music, eWindScreen™, speech and noise management and feedback cancellation.	●●●●●	●●●●●	●●●●
Speech Quality	Speech Quality delivers improved understanding of other people's voices in noise by using an exclusive technology called Ultra HD e2e. Features which utilize Ultra HD e2e include narrow directionality, Spatial SpeechFocus™ and TwinPhone.	●●●●●	●●●●●	●●●●
Connectivity				
Direct Streaming	Signia Nx hearing aids can be connected directly to Apple devices for phone calls and audio streaming – without additional accessories. With the StreamLine Mic accessory, direct streaming is possible also from Android devices and other Bluetooth mobile devices. With the StreamLine TV accessory, high-quality stereo sound can be streamed directly to Nx hearing aids. Available for Bluetooth-equipped devices.	●●●●●	●●●●●	●●●●●
You Adjust	Technically inclined wearers often appreciate having extra control over their hearing aids. Signia Nx hearing aids offer these wearers a comprehensive range of adjustment options via the myControl App™ or touchControl App™.	●●●●●	●●●●●	●●●●●
Remote hearing care				
TeleCare 3.0	TeleCare provides a structured communication approach for hearing care professionals and their clients, which is especially important during the initial stages of a fitting. In addition, TeleCare provides full live remote tuning to assist with rapid and effective troubleshooting. Available for Bluetooth-equipped devices.	●●●●●	●●●●●	●●●●●

●●●●● = highest feature performance

primax feature overview.

Redefining the ease of listening, the outstanding primax™ star platform offers effortless hearing in every situation throughout the day.*

Features		Information	Performance levels				
Audiology			7px	5px	3px	2px	1px
Signal processing and Gain/MPO	Channels/Handles in Connexx™		48/20	32/16	24/12	16/8	16/8
Hearing programs	Maximum available		6	6	6	4	4
SpeechMaster			●	●	●	●	●
HD Music	Number of preset programs		3	1	1	-	-
TwinPhone ^{A,E}	Req. bilateral fitting & e2e wireless™ 3.0		●	●	●	-	-
EchoShield			●	-	-	-	-
Wireless CROS/BiCROS ^{A,E}	Req. bilateral fitting & e2e wireless 3.0		●	●	●	-	-
Directionality	Resolution in channels		48	32	24	16	16
Binaural OneMic directionality ^B	Req. bilateral fitting & e2e wireless 3.0		●●●●●	-	-	-	-
Narrow directionality ^{A,C}	Req. bilateral fitting & e2e wireless 3.0		●●●●●	●●●●	●●●	-	-
Directional microphones ^C			●●●●●	●●●●	●●●	●●	●
Spatial SpeechFocus ^{A,C}	Req. bilateral fitting & e2e wireless 3.0		●●●●●	-	-	-	-
SpeechFocus ^C			●●●●	●●●●	-	-	-
TruEar™	Req. BTE or RIC hearing aids		●●●●●	●●●●	●●●	-	-
Frequency compression			●	●	●	●	●
Extended bandwidth	Up to 12 kHz		●	-	-	-	-
Feedback cancellation			●	●	●	●	●
eWindScreen binaural ^{A,C,D}	Req. bilateral fitting & e2e wireless 3.0		●	●	-	-	-
eWindScreen	Adjustment steps		3	3	1	1	-
Noise reduction	Resolution in channels/Combined adjustment steps		48/5	32/5	24/3	16/3	16/1
Speech and noise management	Individual adjustment steps		7	5	3	3	1
SoundSmoothing™	Individual adjustment steps		3	3	1	1	-
Directional speech enhancement ^C	Individual adjustment steps		3	1	-	-	-
Adaptive streaming volume ^A	Streaming only, req. easyTek & e2e wireless 3.0		●	-	-	-	-
SoundBrilliance™ ^A	Streaming only, req. easyTek & e2e wireless 3.0		●	●	-	-	-
Sound equalizer	Number of classes		6	3	-	-	-
User controls, remote controls, streaming							
touchControl App ^F	iOS and Android devices only		●	●	●	●	●
easyTek/easyPocket ^A	Req. e2e wireless 3.0		●	●	●	●	●
Spatial Configurator ^{A,C}	Req. bilateral fitting & e2e wireless 3.0		●	●	-	-	-
Span	Req. easyTek & easyTek App, touchControl App or rocker switch		●	●	-	-	-
Direction	Req. easyTek & easyTek App		●	●	-	-	-
SoundBalance™			●	●	●	●	●
Fitting							
InSituGram™			●	●	●	●	●
Learning/Data logging	Number of classes		6/6	3/6	1/6	-/6	-/6
Acclimatization manager			●	●	●	●	●
Tinnitus							
Tinnitus therapy							
Static therapy signal (5 presets)	Handles in Connexx		20	16	12	4 ^F	-
Ocean wave therapy signal (4 presets)			●	●	●	-	-
Notch therapy			●	●	●	-	-

A) Not for Ace primax. B) Only for hearing aids with one microphone each. C) Only for hearing aids with two microphones each. D) Not available in the universal program on 5px. E) Not for Motion S primax. F) No presets.

●●●●● = highest feature performance

* See legal note on final page in this document.

Essential segment feature overview.

Proven technology, outstanding comfort: In a broad range of hearing aid models, the Essential line offers a reliable hearing experience, a high level of wearing comfort and easy handling for almost every customer.

Features

Performance levels

Essential

Better hearing

Intuis 3

Frequency channels/Compression controls	12/6
Number of programs	4
Directional microphones ^A	automatic ^B
Frequency compression	●

Sound comfort

Feedback cancellation	●
Speech and noise management	on/off
eWindScreen	-
e2e wireless 3.0	-
e2e wireless 2.0	-
SoundSmoothing	-

Fitting/Individuality

Data logging	●
Acclimatization manager	-
Tinnitus therapy	-

Audio streamer/Hearing aid control

touchControl App	●
easyTek and easyTek App	-
easyPocket	-
ConnexxAir	-

A) Not for CIC. B) Only BTE.

Lotus segment feature overview.

Available in a wide range of fittings and suitable for all levels of hearing loss from mild to profound, the Lotus range offers a good hearing experience for every customer.

Features	Prompt			Click ITEs	CIC	ITC	ITE
	S	P	SP				
Matrix	LifeTube 45/124 dB Earhook 55/124 dB	70/134 dB	80/140 dB	Click CIC 50/113 dB Click ITC 55/116 dB	65/124 dB ¹	65/124 dB ¹	65/124 dB ¹
Channels/Controls	8/4	8/4	8/4	8/4	8/4	8/4	8/4
Programs	4	4	4	Click CIC: 1 Click ITC: 4	4	4	4
Directional microphones	●	●	●	-	-	-	-
T-coil	●	●	●	-	-	○	○
Audio input	-	●	●	-	-	-	-
Noise reduction	situation based	situation based	situation based	situation based	situation based	situation based	situation based
SoundBalance™	●	●	●	●	●	●	●
Feedback cancellation	●	●	●	●	●	●	●

Features	Run			Click ITEs	Fun		Fast
	S	P	SP		P	SP	P
Matrix	LifeTube 45/124 dB Earhook 55/124 dB	70/132 dB	79/140 dB	Click CIC 50/113 dB Click ITC 55/116 dB	70/132 dB	79/140 dB	62/130 dB
Channels	8	8	8	8	6	6	4
Programs	6 preconfigured sound profiles	6 preconfigured sound profiles	6 preconfigured sound profiles	6 preconfigured sound profiles	4 preconfigured sound profiles	4 preconfigured sound profiles	3 preconfigured sound profiles
Directional microphones	●	●	-	-	-	-	-
T-coil	-	-	-	-	-	-	-
Audio input	-	-	-	-	-	-	-
Noise reduction	situation based	situation based	situation based	situation based	always on	always on	always on
SoundBalance	●	●	●	●	-	-	-
Feedback cancellation	●	●	●	●	●	●	●

(1) Max. values



Introducing Signia Styletto.

Too many people who need hearing aids don't want to wear them simply because the look and feel does not match their image of themselves. At last you can attract new and existing clients of all ages and convert them into satisfied wearers.

Styletto is more than a hearing aid. It is hearwear with a revolutionary design: the world's first SLIM-RIC.

It's not just a hearing aid.

It's embracing life with stylish high-tech.

By thinking differently, we have produced a hearing solution that transforms the image of hearing aids and is as easy to sell as traditional RICs. Styletto combines all the audiological advantages of our acclaimed Signia Nx platform with a stylishly designed new form factor.

It provides the latest advancements in terms of energy-efficient chip design, unique lithium-ion batteries, and the most natural sound. All in a refined slim-line styling that is preferred by more than 8 out of 10 people compared to conventional hearing aids.*

Experience the Styletto hearwear revolution.

* Source: Signia consumer study 2018.

Performance levels

7

5

3

sDemo

Powered by Signia Nx

New!

OWN
VOICE
PROCESSING
OVP



Design
The unique, beautifully designed Styletto delivers strong and proven consumer preference.

A new charging paradigm
Wearers enjoy full-day wearing and **charging-on-the-go** via our ultra-slim portable charging case.



DARK GRANITE/SILVER



SNOW WHITE/ROSE GOLD





COSMIC BLUE/ROSE GOLD


Colors
Wearers express their individual style with three contemporary color combinations.

Key features

- Stylishly designed new form factor
- Most natural own voice and best speech understanding in noise thanks to Nx chip
- Portable lithium-ion rechargeability

 Compatible with touchControl App.

 TeleCare 2.0 ready.

 Compatible with miniPocket.

SLIM-RIC

Receiver type miniReceiver 2.0 M

Matrix (2 ccm) 60/119 dB
Battery type Lithium-ion power cell
Power consumption 1.4 mA



touchControl App and miniPocket
To enhance its convenience, Styletto is remotely controllable via the touchControl App and the miniPocket.

The most natural own voice with convenient wireless rechargeability.



Performance levels

7Nx

5Nx

3Nx

sDemo

Pure Charge&Go Nx.

The ultimate RIC, Pure™ Charge&Go Nx combines all the advantages of Signia Nx with convenient wireless rechargeability. It provides the most natural own voice and unrivaled speech intelligibility, top Bluetooth connectivity, full live remote support with Signia TeleCare 3.0, and the most high-tech wireless charging in a remarkably small design. Wearers no longer worry that they might struggle to change disposable batteries. Simply Charge&Go!

Key features

- Most natural own voice and best speech understanding in noise thanks to Nx chip
- Long-lasting high quality audio streaming
- Ultra HD e2e and Bluetooth connectivity



Compatible with myControl App.



TeleCare 3.0 ready.



Compatible with CROS device for people with unaidable hearing loss on one side.



StreamLine TV

With StreamLine TV, the TV signal is directly streamed into Pure Charge&Go Nx hearing aids.



StreamLine Mic

The StreamLine Mic enables hands-free audio streaming (phone calls, music, etc.) from any device with Bluetooth connectivity.

RIC

Model				
Receiver type	miniReceiver 2.0 S	miniReceiver 2.0 M	miniReceiver 2.0 P	miniReceiver 2.0 HP
Performance				
Matrix (2 ccm)	45/108 dB	60/119 dB	70/124 dB	75/130 dB
Battery				
Type	Lithium-ion power cell	Lithium-ion power cell	Lithium-ion power cell	Lithium-ion power cell
Power consumption	1.3 mA	1.5 mA	1.4 mA	1.4 mA

The smallest fully featured hearing aid with direct streaming.

Key features

- Most natural own voice and best speech understanding in noise thanks to Nx chip
- Small, discreet design
- Ultra HD e2e and Bluetooth connectivity



Compatible with myControl App.



TeleCare 3.0 ready.



Compatible with CROS device for people with unaidable hearing loss on one side.



Performance levels

7Nx

5Nx

3Nx

sDemo

Pure 312 Nx.

Equipped with all benefits of the Signia Nx platform, Pure 312 Nx keeps the wearer's own voice and remaining soundscape in perfect balance. This way, it offers the most natural own voice and best speech understanding in noise in a small, elegant housing. Despite its discreet design, Pure 312 Nx provides supreme connectivity via Bluetooth-based direct streaming and the myControl App*. Plus, TeleCare 3.0 lets your customers enjoy full live remote support and fast, smooth trials.

* For Android and iOS.







StreamLine TV
With StreamLine TV, the TV signal is directly streamed into Pure 312 Nx hearing aids.



StreamLine Mic
The StreamLine Mic enables hands-free audio streaming (phone calls, music, etc.) from any device with Bluetooth connectivity.

RIC

Model	RIC			
Receiver type	miniReceiver 2.0 S	miniReceiver 2.0 M	miniReceiver 2.0 P	miniReceiver 2.0 HP
				
Performance				
Matrix (2 ccm)	45/108 dB	60/119dB	70/124 dB	75/130 dB
Battery				
Size	312	312	312	312
Power consumption	1.3 mA	1.5 mA	1.4 mA	1.4 mA

Superior connectivity combined with the longest streaming time.



Performance levels

7Nx

5Nx

3Nx

sDemo

Pure 13 Nx.

By separately processing the wearer's own voice and integrating it smoothly into the soundscape, Signia's Pure 13 Nx provides the most natural own voice and best speech understanding in noise. Thanks to Bluetooth, it comes with superior connectivity with direct streaming and the myControl App*. Users enjoy the longest wearing time while streaming. Full live remote support via TeleCare 3.0 quickly familiarizes them with their Pure 13 Nx hearing aids.

* For Android and iOS.

Key features

- Most natural own voice and best speech understanding in noise thanks to Nx chip
- Longest wearing time while streaming
- Ultra HD e2e and Bluetooth connectivity



Compatible with myControl App.



TeleCare 3.0 ready.




Standard battery door can be exchanged for one with an integrated T-Coil.



Compatible with CROS device for people with unaidable hearing loss on one side.



StreamLine TV
With StreamLine TV, the TV signal is directly streamed into Pure 13 Nx hearing aids.



StreamLine Mic
The StreamLine Mic enables hands-free audio streaming (phone calls, music, etc.) from any device with Bluetooth connectivity.

RIC

Model				
Receiver type	miniReceiver 2.0 S	miniReceiver 2.0 M	miniReceiver 2.0 P	miniReceiver 2.0 HP
Performance				
Matrix (2 ccm)	45/108 dB	60/119 dB	70/124 dB	75/130 dB
Battery				
Size	13	13	13	13
Power consumption	1.3 mA	1.5 mA	1.4 mA	1.4 mA

Simplicity with a naturally clear sound.

Key features

- Ultra-small design
- Outstanding speech understanding in noise
- Remote control option



Compatible with touchControl App.



TeleCare 2.0 ready.



touchControl App and miniPocket

To enhance its discretion, Pure 10 Nx is remotely controllable via the touchControl App and the miniPocket.

Performance levels

7Nx

5Nx


3Nx

sDemo

Pure 10 Nx.

The ultra-small Pure 10 Nx expertly combines a highly discreet design with a more natural sound. Its sleek appearance is complemented by a wealth of audiological innovations that produce a clear sound and excellent speech understanding in noise. The non-wireless Pure 10 Nx is the ideal solution for those who wish their hearing aids to be simple and discreet.

RIC

Model				
Receiver type	miniReceiver 2.0 S	miniReceiver 2.0 M	miniReceiver 2.0 P	miniReceiver 2.0 HP
				
Performance				
Matrix (2 ccm)	45/108 dB	55/113 dB	60/118 dB	65/123 dB
Battery				
Size	10	10	10	10
Power consumption	1.2 mA	1.3 mA	1.2 mA	1.2 mA

The ultimate rechargeable BTE for the most natural own voice.

OWN VOICE PROCESSING
OVP

NEW



Key features

- Ultra HD e2e and Bluetooth connectivity
- Rechargeable Lithium-Ion power cell for durable performance even when streaming
- Robust housing



Compatible with myControl App.



TeleCare 3.0 ready.



Compatible with CROS device for people with unaidable hearing loss on one side.



StreamLine TV
With StreamLine TV, the TV signal is directly streamed into Motion Charge&Go Nx hearing aids.




StreamLine Mic
The StreamLine Mic enables hands-free audio streaming (phone calls, music, etc.) from any device with Bluetooth connectivity.

Performance levels

- 7Nx
- 5Nx
- 3Nx
- sDemo

Motion Charge&Go Nx.

The incredibly versatile Motion™ Charge&Go Nx combines all the advantages of Signia Nx with Li-ion powered wireless rechargeability and seamless Bluetooth connectivity. It is intuitive to handle because wireless charging is integrated into its dependable design. Ideal for the most varied fitting requirements, Motion Charge&Go Nx provides the most natural own voice for highest acceptance, convenient handling, and effortless enjoyment of TV, music and phone calls.

BTE		
Model		
		
Performance		
Matrix (2 ccm)	ThinTube 56/126 dB	Earhook 63/134 dB
Battery		
Type	Lithium-ion power cell	Lithium-ion power cell
Power consumption	1.3 mA	1.5 mA

Direct streaming with the most versatile fitting options.



Performance levels

7Nx

5Nx

3Nx

sDemo

Motion 13 Nx.

Available in a small S model and a sleek P model, the incredibly versatile Motion 13 Nx BTE is ideal for wearers with individual fitting requirements and moderate to severe or even profound hearing loss. It harmonizes the wearer's voice and surrounding sounds, providing the most natural own voice and best speech understanding in noise. Motion 13 Nx and Motion 13P Nx provide full Bluetooth connectivity and Signia TeleCare 3.0 with an optional modular T-Coil integrated in the battery door.

Key features

- Most natural own voice and best speech understanding in noise thanks to Nx chip
- For individual handling and fitting requirements
- Ultra HD e2e and Bluetooth connectivity



Compatible with myControl App.



TeleCare 3.0 ready.




Standard battery door can be exchanged for one with an integrated T-Coil.



Compatible with CROS device for people with unaidable hearing loss on one side.



StreamLine TV
With StreamLine TV, the TV signal is directly streamed into Motion 13 Nx hearing aids.



StreamLine Mic
The StreamLine Mic enables hands-free audio streaming (phone calls, music, etc.) from any device with Bluetooth connectivity.

BTE

Model		
	<p>Motion 13 Nx</p>	<p>Motion 13P Nx</p>
Performance	<p>ThinTube 56/126 dB Earhook 63/134 dB</p>	<p>ThinTube 65/129 dB Earhook 77/135 dB</p>
Battery		
Size	13	13
Power consumption	1.3 - 1.5 mA	1.3 - 1.5 mA

A practically invisible natural experience. Simply fit and go.

Key features

- Best speech understanding in noise and a most natural sound thanks to Nx chip
- Click Sleeves for a quick and easy fit
- Even smaller than its predecessor
- Invisible and discreet
- Binaural OneMic directionality
- Ultra HD e2e for superior speech understanding and a most natural sound




Compatible with touchControl App.



TeleCare 2.0 ready.



Compatible with CROS device for people with unaidable hearing loss on one side.



One click and you're done.

Available in four sizes, Click Sleeves simply click on to Silk Nx hearing aids and sit securely in the ear.

Performance levels

7Nx

5Nx

3Nx


sDemo

Silk Nx.

The Silk™ Nx delivers a natural hearing experience in its discreet size, instant fit, and exceptional sound quality. The only ready-to-wear CIC with Ultra HD e2e for superior wireless binaural processing is also 20% smaller than Silk primax™ for ultimate discretion. Thanks to its flexible silicone Click Sleeves, Silk Nx offers a precise fit without waiting times.

And for wearers with unaidable hearing loss in one ear, CROS Silk Nx provides the most discreet solution – the world's first CROS CIC hearing aid.

Click ITE

Model	
	Click CIC
	
Performance	
Matrix (2 ccm)	50/113 dB
Battery	
Size	10
Power consumption	1.3 mA

The most discreet natural hearing experience with full connectivity.

NEW


Performance levels

7Nx

5Nx

3Nx

Insio Nx.

Signia's new Insio™ Nx ITC and ITE models provide top Bluetooth connectivity with best-in-class sound quality. Tailor-made for the most discerning wearers, Insio Nx delivers a more natural sound with maximum discretion. Equipped with Ultra HD e2e for superior wireless binaural processing including as a small CIC with a single microphone, Insio offers crystal clear sound and unrivalled speech understanding even in very noisy environments.

Key features

- Best speech understanding in noise and a most natural sound thanks to Nx chip
- Small size up to high power matrix
- Full wireless and remote control options for all models
- Ultra HD e2e for superior speech understanding and a most natural sound



Compatible with myControl App.*



TeleCare 3.0 ready.**



Compatible with CROS device for people with unaidable hearing loss on one side.***







StreamLine TV
With StreamLine TV, the TV signal is directly streamed into Insio Nx ITC and ITE hearing aids.



StreamLine Mic
The StreamLine Mic enables hands-free audio streaming (phone calls, music, etc.) from any device with Bluetooth connectivity for Insio Nx ITCs and ITEs.

ITE

Model	IIC	CIC	ITC	ITE
				
Performance				
Matrix (2 ccm)	50/113 dB	S 50/113 dB M 55/118 dB P 65/124 dB	S 50/113 dB M 55/118 dB P 65/124 dB	M 55/118 dB P 65/124 dB
Battery				
Size	10	10	312	312
Power consumption	1.3 mA	1.3 mA ¹	1.3 mA ¹	1.3 mA ¹

(1) Max. values

* Insio Nx ITC and ITE models, only. Insio Nx IIC and CIC models are compatible with touchControl App.

** Insio Nx ITC and ITE models, only. Insio Nx IIC and CIC models are compatible with TeleCare 2.0.

*** Insio Nx ITC and ITE models are compatible with CROS Pure 312 Nx. Insio Nx IIC and CIC models are compatible with CROS Silk Nx.

Discreet.
Powerful.
Effortless.

Key features

- Ultra-small design
- HD Binaural Hearing
- Clinically proven to reduce hearing effort*
- Battery size 10



Compatible with touchControl App.



TeleCare 2.0 ready.



With Notch Therapy for tinnitus treatment.

* See legal note on final page of this document.



Performance levels

7px

5px

3px

2px

1px

Ace primax.

Signia's Ace primax™ skillfully combines the three main elements of discretion: tiny size, automatic functionality and remote control options. Equipped with sophisticated primax star technology, it offers outstanding speech understanding and is an ideal solution for those who wish for a high level of discretion in a hearing aid.

sDemo

RIC

Model				
Receiver type	miniReceiver 2.0 S	miniReceiver 2.0 M	miniReceiver 2.0 P	miniReceiver 2.0 HP
				
Performance				
Matrix (2 ccm)	45/108 dB	55/113 dB	60/118 dB	65/123 dB
Battery				
Size	10	10	10	10
Power consumption	0.9 mA	1.0 mA	0.9 mA	0.9 mA

Fully featured. Binaural. Effortless.



Performance levels

7px

5px

3px

2px

1px

Pure primax.

Equipped with the outstanding primax star platform, Signia's Pure primax™ provides outstanding sound comfort and sophisticated features. Its elegant, discreet housing design also allows for a quick and easy exchange of batteries. Pure primax is rechargeable via the eCharger*, remotely controllable via touchControl App and offers wearers full connectivity for audio streaming via easyTek™.

Key features

- Rechargeable
- Wireless connectivity
- Battery size 312
- HD Binaural Hearing
- Clinically proven to reduce hearing effort*



TeleCare 2.0 ready.



Audio streaming via easyTek and easyTek App.



Compatible with CROS device for people with unaidable hearing loss on one side.

* See legal note on final page of this document.



eCharger

The eCharger recharges and dehumidifies Pure primax hearing aids overnight, so your clients don't need to worry about handling batteries.

sDemo

* Available with performance levels 7px, 5px and 3px only.

RIC

Model				
Receiver type	miniReceiver 2.0 S	miniReceiver 2.0 M	miniReceiver 2.0 P	miniReceiver 2.0 HP
				
Performance				
Matrix (2 ccm)	45/108 dB	60/119 dB	70/124 dB	75/130 dB
Battery				
Size	312	312	312	312
Power consumption	0.9 mA	1.0 mA	1.3 mA	1.3 mA

Revolutionary battery-free hearing.



Key features

- Rechargeable lithium-ion power cell
- 24 hours continuous use with unlimited streaming
- Contactless, inductive charging
- HD Binaural Hearing
- Clinically proven to reduce hearing effort*



TeleCare 2.0 ready.



Audio streaming via easyTek and easyTek App.



Equipped with an integrated T-Coil.

* See legal note on final page of this document.



Simply smart.
When placed in the wells, the Cellion primax hearing aids start charging automatically.

Performance levels

7px

5px





3px

2px

Cellion primax.

Signia's groundbreaking Cellion primax™ offers revolutionary lithium-ion technology in a unique, contact-free OneShell design with no battery door. Once fully charged via the inductive charger, the high-capacity power cell makes Cellion primax the only hearing aids that last all day while offering 24-hour continuous streaming.

sDemo

RIC				
Model				
Receiver type	miniReceiver 2.0 S	miniReceiver 2.0 M	miniReceiver 2.0 P	miniReceiver 2.0 HP
				
Performance				
Matrix (2 ccm)	45/108 dB	60/118 dB	70/124 dB	75/130 dB
Battery				
Type	Lithium-ion power cell	Lithium-ion power cell	Lithium-ion power cell	Lithium-ion power cell
Power consumption	0.9 mA	1.0 mA	1.0 mA	1.1 mA

Power meets versatility.



Key features

- Rechargeable
- Battery size 13
- Wireless connectivity
- HD Binaural Hearing
- Clinically proven to reduce hearing effort*



Compatible with touchControl App.



TeleCare 2.0 ready.



Audio streaming via easyTek and easyTek App.



With Notch Therapy for tinnitus treatment.

* See legal note on final page of this document.

Performance levels

7px

5px





3px

Carat primax.

Signia's Carat primax™ was designed to deliver optimum amplification for almost all levels of hearing loss. Its optional rechargeable battery packs longer lasting power, so wearers can enjoy highly convenient handling. Carat primax' remote control and wireless streaming options make these hearing aids a comfortable and discreet solution.

sDemo

RIC

RIC				
Model				
Receiver type	<div>miniReceiver 2.0 S</div> 	<div>miniReceiver 2.0 M</div> 	<div>miniReceiver 2.0 P</div> 	<div>miniReceiver 2.0 HP</div> 
Performance				
Matrix (2 ccm)	45/108 dB	60/119 dB	70/124 dB	75/130 dB
Battery				
Size	13	13	13	13
Power consumption	0.9 mA	1.1 mA	1.0 mA	1.1 mA

Sleek. Convenient. Effortless.



Key features

- Sleek and compact housing
- Easy handling
- HD Binaural Hearing
- Clinically proven to reduce hearing effort*
- Option for Motion SA, P and SP primax models: Direct audio input via audio shoe



Compatible with touchControl App.



TeleCare 2.0 ready.



Audio streaming via easyTek and easyTek App.



With Notch Therapy for tinnitus treatment.



Direct Audio Input

* See legal note on final page of this document.

Performance levels

7px

5px

3px

2px





1px

sDemo

Motion primax.

Motion primax™ from Signia comes in a broad and versatile BTE range. From the small S model to the sleek but powerful P model, from the SA model that can be equipped with an audio shoe to the rechargeable SX to the super power SP model, Motion primax provides a solution for almost every ear and every hearing loss. For a maximum level of convenience, all models of the Motion family are remotely controllable and offer full connectivity for audio streaming.

BTE

Model	Motion S primax	Motion SX/SA primax	Motion P primax	Motion SP primax
				
Performance	LifeTube 45/124 dB Earhook 55/124 dB	ThinTube 53/126 dB Earhook 60/130 dB	ThinTube 65/129 dB Earhook 77/135 dB	Earhook 82/140 dB
Battery				
Size	312	13	13	675
Power consumption	0.8 mA	1.0 mA	1.1 - 1.3 mA	2.1 mA

Click.
Fit.
Go.

Key features

- Click Sleeves for a swift fit
- 3 sleeve sizes to fit every ear
- Invisible and discreet
- Binaural OneMic directionality
- Clinically proven to reduce hearing effort*



TeleCare 2.0 ready.



Audio streaming via easyTek and easyTek App.



With Notch Therapy for tinnitus treatment.

* See legal note on final page of this document.



Performance levels

7px

5px


3px

2px

1px

Silk primax.


Delivering the full convenience of Click CICs, the nearly invisible Signia Silk primax lets you provide the smoothest, fastest ITE fittings. Thanks to the flexible, silicone Click Sleeves, Silk primax Click models offer a precise fit, binaural OneMic directionality and full wireless connectivity. This makes them ideal for first-time wearers looking for a discreet on-the-spot solution without any fuss.



One click and you're done.
Click Sleeves simply click on to Silk primax hearing aids and sit securely in the ear.

sDemo

Click ITE

Model	
	Click CIC
	
Performance	
Matrix (2 ccm)	50/113 dB
Battery	
Size	10
Power consumption	1.1 mA

Ultra-discreet. Custom fit. Effortless.

Key features

- Binaural OneMic directionality
- Small size up to high power matrix
- Full wireless and remote control options for all models
- Power matrix in all models, from IICs to ITEs
- Clinically proven to reduce hearing effort*



TeleCare 2.0 ready.



Audio streaming via easyTek and easyTek App.



With Notch Therapy for tinnitus treatment.

* See legal note on final page of this document.



Performance levels

7px

5px


3px

2px


1px

Insio primax.

Tailor-made for demanding wearers, Insio primax™ from Signia offers the ultimate in discretion and performance. Thanks to their binaural OneMic directionality function, even the small CIC models with a single microphone deliver superior directionality to the front. All models can be controlled remotely for a maximum of discretion.



CanalGrip®
The CanalGrip® is a durable pattern for retention in difficult ear canals.



Concha Lock
The Concha Lock provides a precise and stable retention. The transparent color supports high discretion.

ITE

Model	IIC	CIC	ITC	ITE
Performance				
Matrix (2 ccm)	65/124 dB ¹	65/124 dB ¹	65/124 dB ¹	65/124 dB ¹
Battery				
Size	10	10	312/10	13/312
Power consumption	1.1 mA ¹	1.1 mA ¹	1.1 mA ¹	1.1 mA ¹

(1) Max. values

Proven technology, outstanding comfort.



Key features

- For all levels of hearing loss
- Ergonomic design for a high level of wearing comfort
- Remotely controllable for a maximum of convenience and discretion



Compatible with touchControl App.



TeleCare 2.0 ready.







Direct Audio Input

Essential

Intuis 3 BTEs.

All Signia’s Intuis 3 BTE models offer proven technology with 6 channel MPO for the most advanced fitting options and for more individual adjustments. Suitable for all levels of hearing loss, the easy to handle Intuis 3 BTEs are equipped with automatic directional microphones for a precise hearing experience, offer remote control options, and are TeleCare ready for smoother, faster trials and higher customer satisfaction.

BTE				
Model	Intuis 3 S	Intuis 3 M	Intuis 3 P	Intuis 3 SP
Performance				
Matrix (2 ccm)	LifeTube 45/124 dB Earhook 55/124 dB	LifeTube 53/125 dB Earhook 60/130 dB	Earhook 70/134 dB	Earhook 80/140 dB
Battery	312	13	13	675
Size	0.9 mA	1.0 mA	1.4 mA	2.4 mA
Power consumption				

Key features

- For mild to severe hearing loss
- 6 channel MPO for the most advanced fitting options and more individual adjustments
- Remotely controllable for a maximum of discretion and convenience



Compatible with touchControl App.



TeleCare 2.0 ready.



Essential

Intuis 3 Click ITEs and ITEs.

With Intuis 3 Click ITEs and tailor-made ITEs, you can offer your customers the full range of in-the-ear models: No matter if they wish for a solution that fits right away or prefer a customized model, Intuis 3 is the perfect match. Equipped with proven technology, all models offer a precise hearing experience, are remotely controllable for a maximum of discretion and TeleCare ready for successful trials.

Click ITE/ITE		
Model	<div>Intuis 3 Click CIC, Click ITC</div> <div></div>	<div>Intuis 3 CIC, ITC, ITE</div> <div></div>
Performance		
Matrix (2 ccm)	Click CIC 50/113 dB Click ITC 55/116 dB	CIC 65/124 dB ¹ ITC 65/124 dB ¹ ITE 65/124 dB ¹
Battery		
Size	Click CIC 10, Click ITC 312	CIC 10, ITC 312, ITE 13
Power consumption	1.1 mA	1.1 mA ¹

(1) Max. values

Make your customers the center of attention.

Key features

- A full range of BTE, Click ITE and ITE models
- For all levels of hearing loss
- Remotely controllable for a maximum of convenience and discretion



Compatible with touchControl App.



TeleCare ready.




Lotus

The Prompt family.



Signia's Prompt family offers outstanding wearing comfort and ease of use for every form of hearing loss. All models offer remote control options via the touchControl App and the miniPocket™ and are a convenient solution for almost all your customers.



BTE

Model	Prompt S	Prompt P	Prompt SP
			
Performance			
Matrix (2 ccm)	ThinTube 45/124 dB Earhook 55/124 dB	Earhook 70/134 dB	Earhook 80/140 dB
Battery			
Size	312	13	675
Power consumption	0.9 mA	1.4 mA	2.4 mA

Click ITE/ITE

Model	Prompt Click CIC, Click ITC	Prompt CIC, ITC, ITE
		
Performance		
Matrix (2 ccm)	Click CIC 50/113 dB Click ITC 55/116 dB	CIC 65/124 dB ¹ ITC 65/124 dB ¹ ITE 65/124 dB ¹
Battery		
Size	Click CIC 10, Click ITC 312	CIC 10, ITC 312/10, ITE 13/312
Power consumption	1.1 mA	1.1 mA ¹

(1) Max. values

Take hearing to a new level.

Key features

- Quick and easy fittings everywhere via the MobileFitting App
- For mild to profound hearing loss
- Ergonomic shape for a high level of wearing comfort



Compatible with touchControl App.



Compatible with MobileFitting App.




Lotus

Run.



Easy to fit in just a few minutes with the convenient MobileFitting App, Run from Signia is ideal for on-the-spot fittings. Incredibly comfortable thanks to its ergonomic shape, the convenient and easy-to-use Run delivers an outstanding hearing experience.



BTE

Model	Run S	Run P	Run SP
			
Performance	ThinTube 45/124 dB Earhook 55/124 dB	Earhook 70/132 dB	Earhook 79/140 dB
Battery			
Size	312	13	675
Power consumption	0.9 mA	1.3 mA	2.2 mA

Click ITE

Model	Run Click CIC	Run Click ITC
		
Performance		
Matrix (2 ccm)	50/113 dB	55/116 dB
Battery		
Size	10	312
Power consumption	1.1 mA	1.1 mA

Provide an outstanding hearing experience.

- Key features
- A high level of amplification for severe to profound hearing loss
 - Ergonomic shape for a high level of wearing comfort
 - Simple fittings everywhere via the rocker switch

Lotus

Fun.

Equipped with Signia’s proven hearing aid technology, Fun covers all essential customer requirements. It offers superb sound quality for severe to profound hearing loss, excellent wearing comfort and is simply fitted by pushing the handy rocker switch.



BTE		
Model		
	<div>Fun P</div> 	<div>Fun SP</div> 
Performance		
Matrix (2 ccm)	Earhook 70/132 dB	Earhook 79/140 dB
Battery		
Size	13	675
Power consumption	1.2 mA	2.2 mA

Deliver a powerful listening experience.

Key features


- Simple fittings everywhere via the rocker switch
- For mild to severe hearing loss
- Easy handling

Lotus

Fast.

Fast BTE hearing aids offer a high level of performance and deliver just the right amount of amplification for customers with severe hearing loss. They are easily adjusted via the rocker switch for simplified fitting. Sleek and comfortable with a range of features, Signia's Fast ensures an outstanding hearing experience suitable for every customer.



BTE	
Model	Fast P
	
Performance	
Matrix (2 ccm)	Earhook 62/130 dB
Battery	
Size	13
Power consumption	1.1 mA


Sharpening your customers' senses.

SecureEar.

SecureEar™ is an active digital hearing protection that prevents the wearer from excessive noise and loud sounds. Customized for a perfect fit, it ideally combines hearing protection with the ability to communicate effectively in all environments. With SecureEar, the wearer can also receive cell phone or CB calls at any time thanks to the integrated T-coil.







SecureEar with T-coil: Protects hearing and also allows the user to communicate without any constraints.

ITE	
Model	SecureEar
	
Performance	
Matrix (2 ccm)	85 dB
Battery	
Size	312
Battery life	~ 230 h

Apps feature overview.

Remote support, audio streaming and remote control at the tip of your clients' fingers: Most Signia hearing aids are compatible with useful apps for Android and iOS smartphones. This overview shows you the features of each app in detail.

	Remote support	Audio streaming and remote control	Remote control	Audio streaming and remote control
	 myHearing App	 myControl App	 touchControl App	 easyTek App ^A
Availability	Apple App Store/ Google Play Store	Apple App Store/ Google Play Store	Apple App Store/ Google Play Store	Apple App Store/ Google Play Store
Use case & compatibility	Signia TeleCare	Nx hearing aids with Bluetooth connectivity	Styletto, Nx hearing aids without Bluetooth connectivity, all primax, Intuis 3, Prompt and Run hearing aids	easyTek required
Remote control features				
Transmission signal	Bluetooth High frequency coded signal	Bluetooth High frequency coded signal	High frequency coded signal	Bluetooth High frequency coded signal
Volume adjustment	●	●	●	●
Independent volume adjustment (R/L)	-	● ^B	-	-
Program change	●	●	●	●
Mute	●	-	●	● ^C
Remote microphone	-	● ^D	-	-
SoundBalance adjustment	●	●	●	●
Adjustable system sounds	-	●	-	-
Spatial Configurator adjustment	● ^E	● ^E	● ^E	● ^E
Automatic listening situations	-	● ^F	-	-
Power-on-Delay configuration	-	●	-	-
Battery status read-out	-	●	● ^G	●
Low battery level	-	●	-	-
Voice activity tracking	-	●	-	-
Sound exposure tracking	● ^H	●	-	-
Wearing time tracking	● ^H	●	-	-
User guide	●	-	-	-
Adaptive streaming volume	-	● ^D	-	-
Smartphone motion sensors usage	● ^H	●	-	-
Hearing aid airplane mode	-	●	-	-
Signia TeleCare features				
Usage statistics	● ^H	-	-	-
Daily satisfaction rating	● ^H	-	-	-
Hearing lessons	●	-	-	-
CareChat	●	-	-	-
Private label option	● ^I	-	-	-

A) Programmable via Connexx.

B) Configurable via Connexx.

C) Hearing aid microphones can be muted while streaming via easyTek device.

D) Not available in the app version on Google Play Store.

E) Available with hearing aids with two microphones each and at performance levels 7 and 5.

F) Available in binaural fittings only.

G) Output via the number of signal tones in the hearing aid.

H) Available with Nx hearing aids with Bluetooth connectivity.

I) After activation via customer service.

App overview.

Remote support



Via the **myHearing App**, your customers can experience all benefits of TeleCare to familiarize themselves with their new hearing aids in as little time as possible. They can follow special listening exercises, consult the user guide and connect to you on the go.

Audio streaming and remote control



The **myControl App** enables direct streaming and remote control for hearing aids using Bluetooth. It also optimizes hearing aid settings for each individual listening environment and lets your clients track their communication time and hearing performance.

Remote control



The **touchControl App** for hearing aids without Bluetooth enables discreet and convenient changes to hearing programs, volume, bass, and treble via the wearer's smartphone - without any extra hardware.

Audio streaming and remote control















With the **easyTek App™**, hearing aid wearers can discreetly control the easyTek via their smartphones for program or volume adjustments and audio streaming.



All apps are compatible with Android and iOS. Their download on the Apple App Store or Google Play Store is free of charge.

Apps, accessories and hearing aid compatibility.

This overview shows you which apps and accessories work with our hearing aids. Depending on the solution your customers choose, they can even enhance their hearing experience and convenience with this well-matching portfolio.

Hearing aid	Apps				CROS/ BiCROS	Streamers and transmitters				Remote controls		Charging stations
	myHearing App	myControl App	touchControl App	easyTek App	CROS	StreamLine TV	StreamLine Mic	easyTek	TV Transmitter	miniPocket	easyPocket	Charger
												
Stylecto New!	●	-	●	-	-	-	-	-	-	●	-	● ^A
Pure 10 Nx	●	-	●	-	-	-	-	-	-	●	-	-
Pure 312 Nx	●	●	-	-	● ^B	●	●	-	-	●	-	-
Pure 13 Nx	●	●	-	-	● ^B	●	●	-	-	●	-	-
Pure Charge&Go Nx	●	●	-	-	● ^B	●	●	-	-	●	-	● ^C
Motion 13 Nx	●	●	-	-	● ^B	●	●	-	-	●	-	-
Motion 13P Nx	●	●	-	-	● ^B	●	●	-	-	●	-	-
Motion Charge&Go Nx	●	●	-	-	● ^B	●	●	-	-	●	-	● ^C
Silk Nx	●	-	●	-	● ^D	-	-	-	-	●	-	-
Insio Nx	●	● ^E	● ^F	-	● ^{G,H}	● ^E	● ^E	-	-	●	-	-
Ace primax	●	-	●	-	-	-	-	-	-	●	-	-
Pure primax	●	-	●	●	● ^{I,J}	-	-	●	●	●	●	● ^{I,K}
Cellion primax	●	-	●	●	● ^{J,L}	-	-	●	●	●	●	● ^C
Carat primax	●	-	●	●	● ^J	-	-	●	●	●	●	● ^K
Motion SP primax	●	-	●	●	● ^{I,J}	-	-	●	●	●	●	-
Motion P primax	●	-	●	●	● ^{I,J}	-	-	●	●	●	●	-
Motion SX primax	●	-	●	●	● ^{I,J}	-	-	●	●	●	●	● ^{I,K}
Motion SA primax	●	-	●	●	● ^{I,J}	-	-	●	●	●	●	-
Motion S primax	●	-	●	●	-	-	-	●	●	●	●	-
Silk primax	●	-	●	●	● ^{I,J}	-	-	●	●	●	●	-
Insio primax	●	-	●	●	● ^{I,J}	-	-	●	●	●	●	-
Intuis 3	●	-	●	-	-	-	-	-	-	●	-	-
Prompt	●	-	●	-	-	-	-	-	-	●	-	-
Run	-	-	●	-	-	-	-	-	-	●	-	-

A) Compatible with Stylecto charging case.

D) Compatible with CROS Silk Nx.

G) Insio Nx IIC and CIC are compatible with CROS Silk Nx.

J) Compatible with CROS Pure.

B) Compatible with CROS Pure 312 Nx.

E) Compatible with Insio Nx ITC/ITE.

H) Insio Nx ITC and ITE are compatible with CROS Pure 312 Nx.

K) Compatible with eCharger.

C) Compatible with the Cellion/Pure Charge&Go/Motion Charge&Go charger.

F) Compatible with Insio Nx IIC/CIC.

I) Not available with 2px/1px.

L) Not available with 2px.

Accessory overview.

Streamers



Via **StreamLine TV**, TV sound can be fed directly into Bluetooth-equipped Signia Nx hearing aids. Wearers can even adjust the TV volume individually via the myControl App.



The **StreamLine Mic** transforms Signia Nx devices with Bluetooth into more than state-of-the-art hearing aids. They also become hands-free mobile headsets and wireless stereo headphones that stream highest quality audio into both ears from any Bluetooth device. Wearers can even adjust the streaming volume individually via the myControl App.



The **easyTek** wirelessly connects hearing aids to different Bluetooth enabled devices so sound is streamed directly into the hearing aids. Plus, it also acts as remote control for program or volume adjustments and is controllable via the easyTek App.

CROS/BiCROS



For people with unaidable hearing loss on one side. **CROS Pure Nx**, **CROS Silk Nx**, and **CROS Pure** receive sound on that side and transmit it to a hearing aid on the other side. **CROS Silk Nx** is the world's only **CROS CIC**.

Transmitter



The **TV Transmitter** connects TVs to the easyTek. This way, sound can be wirelessly streamed directly into the hearing aids.

Charging stations



The **chargers** recharge hearing aids so your clients don't need to worry about handling batteries.

Remote controls



Small enough to carry on a keyring, the **miniPocket** is the ideal solution for anyone who wants discreet, convenient hearing aid control without the need for a large remote control or smartphone.



The sleek **easyPocket™** with large buttons and an easy-to-read display even gives users with dexterity issues easy access to their hearing aids.







Color overview SLIM-RIC, RICs and BTEs.

Products	Colors									
	Beige	Granite	Grey	Brown	Black	Silver	Pearl White	Golden Blonde	Sandy Brown	Dark Champagne
	1	2	3	4	5	6	7	30	31	37
Pure 10 Nx	●	●	●		●	●	●	●	●	●
Pure 312 Nx	●	●	●		●	●	●	●	●	●
Pure 13 Nx	●	●	●		●	●	●	●	●	●
Pure Charge&Go Nx	●	●	●		●	●	●	●	●	●
Motion 13 Nx	●	●	●		●	●	●	●	●	●
Motion 13P Nx	●	●	●		●	●	●	●	●	●
Motion Charge&Go Nx	●	●	●		●	●	●	●	●	●
Ace primax	●	●	●		●	●	●	●	●	
Pure primax	●	●	●		●	●	●	●	●	●
Cellion primax	●	●	●		●	●	●	●	●	●
Carat primax	●	●	●		●	●	●	●	●	●
Motion primax	●	● ⁴	●		●	● ⁵		● ⁶	● ⁶	● ⁶
Intuis 3	●		●	●	●	●				
Prompt	●		●							
Run	●		●							
Fun	●		●							
Fast	●									

(2) Not available with 2px/1px
 (3) Not available with 2px
 (4) Not for Motion S primax
 (5) Not for Motion SP primax
 (6) Only for Motion SX, SA and P primax
 (7) Only for Motion SA and P primax
 (8) Only for Motion SX and SA primax
 (9) Only for Orion 2 RIC312

Color overview Click ITEs.

Standard colors for Click ITEs

Description	Type	Silk Nx	Silk primax	Intuis 3	Prompt	Run
 Blue/Mocha	left Click CIC	●	●			
 Red/Mocha	right Click CIC	●	●			
 Blue/Black	left Click CIC	●				
 Red/Black	right Click CIC	●				
 Mocha/Beige	Click CIC			●	●	●
 Tan/Beige	Click ITC			●	●	●

Further colors available on request.

© 2006 The Authors
Journal compilation © 2006 Blackwell Publishing Ltd

Candy Pink	Galactic Blue	Dark Granite	Deep Brown	Deep Red	Spirit	Elegance	Dark Champagne/Black	Sandy Brown/Black
60	61	66	72	73	C5	C12	HC2	HC4
		●	●					
		●	●					
		●	●					
		●	●				●	
		●	●					●
		●	●					●
		●	●				●	
		●		●	●			
		●	●	● ²	● ²			
		●	●	● ³				
		●	●	●	●			
● ^{2,7}	● ^{2,7}	●	●	● ^{2,4,5}	● ^{2,8}	● ^{2,8}		
		●						
			●					
			●					
			●					

Product	Colors		
	Dark Granite/ Silver	Cosmic Blue/ Rose Gold	Snow White/ Rose Gold
			
Styletto			










Color Flipper

Presenting RIC, BTE and ITE colors, and product-specific overview of available colors per product.

Item No. 03287-99C1-7600

Color overview ITEs.

Standard colors for custom-made ITEs

Description	Type	Insio Nx	Insio Nx ITC/ITE	Insio primax	Intuis 3	Prompt
 Beige	ITC, ITE		●	●	●	●
 Tan	ITC, ITE		●	●	●	●
 Mocha	CIC, ITC, ITE	●	●	●	●	●
 Brown	CIC, ITC, ITE	●	●	●	●	●
 Dark Brown	IIC, ITC, ITE	●	●	●	●	●
 Black	IIC, CIC	●				
 Deep Red	IIC, ITC, ITE			●*		

Further colors available on request.

* Insio primax 7px/5px/3px only.



Meeting all your needs:

The Signia hearing care professional support.

Everyday you try to meet and exceed your customers' hearing needs and wishes. In Signia you have a partner who fully supports you. Our broad range of hearing aid components, counseling and trial support as well as fitting tools enhances your workflow and provides you with all you need to find the best individual solution for your customers.

For detailed information, see our website signia-pro.com

Perfection is in the detail:

Components for hearing aids.



miniReceivers



Click Sleeves 2.0



Click Domes
and LifeTubes



Click Molds

From RIC miniReceivers™ to BTE earpieces and coupling options to Click Sleeves 2.0 for Click ITEs and RICs, we supply you with all necessary items so you can tailor the perfect hearing solution for your customers.

Hearing aid components

RICs*: miniReceivers S, M, P, HP
Click Sleeves 2.0 XS, S, M, L

BTEs: Click Domes™
Click Molds™
ThinTubes
LifeTubes
S-LifeTubes
LifeTips
Lotus Coupling
Earmolds

Click ITEs: Click Sleeves 2.0 XS, S, M, L

* Styletto is compatible to M-size receivers and Click Sleeves, only.

Reducing effort, increasing satisfaction:

Support for counseling, fittings and trials.



Connex 8
fitting software



Unity 3
mobile fitting unit system



SoundDemo Toolkit
incl. sDemo hearing aids



Counseling Suite
advisory tool

Signia offers you a broad range of counseling and fitting support – from programming software and tools to mobile fitting units to hearing aids for demonstration and advisory tools – up to comprehensive remote hearing care. This enhances your workflow and makes trials and fittings smoother and more pleasant for your customers.

Tools and support for counseling, fittings and trials

Connex (Update 8.5)
iScan™ II
Unity™ 3
ConnexAir
Noahlink
Noahlink Wireless
HI-PRO 2 USB
TeleCare
Counseling Suite
soundDemo Toolkit





Detailed feature overview **SLIM-RIC** and **RICs**.

Type	Styletto	Pure 10 Nx	Pure 312 Nx	Pure 13 Nx	Pure Charge&Go Nx
Performance levels	7/5/3	7Nx/5Nx/3Nx			
Picture	Powered by Signia Nx 				
Product details					
Matrix (2 ccm)	- M 60/119 dB - -	S 45/108 dB M 55/113 dB P 60/118 dB HP 65/123 dB	S 45/108 dB M 60/119 dB P 70/124 dB HP 75/130 dB	S 45/108 dB M 60/119 dB P 70/124 dB HP 75/130 dB	S 45/108 dB M 60/119 dB P 70/124 dB HP 75/130 dB
Battery size	Lithium-ion power cell	10	312	13	Lithium-ion power cell
Power consumption	1.4 mA	1.2 - 1.3 mA	1.3 - 1.5 mA	1.3 - 1.5 mA	1.3 - 1.5 mA
Features					
Frequency channels/ Compression controls	7 - 48/20 5 - 32/16 3 - 24/12 - -	7Nx - 48/20 5Nx - 32/16 3Nx - 24/12 - -	7Nx - 48/20 5Nx - 32/16 3Nx - 24/12 - -	7Nx - 48/20 5Nx - 32/16 3Nx - 24/12 - -	7Nx - 48/20 5Nx - 32/16 3Nx - 24/12 - -
Hearing programs	6	6	6	6	6
e2e wireless	● ¹³	-	● ¹³	● ¹³	● ¹³
Direct streaming	-	-	●	●	●
TwinPhone	● ¹⁵	-	● ¹⁶	● ¹⁶	● ¹⁷
TeleCare	● ¹⁸	● ¹⁸	● ¹⁹	● ¹⁹	● ¹⁹
Tinnitus therapies	●	●	●	●	●
CROS/BiCROS	-	-	● ²¹	● ²¹	● ²¹
Push button	-	○	-	-	●
Volume control	-	-	● ²⁴	● ²⁴	● ²⁴
Rechargeable	●	-	-	-	●
T-coil	-	-	-	○	-
Battery door on/off function	-	●	●	●	●
Audio input	-	-	● ²⁸	● ²⁸	● ²⁸
Modifications					
Color conversion kit	-	○	○	○	○
App/Accessory compatibility					
myControl App	-	-	●	●	●
myHearing App	●	●	●	●	●
touchControl App	●	●	-	-	-
easyTek/TV Transmitter	-	-	-	-	-
StreamLine TV/StreamLine Mic	-	-	●	●	●
miniPocket	●	●	●	●	●
easyPocket	-	-	-	-	-
ConnexxAir	●	-	-	-	-
Noahlink Wireless	-	-	●	●	●
Programming adapter/cable	-	size 10	size 312	size 13	-
Fitting set	-	miniReceiver 2.0	miniReceiver 2.0	miniReceiver 2.0	miniReceiver 2.0

- (2) Not available with 2px/1px
 (3) Not available with 2px
 (10) 6 hearing programs for 7px/5px/3px
 (11) 4 hearing programs for 2px/1px
 (12) 4 hearing programs for 2px
 (13) Ultra HD e2e
 with binaural audio signal processing
 (14) HD Binaural Hearing
 with binaural audio signal processing
 (15) To be activated via touchControl App
 or miniPocket

- (16) To be activated via rocker switch,
 myControl App or miniPocket
 (17) To be activated via push button,
 myControl App or miniPocket
 (18) TeleCare 2.0
 (19) TeleCare 3.0
 (20) Not for 1px
 (21) With CROS Pure 312 Nx
 (22) With CROS Pure

Overview SLIM-RIC/RICs

Ace primax	Pure primax	Cellion primax	Carat primax
7px/5px/3px/2px/1px		7px/5px/3px/2px	7px/5px/3px
			
S 45/108 dB M 55/113 dB P 60/118 dB HP 65/123 dB	S 45/108 dB M 60/119 dB P 70/124 dB HP 75/130 dB	S 45/108 dB M 60/118 dB P 70/124 dB HP 75/130 dB	S 45/108 dB M 60/119 dB P 70/124 dB HP 75/130 dB
10	312	Lithium-ion power cell	13
0.9 - 1.0 mA	0.9 - 1.3 mA	0.9 - 1.1 mA	0.9 - 1.1 mA
7px - 48/20 5px - 32/16 3px - 24/12 2px - 16/8 1px - 16/8	7px - 48/20 5px - 32/16 3px - 24/12 2px - 16/8 1px - 16/8	7px - 48/20 5px - 32/16 3px - 24/12 2px - 16/8 -	7px - 48/20 5px - 32/16 3px - 24/12 -
6 ^{10/4} ¹¹	6 ^{10/4} ¹¹	6 ^{10/4} ¹²	6
-	● ¹⁴	● ¹⁴	● ¹⁴
-	-	-	-
-	● ²	● ³	●
● ¹⁸	● ¹⁸	● ¹⁸	● ¹⁸
● ²⁰	● ²⁰	●	●
-	● ^{2, 22}	● ^{3, 22}	● ²²
●	● ²³	●	● ²⁴
● ²⁵	● ²³	● ²⁵	● ²⁴
-	● ^{2, 26}	●	● ²⁷
-	●	●	●
●	●	-	●
-	○ ²⁹	○ ²⁹	○ ²⁹
○	○	○	○
-	-	-	-
●	●	●	●
●	●	●	●
-	●	●	●
-	-	-	-
●	●	●	●
-	●	●	●
-	●	●	●
-	-	-	-
size 10	size 312	CS 44 cable	size 13
miniReceiver 2.0	miniReceiver 2.0	miniReceiver 2.0	miniReceiver 2.0

(23) Exchangeable controls – programmable via Connexx

(24) Programmable via Connexx

(25) With touchControl App/miniPocket

(26) Only recommended with miniReceiver S

(27) Reduced battery life

(28) Via direct streaming

(29) With easyTek

● Standard

○ Optional







Detailed feature overview **BTEs**.

Type	Motion 13 Nx	Motion 13P Nx	Motion Charge&Go Nx	Motion SP primax	Motion P primax
Performance levels	7Nx/5Nx/3Nx			7px/5px/3px/2px/1px	
Picture					
Product details					
Matrix (2 ccm) ThinTube/ LifeTube/ Earhook	56/126 dB - 63/134 dB	65/129 dB - 77/135 dB	56/126 dB - 63/134 dB	- 82/140 dB	65/129 dB - 77/135 dB
Battery size	13	13	Lithium-ion power cell	675	13
Power consumption	1.3 - 1.5 mA	1.3 - 1.5 mA	1.3 - 1.5 mA	2.1 mA	1.1 - 1.3 mA
Features					
Frequency channels/ Compression controls	7Nx - 48/20 5Nx - 32/16 3Nx - 24/12 -	7Nx - 48/20 5Nx - 32/16 3Nx - 24/12 -	7Nx - 48/20 5Nx - 32/16 3Nx - 24/12 -	7px - 48/20 5px - 24/16 3px - 24/12 2px - 16/8 1px - 16/8	7px - 48/20 5px - 32/16 3px - 24/12 2px - 16/8 1px - 16/8
Hearing programs	6	6	6	6 ^{10/4 11}	6 ^{10/4 11}
e2e wireless	● ¹³	● ¹³	● ¹³	● ¹⁴	● ¹⁴
Direct streaming	●	●	●	-	-
TwinPhone	● ¹⁶	● ¹⁶	● ¹⁷	● ²	● ²
TeleCare	● ¹⁹	● ¹⁹	● ¹⁹	● ¹⁸	● ¹⁸
Tinnitus therapies	●	●	●	● ²⁰	● ²⁰
CROS/BiCROS	● ²¹	● ²¹	● ²¹	● ^{2, 22}	● ^{2, 22}
Push button	-	-	●	●	-
Volume control	● ²⁴	● ²⁴	● ²⁴	● ²⁴	● ²⁴
Rechargeable	-	-	●	-	-
T-coil	○	○	-	●	●
Battery door on/off function	●	●	-	●	●
Audio input	● ²⁸	● ²⁸	● ²⁸	○ ³²	○ ³²
Modifications					
Color conversion kit	○	○	○	-	○
Eyeglass adapter	-	-	-	-	-
App/Accessory compatibility					
myControl App	●	●	●	-	-
myHearing App	●	●	●	●	●
touchControl App	-	-	-	●	●
easyTek/TV Transmitter	-	-	-	●	●
StreamLine TV/StreamLine Mic	●	●	●	-	-
miniPocket	●	●	●	●	●
easyPocket	-	-	-	●	●
ConnexxAir	-	-	-	●	●
Noahlink Wireless	●	●	●	-	-
Programming adapter/cable	size 13	size 13	-	CS 44 cable	size 13
Fitting set	ThinTube	ThinTube	ThinTube	-	ThinTube

- (2) Not available with 2px/1px
 (10) 6 hearing programs for 7px/5px/3px
 (11) 4 hearing programs for 2px/1px
 (13) Ultra HD e2e
 with binaural audio signal processing
 (14) HD Binaural Hearing
 with binaural audio signal processing
 (16) To be activated via rocker switch,
 myControl App or miniPocket

- (17) To be activated via push button,
 myControl App or miniPocket
 (18) TeleCare 2.0
 (19) TeleCare 3.0
 (20) Not for 1px
 (21) With CROS Pure 312 Nx
 (22) With CROS Pure

Overview BTEs

Motion SX/SA primax	Motion S primax	Intuis 3 S	Intuis 3 M	Intuis 3 P	Intuis 3 SP
7px/5px/3px/2px/1px	3px/2px/1px	Essential			
					
SX SA					
53/126 dB - 60/130 dB	- 45/124 dB 55/124 dB	- 45/124 dB 55/124 dB	- 53/125 dB 60/130 dB	- 70/134 dB	- 80/140 dB
13	312	312	13	13	675
1.0 mA	0.8 mA	0.9 mA	1.0 mA	1.4 mA	2.4 mA
7px - 48/20 5px - 32/16 3px - 24/12 2px - 16/8 ³⁰ 1px - 16/8 ³⁰	- - 3px - 24/12 2px - 16/8 1px - 16/8	12/6	12/6	12/6	12/6
6 ^{10/4} 11	6 ^{31/4} 11	4	4	4	4
● ¹⁴	● ¹⁴	-	-	-	-
-	-	-	-	-	-
● ²	-	-	-	-	-
● ¹⁸	● ¹⁸	● ¹⁸	● ¹⁸	● ¹⁸	● ¹⁸
● ²⁰	● ²⁰	-	-	-	-
● ^{2, 22}	-	-	-	-	-
○ ²⁴	● ²⁴	● ²⁴	○ ²⁴	●	●
● ²⁴	○ ²⁴	-	● ²⁴	● ²⁴	● ²⁴
● ²	-	-	-	-	-
●	-	-	●	●	●
●	●	●	●	●	●
○ ²⁹ ○ ³²	○ ²⁹	-	○ ³³	○ ³³	○ ³³
○	○	○	○	○	-
-	-	-	-	-	-
-	-	-	-	-	-
●	●	●	●	●	●
●	●	●	●	●	●
●	●	-	-	-	-
-	-	-	-	-	-
●	●	●	●	●	●
●	●	-	-	-	-
●	●	-	-	-	-
size 13	size 312	size 312	size 13	Programming shoe/ CS 44 cable	Programming shoe/ CS 44 cable
ThinTube	LifeTube	LifeTube	LifeTube	-	-

(24) Programmable via Connexx

(28) Via direct streaming

(29) With easyTek

(30) For Motion SA primax only

(31) 6 hearing programs for 3px

(32) With battery compartment modification and easyTek

(33) With battery compartment modification

● Standard

○ Optional

Detailed feature overview **BTEs**.

Type	Prompt S	Prompt P	Prompt SP	Run S
Performance level	Lotus			
Picture				
Product details				
Matrix (2 ccm)	45/124 dB - 55/124 dB	- - 70/134 dB	- - 80/140 dB	45/124 dB - 55/124 dB
Battery size	312	13	675	312
Power consumption	0.9 mA	1.4 mA	2.4 mA	0.9 mA
Features				
Frequency channels/ Compression controls	8/4	8/4	8/4	8
Hearing programs	4	4	4	6 preconfigured sound profiles
e2e wireless	-	-	-	-
Direct streaming	-	-	-	-
Push button	●	●	●	-
Volume control	-	● ⁽²⁴⁾	● ⁽²⁴⁾	●
T-coil	-	●	●	-
Battery door on/off function	●	●	●	●
Audio input	-	○ ⁽³³⁾	○ ⁽³³⁾	-
App/Accessory compatibility				
touchControl App	●	●	●	●
easyTek/TV Transmitter	-	-	-	-
miniPocket	●	●	●	●
ConnexxAir	-	-	-	-
Noahlink Wireless	-	-	-	-
Programming adapter/cable	size 312	Programming shoe/ CS 44 cable	Programming shoe/ CS 44 cable	MobileFitting App
Fitting set	LifeTube	Comfort acoustic coupling	Comfort acoustic coupling	Comfort acoustic coupling






(24) Programmable via Connexx

(33) With battery compartment modification

● Standard

○ Optional

Overview BTEs

Run P	Run SP	Fun P	Fun SP	Fast P
Lotus				
				
-	-	-	-	-
70/132 dB	79/140 dB	70/132 dB	79/140 dB	62/130 dB
13	675	13	675	13
1.3 mA	2.2 mA	1.2 mA	2.2 mA	1.1 mA
8	8	6	6	4
6 preconfigured sound profiles	6 preconfigured sound profiles	4 preconfigured sound profiles	4 preconfigured sound profiles	3 preconfigured sound profiles
-	-	-	-	-
-	-	-	-	-
-	-	-	-	-
●	●	●	●	●
-	-	-	-	-
●	●	●	●	●
-	-	-	-	-
●	●	-	-	-
-	-	-	-	-
●	●	-	-	-
-	-	-	-	-
-	-	-	-	-
MobileFitting App	MobileFitting App	Fitting via rocker switch	Fitting via rocker switch	Fitting via rocker switch
Comfort acoustic coupling	Comfort acoustic coupling	Comfort acoustic coupling	Comfort acoustic coupling	Simple acoustic coupling

Detailed feature overview Click ITEs and ITEs.

Type	Silk Nx Click CIC	Insio Nx IIC/CIC	Insio Nx ITC	Insio Nx ITE	Silk primax Click CIC	Insio primax IIC/CIC
Performance levels	7Nx/5Nx/3Nx				7px/5px/3px/2px/1px	
Picture			 New!	 New!		
Product details						
Matrix (2 ccm)	50/113 dB	S 50/113 dB M 55/118 dB P 65/124 dB	S 50/113 dB M 55/118 dB P 65/124 dB	M 55/118 dB P 65/124 dB	50/113 dB	65/124 dB ¹
Battery size	10	10	312	312	10	10
Power consumption	1.3 mA	1.3 mA ¹	1.3 mA ¹	1.3 mA ¹	1.1 mA	1.1 mA ¹
Features						
Frequency channels/ Compression controls	7Nx - 48/20 5Nx - 32/16 3Nx - 24/12 - -	7Nx - 48/20 5Nx - 32/16 3Nx - 24/12 - -	7Nx - 48/20 5Nx - 32/16 3Nx - 24/12 - -	7Nx - 48/20 5Nx - 32/16 3Nx - 24/12 - -	7px - 48/20 5px - 32/16 3px - 24/12 2px - 16/8 1px - 16/8	7px - 48/20 5px - 32/16 3px - 24/12 2px - 16/8 1px - 16/8
Hearing programs	6	6	6	6	6 ^{10/4} ¹¹	6 ^{10/4} ¹¹
e2e wireless	○ ¹³	○ ¹³	○ ¹³	○ ¹³	○ ¹⁴	○ ¹⁴
Direct streaming	-	-	●	●	-	-
TwinPhone	● ¹⁵	● ³⁴	● ³⁵	● ³⁵	-	-
TeleCare	● ¹⁸	● ¹⁸	● ¹⁹	● ¹⁹	● ¹⁸	● ¹⁸
Tinnitus therapies	●	●	●	●	● ²⁰	● ²⁰
CROS/BiCROS	● ³⁶	● ³⁶	● ²¹	● ²¹	● ^{2, 22}	● ^{2, 22}
Push button	-	○ ³⁷	○	○	-	○ ³⁴
Volume control	○ ²⁵	○ ²⁵	○ ³⁸	○ ³⁸	○ ^{25, 29}	○ ^{25, 29}
T-coil	-	-	-	-	-	-
Battery door on/off function	●	●	●	●	●	●
Audio input	-	-	○ ²⁸	○ ²⁸	○ ²⁹	○ ²⁹
Optivent	-	○	○	○	-	○
App/Accessory compatibility						
myControl App	-	-	●	●	-	-
myHearing App	●	●	●	●	●	●
touchControl App	●	●	-	-	●	●
easyTek/TV Transmitter	-	-	-	-	●	●
StreamLine TV/StreamLine Mic	-	-	●	●	-	-
miniPocket	●	●	●	●	●	●
easyPocket	-	-	-	-	●	●
ConnexxAir	● ³⁹	● ³⁹	-	-	● ³⁹	● ³⁹
Noahlink Wireless	-	-	●	●	-	-
Programming adapter/cable	Flex connector	Flex connector ³⁷	Flex connector	Flex connector	Flex connector	Flex connector

(1) Max. values

(2) Not available with 2px/1px

(10) 6 hearing programs for 7px/5px/3px

(11) 4 hearing programs for 2px/1px

(13) Ultra HD e2e
with binaural audio signal processing(14) HD Binaural Hearing
with binaural audio signal processing(15) To be activated via
touchControl App or miniPocket

(18) TeleCare 2.0

(20) Not for 1px

(21) With CROS Pure 312 Nx

(22) With CROS Pure








(25) With touchControl App/
miniPocket

(28) Via direct streaming

(29) With easyTek

(34) To be activated via
touchControl App, miniPocket
or optional push button

Overview Click ITEs/ITEs

Insio primax ITC	Insio primax ITE	Intuis 3 Click CIC	Intuis 3 Click ITC	Intuis 3 CIC	Intuis 3 ITC	Intuis 3 ITE
7px/5px/3px/2px/1px		Essential				
						
65/124 dB ¹	65/124 dB ¹	50/113 dB	55/116 dB	65/124 dB ¹	65/124 dB ¹	65/124 dB ¹
312/10	13/312	10	312	10	312	13
1.1 mA ¹	1.1 mA ¹	1.1 mA	1.1 mA	1.1 mA ¹	1.1 mA ¹	1.1 mA ¹
7px - 48/20 5px - 32/16 3px - 24/12 2px - 16/8 1px - 16/8	7px - 48/20 5px - 32/16 3px - 24/12 2px - 16/8 1px - 16/8	12/6	12/6	12/6	12/6	12/6
6 ¹⁰ /4 ¹¹	6 ¹⁰ /4 ¹¹	4	4	4	4	4
○ ¹⁴	○ ¹⁴	-	-	-	-	-
-	-	-	-	-	-	-
○ ²	○ ²	-	-	-	-	-
● ¹⁸	● ¹⁸	● ¹⁸	● ¹⁸	● ¹⁸	● ¹⁸	● ¹⁸
● ²⁰	● ²⁰	-	-	-	-	-
● ^{2, 22}	● ^{2, 22}	-	-	-	-	-
○	○	-	●	○	○	○
○ ^{25, 29}	○ ^{25, 29}	○ ²⁵	○ ²⁵	○ ²⁵	○ ²⁵	○ ²⁵
○	○	-	-	-	○	○
●	●	●	●	●	●	●
○ ²⁹	○ ²⁹	-	-	-	-	-
○	○	-	-	○	○	○
-	-	-	-	-	-	-
●	●	●	●	●	●	●
●	●	●	●	●	●	●
●	●	-	-	-	-	-
-	-	-	-	-	-	-
●	●	●	●	●	●	●
●	●	-	-	-	-	-
● ³⁹	● ³⁹	-	-	-	-	-
-	-	-	-	-	-	-
Flex connector	Flex connector	Flex connector	Flex connector	Flex connector	Flex connector	Flex connector

(35) To be activated via myControl App, miniPocket or optional push button

(36) With CROS Silk Nx






(37) Not available for IIC

(38) With myControl App or miniPocket

(39) Only if e2e wireless is available

- Standard
- Optional

Detailed feature overview **Click ITEs** and **ITEs**.

Type	Prompt Click CIC	Prompt Click ITC	Prompt CIC	Prompt ITC	Prompt ITE
Performance levels	Lotus				
Picture					
Product details					
Matrix (2 ccm)	50/113 dB	55/116 dB	65/124 dB ⁽¹⁾	65/124 dB ⁽¹⁾	65/124 dB ⁽¹⁾
Battery size	10	312	10	312/10	13/312
Power consumption	1.1 mA	1.1 mA	1.1 mA ⁽¹⁾	1.1 mA ⁽¹⁾	1.1 mA ⁽¹⁾
Features					
Frequency channels/ Compression controls	8/4	8/4	8/4	8/4	8/4
Hearing programs	1	4	4	4	4
e2e wireless	-	-	-	-	-
Direct streaming	-	-	-	-	-
TwinPhone	-	-	-	-	-
TeleCare	-	-	-	-	-
Tinnitus therapies	-	-	-	-	-
CROS/BiCROS	-	-	-	-	-
Push button	-	●	○	○	○
Volume control	-	-	○ ⁽²⁵⁾	○ ⁽²⁵⁾	○ ⁽²⁵⁾
Rechargeable	-	-	-	-	-
T-coil	-	-	-	○	○
Battery door on/off function	●	●	●	●	●
Audio input	-	-	-	-	-
Optivent	-	-	○	○	○
App/Accessory compatibility					
touchControl App	●	●	●	●	●
miniPocket	●	●	●	●	●
ConnexxAir	-	-	-	-	-
Noahlink Wireless	-	-	-	-	-
Programming adapter/cable	Flex connector	Flex connector	Flex connector	Flex connector	Flex connector



(1) Max. values

(25) With touchControl App/miniPocket

● Standard

○ Optional

Overview Click ITEs/ITEs

Run Click CIC	Run Click ITC
Lotus	
	
50/113 dB	55/116 dB
10	312
1.1 mA	1.1 mA
8	8
6 preconfigured sound profiles	6 preconfigured sound profiles
-	-
-	-
-	-
-	-
-	-
-	-
-	-
-	-
-	-
●	●
-	-
-	-
●	●
●	●
-	-
-	-
MobileFitting App	MobileFitting App

The information in this document contains general descriptions of the technical options available, which do not always have to be present in individual cases and are subject to change without prior notice.

The Bluetooth word mark and logos are owned by the Bluetooth SIG, Inc. and any use of such marks by Sivantos is under license. Other trademarks and trade names are those of their respective owners.

Android and Google Play are trademarks of Google Inc.
Apple App Store is a trademark of Apple Inc.

iPhone is a trademark of Apple Inc., registered in the U.S. and other countries.

* Study conducted at the University of Northern Colorado, 2015, examined the effectiveness of the new features of primax by collecting and analyzing ongoing EEG data while subjects performed speech testing. For both primax features SpeechMaster and EchoShield, the objective brain behavior measures revealed a significant reduction in listening effort when the feature was activated.

Legal Manufacturer

For Siemens products:
Sivantos GmbH
Henri-Dunant-Strasse 100
91058 Erlangen
Germany

For Signia products:
Signia GmbH
Henri-Dunant-Strasse 100
91058 Erlangen
Germany

Sivantos GmbH is a Trademark
Licensee of Siemens AG.

Signia GmbH is a Trademark
Licensee of Siemens AG.

Order No. 03635-99C1-7600 | Printed in Germany
©08.2018 Signia GmbH

