

Phonak Audéo B is the comprehensive RIC portfolio for mild to severe hearing losses. Phonak Audéo B includes six designs, three external receivers and four performance levels, including a noise generator for tinnitus management. Powered by the innovative technology of the Phonak Belong platform, Phonak Audéo B hearing aids are the ideal solution for those clients that prefer best hearing performance in flexible and discreet RIC housings.

For more information visit www.phonakpro-us.com



Product description

Troduct acse						9		
			Audéo B-10 B90/B70/B50/B30	Audéo B-312 B90/B70/B50/B30	Audéo B-312T B90/B70/B50/B30	Audéo B-13 B90/B70/B50/B30	Audéo B-Direct B90/B70/B50	Audéo B-R B90/B70/B50
Battery size			10	312	312	13	13	
Push button (Program/Volume)			•	•	•	•	•	•
Volume control					-	•		
Telecoil					•	•		
EasyPhone			•	•	•	•		•
Bluetooth®					-		•	
Nano coating			•	•	•	•	•	•
IP rating			IP681	IP681	IP681	IP681	IP681	IP681
Dimensions	LxWxD			26.6x 10.8x 6.9 mm (1.05x 0.43 x 0.27")	28.7 x 12.3 x 6.8 mm (1.13 x 0.48 x 0.27")		27.2 x 10.9 x 7.3 mm (1.07 x 0.43 x 0.3")	29.7 x 12.9 x 7.9 mm (1.17 x 0.51 x 0.31")
	Weight		1.0 g (0.04 oz)	1.3 g (0.05 oz)	1.5 g (0.05 oz)	2.2 g (0.08 oz)	1.2 g (0.04 oz)	2.4 g (0.09 oz)
Battery life ²	No audio/media streaming							24 h
	5 h of audio/media streaming							20 h
	10 h of audio/media steaming							16 h
Charging times ³	0% to 100%				-			3 h
	0% to 80%							1.5 h
	30 min. fitting							10 min.
Standard xReceiver (xS)								
Max. Power Output (dB SPL)		2 cc coupler	111	111	111	111	111	111
	-	Ear simulator	122	122	122	122	122	122
Max. gain (dB)		2 cc coupler	46	46	46	46	46	46
		Ear simulator	56	56	56	56	56	56
Frequency range – Ear simulator	(Hz)		<100 - 9200	<100 - 9200	<100 - 9200	<100 - 9200	<100 - 9200	<100 - 9200
Working current (mA)			1.2	1.2	1.2	1.2	1.8	
Power xReceiver (xP)								
Maximum Power Output (dB SPL)	2 cc coupler	124	124	124	124	124	124
		Ear simulator	132	132	132	132	132	132
Max. gain (dB)		2 cc coupler	57	57	57	57	57	57
		Ear simulator	65	65	65	65	65	65
Frequency range – Ear simulator	(Hz)		<100 - 6400	<100 - 6400	<100 - 6400	<100 - 6400	<100 - 6400	<100 - 6400
Working current (mA)			1.2	1.2	1.2	1.2	1.8	

UltraPower xReceiver (xUP)

Working current (mA)

Max. Power Output (dB SPL)	2 cc coupler	130	130	130	130	130
	Ear simulator	136	136	136	136	136
Max. gain (dB)	2 cc coupler	66	66	66	66	66
	Ear simulator	72	72	72	72	72
Frequency range – Ear simulator (Hz)		<100 - 5800	<100 - 5800	<100 - 5800	<100 - 5800	<100 - 5800

1.3

1.3

1.3

2.0

IP68 indicates that the hearing aid is water resistant and dust tight. It survived continuous immersion in 1 meter of water for 60 minutes and 8 hours in a dust chamber as per the IEC60529 standard. No traces of dust were evident within the housing.
Actual battery life also depends on hearing loss, battery age and sound environment.
Actual charge time depends on the remaining battery life, but it will not exceed 3 hours.

Housing colors









Sandalwood



Chestnut



Champagne



Silver Gray



Graphite Gray



Velvet Black

T7

Alpine White



Beige


UltraPower xReceiver

•

Acoustic coupling

		0	Standard xReceiver	Power xReceiver	
Length	Available in fou	r different lengths (0, 1, 2, 3), left and right	•	•	
Dome	Open dome	Available in three sizes (S, M, L)	•	•	
	Closed dome	Available in three sizes (S, M, L)	•	•	
	Power dome	Available in three sizes (S, M, L)	•	•	
SlimTip*	Hard material	Hollow shell	•	•	
	Soft material	Solid shell	•		
cShell	Hard material		•	•	
*Available with ex	tra canal length option		dB	dB HL	dB I

^{*}Available with extra canal length option



Phonak Audéo B-R charging options



Phonak Charger Case RIC

The Phonak Charger Case is a charger, drying kit and protective hard case all in one. Additionally, it includes a cleaning tool in the robust and splash-proof housing designed for two Audéo B-R hearing aids.



Phonak Power Pack

hearing loss.

For charging Phonak Audéo B-R on the go, the Phonak Power Pack can be easily attached to the bottom of the Phonak Charger Case. It provides power for 7 full charges of two Audèo B-R hearing aids.



Phonak Mini Charger RIC The Phonak Mini Charger is a compact charger for two Audéo B-R hearing aids.

Please refer to the technical data sheet for the Phonak Charger Case, Phonak Power Pack and Phonak Mini Charger for further information.

Phonak Audéo B-Direct Bluetooth® versions

Direct calls

Bluetooth® 4.2 wireless technology and most older Bluetooth phones

Phonak Remote App

iOS 10.2 or higher and Android 6.0 or higher

Fitting

	Audéo B-10 B90/B70/B50/B30	Audéo B-312 B90/B70/B50/B30	Audéo B-312T B90/B70/B50/B30	Audéo B-13 B90/B70/B50/B30	Audéo B-Direct B90/B70/B50	Audéo B-R B90/B70/B50
Phonak Target	5.0 or higher	5.0 or higher	5.0 or higher	5.0 or higher	5.2 or higher	5.0 or higher
iCube II	•	•	•	•		•
Noahlink Wireless					•	
NOAHLink		•	•	•		
HI-PRO, HI-PRO2		•	•	•		

Wireless communication portfolio

Wireless acces	ssory	Audéo B-10 B90/B70/B50/B30	Audéo B-312 B90/B70/B50/B30	Audéo B-312T B90/B70/B50/B30	Audéo B-13 B90/B70/B50/B30	Audéo B-Direct B90/B70/B50	Audéo B-R B90/B70/B50
1	Phonak CROS B-312	•	•	•	•		
1	Phonak CROS B-13	•	•	•	•		
\(\rightarrow\)	Phonak DECT II	•	•	•	•		•
•	Phonak PilotOne II	•	•	•	•		•
0	EasyCall II	•	•	•	•		•
	Phonak ComPilot II	•	•	•	•		•
	Phonak ComPilot Air II	•	•	•	•		•
	Phonak TVLink II	•1	• 1	•1	•1		•1
	Phonak RemoteMic	•1	• 1	•1	•1		•1
	TV Connector					•	
Roger							
3	Roger 18				•		
\$	Roger X / AS18				•		
#	Roger X / Phonak ComPilot II	•	•	•	•		•
	Roger MyLink			•	•		
Apps							
	Phonak RemoteControl	•	•	•	•		•
	Phonak Remote					•	
						·	

For information on Roger microphones please visit www.phonakpro.com/roger

Performance levels

	Premium B90	Advanced B70	Standard B50	Essential B30*
AutoSense OS	Premium	Advanced	Standard	Essential
Calm Situation	•	•	•	•
Speech in Noise	•	•	•	•
Comfort in Noise	•	•	•	
Music	•	•		
Speech in Loud Noise **	•			
Speech in Car	•			
Comfort in Echo	•			
Additional programs				
Max. additional programs	5	5	5	5
Speech in Wind **	•			
Comfort in Echo	•			
Speech in Loud Noise **	•	•		
Speech in 360° **	•	•		
Speech in Noise	•	•	•	•
Calm Situation	•	•	•	•
Comfort in Noise	•	•	•	•
Music	•	•	•	•
Acoustic phone	•	•	•	•
Custom program	•	•	•	•
Streaming programs				
Audéo B-10, B-312, B-312T, B-13 a	nd B-R			
Max. streaming programs	4	4	4	4
Bluetooth audio + mic	•	•	•	•
Bluetooth phone / DECT + mic	•	•	•	•
RemoteMic / Roger	•	•	•	•
Audio jack	•	•	•	•
Audéo B-Direct				
Max. streaming programs	2	2	2	
TV Connector + mic	•	•	•	
Phone + mic	•	•	•	

Features	Premium B90	Advanced B70	Standard B50	Essential B30*
UltraZoom	Premium	Advanced	Standard	Essential
SNR-Boost	•	•	•	
FlexControl	•	•	•	
FlexVolume	•	•	•	
DuoPhone**	•	•	•	
SoundRecover2	•	•	•	•
User Preference Tuning	•	•	•	•
Real Ear Sound	•	•	•	
Finetuning channels	20	16	12	8
WhistleBlock	•	•	•	•
NoiseBlock	•	•	•	•
WindBlock	•	•		
EchoBlock	•			
SoundRelax	•	•		
QuickSync	•	•	•	•
AOV	•	•	•	•
Tinnitus Balance	•	•	•	•
auto Acclimatization	•	•	•	•
Environmental balance ***	•	•	•	

^{*} Not available for Phonak Audéo B-Direct and Audéo B-R ** Not available with Audéo B-Direct *** Only for Audéo B-Direct

Life is on

At Phonak, we believe that hearing well is essential to living life to the fullest. For more than 70 years, we have remained true to our mission by developing pioneering hearing solutions that change people's lives to thrive socially and emotionally. Life is on.

www.phonakpro-us.com



