







APOMINTON

Dog Wheelchair



-  Fully adjustable length ,width, and height
-  Allows pet to pee and poop while wheelchair is in use
-  For disabled pets due to dog age, conditions, injuries
-  Reabilitations for injured pets when they are healing



APOMINTON

INSTALLATION INSTRUCTIONS


A. Frame

B. Screws



C. Wheels

D. Dog leg fixed belt


E. Chest and back cover




Hip Height



Fix the harness E1 on frame with screw B



Remove the screw B to adjust the best size.
Hammock for the dog leg is adjustable .



Velcro of E2 and Straps of E3 can fix directly on the frame, straps size is adjustable.