

manuthera[®] 242

Master of Spine Techniques



Manuthera 242 isn't just any table –

it's a unique and versatile companion for the most demanding
professionals in manual therapy and patient care


LOJER[®]

manuthera242.com

How to utilize Manuthera 242?


Cervical spine traction

Access lateral aspects of cervical spine joints and ligaments. Treat and relieve tension in the stiff and compressed neck and shoulder region.




Thoracic spine manipulation

Manipulate and adjust thoracic spine and ribs for better range of motion and function. Tilting for gravity assistance, flexion and lateral bend from the Manuthera enables more specific and comfortable treatment. Use the dynamic motion for soft tissue techniques and gaining access to deeper tissues.



Lumbar spine

Smooth flexion distraction with side bend, rotation, and decompression options gives you limitless 3D opportunities for lumbar, intervertebral disc, pelvic, and sacrum area treatment. Add some inversion for a broader treatment angle, a stronger effect, treating both superficial and deeper layers of the tissue.




Discover other possibilities

Access free videos from global experts!

Access the free video Learning Hub platform by scanning the QR-code below and learn how to get the most out of Manuthera's versatile options. Sign up and you get access to hundreds of videos of different practitioners sharing their expertise with Manuthera 242 and success they've had with various patient cases completely free of charge.



Manuthera 242 dimensions



Technical data and features	
Height adjustment	51–106 cm
Weight	170 kg
Tilt Up / Down	Electrical 52°/32°
Tilt Up/Down function	Electronic
Lower body rotation	Yes
Lower body section movement (3D)	+21.5°/max -40° (vertical), 0..18° (horizontal)
Upper thorax -section movement (3D)	+15°/-13° (vertical), 0..18° (horizontal)
Total lock function for upper/lower body	Yes
Castors	Yes
Safe working load (SWL)	150 kg
Drop functions	Traction drop
Safety edge	Yes
Safety lock functions	Automatic (2 h) or Manual / Manually operated
Lying surface	Antimicrobial upholstery
U-shaped face element and breathing hole cover	Yes
Manipulation handle	Yes
Country of Origin	Finland
CE Marked	Yes The device is a CE Marked medical device according to MDR 2017/745 class 1)
Manufacturer	Lojer Oy

FENNIAPRIZE 17

HONOURABLE MENTION

