



REF EU-32041

4 Fr. 1 Lumen 20 cm catheter length .018 inch dia. spring-wire guide

# Pressure Injectable One-Lumen Midline

Contents:

- 1: TaperFree® Catheter: 4 Fr.(1.45 mm OD) x 20 cm Pressure Injectable Radiopaque Polyurethane with T-Port Connector, Extension Line Clamp, Contamination Guard and Stiffening Stylet
- 1: GlideThru™ Peel-Away Sheath: 4 Fr. x 2-3/4" (7 cm) Radiopaque over 4 Fr. Dilator
- 1: Spring-Wire Guide, Nitinol, Marked: .018" (0.46 mm) dia. x 17-3/4" (45 cm) (Straight Soft Tip on One End - Straight Stiff Tip on Other) with Arrow Advancer
- 1: Introducer Needle: Echogenic 21 Ga. x 2-3/4" (7 cm) TW
- 1: Syringe: 10 mL Luer-Lock
- 1: SharpsAway® II Locking Disposal Cup
- 1: Catheter Trimmer
- 1: Safety Scalpel: #11
- 1: Paper Tape Measure
- 1: Tourniquet
- 1: Dressing: GRIP-LOK®<sup>1</sup> Catheter Fixation Device
- 1: Valve: MicroCLAVE®<sup>2</sup> Neutral Displacement Connector

<sup>1</sup>A registered trademark of TIDI.  
<sup>2</sup>A registered trademark of ICU Medical, Inc.

All components are CE 2797 unless otherwise noted.

Warning: Read all package insert warnings, precautions, and instructions prior to use. Failure to do so may result in severe patient injury or death. [www.teleflex.com/IFU](http://www.teleflex.com/IFU)  
Not made with natural rubber latex.  
Fluid path components are non-pyrogenic.

| Lumen           | Priming Volume* (mL) | Gravity Flow Rate† (mL/hr) | Pump Flow Rate†† (mL/hr) | MAX Pressure Injection Flow Rate** (mL/sec) |
|-----------------|----------------------|----------------------------|--------------------------|---------------------------------------------|
| Yellow (18 Ga.) | 0.37                 | 1880                       | 9324                     | 5                                           |

\* Priming volumes are approximate and are done without accessories.  
† Flow rate values are approximate and are determined using deionized water at 100 cm head height.  
†† Pump flow rates are determined at maximum pump pressure of 10 psig and represent approximate flow capabilities.  
\*\* Pressure injection flow rates are determined using media of 11.8 centipoise (cP) viscosity, with 152 cm administration set tubing.

- enbghrcsdanl etfi frde el hu itlv
- ltno plptro sksl es sv tr



EU Authorized Representative and Importer:  
Teleflex Medical  
IDA Business and Technology Park  
Dublin Road, Athlone  
Co. Westmeath, Ireland  
Arrow International LLC  
Subsidiary of Teleflex Incorporated  
3015 Carrington Mill Blvd.  
Morrisville, NC 27560 USA

Teleflex®  
www.teleflex.com

LBL055050 R02 (2022-01)



REF EU-32041

4 Fr. 1 Lumen 20 cm catheter length .018 inch dia. spring-wire guide

# Pressure Injectable One-Lumen Midline