

CLEANING SOLUTION 1

Detergent-based solution for everyday cleaning procedures of heliXcyto and chips.

Dynamic Biosensors GmbH & Inc.
CS-100-1 v1.1



Product Description

Order Number: **CS-100-1**

Table 1. Contents and Storage Information

Material	Amount	Storage
Cleaning solution 1	100 mL	RT

The detergent-based Cleaning Solution 1 (CS 1) is **ready-to-use**.

CS 1 is used for regular cleaning of the **heliX^{cyto}** and **heliX^{cyto}** chip. Please follow the use of CS 1 with RB 1 or another suitable buffer or H₂O.

For research use only.

This product has a limited shelf life, please see expiry date on label.

Contact

Dynamic Biosensors GmbH
Perchtinger Str. 8/10
81379 Munich
Germany

Dynamic Biosensors, Inc.
300 Trade Center, Suite 1400
Woburn, MA 01801
USA

Order Information order@dynamic-biosensors.com
Technical Support support@dynamic-biosensors.com

www.dynamic-biosensors.com

Instruments and chips are engineered and manufactured in Germany.
©2024 Dynamic Biosensors GmbH | Dynamic Biosensors, Inc. All rights reserved.