

Nuance Dragon Medical One EDUCATION

DAX Copilot is interrupted if anything takes over the microphone, such as an incoming call. Do Not Disturb can be used manually or can be set to automatically run any time the PowerMic Mobile app is in use.

Manually Enable Do Not Disturb

- 1. From any screen, **swipe** to open the Control Center.
- 2. Tap Focus.
- 3. Tap **Do Not Disturb** from the Focus modes.
- 4. Do Not Disturb is enabled.
- Follow the steps again to disable Do Not Disturb after the patient encounter is completed.



Note: Calls/messages from favorite contacts may override the Do Not Disturb setting. To limit all interruptions, ensure there are no exceptions added to the Do Not Disturb feature.

Automatically Enable Do Not Disturb

- 1. From any screen, **swipe** to open the Control Center.
- 2. Tap Focus.
- 3. Tap the ellipses (...) to the right of Do Not Disturb.
- 4. Tap Settings.
- 5. Scroll down and tap Add Schedule.



- 6. Select App.
- 7. Find and select the **PowerMic Mobile** app.
- 8. Under **SET A SCHEDULE**, the PowerMic Mobile app is now listed as **On**.
 - Do Not Disturb will automatically enable while the PowerMic Mobile app is in use.
 - Do Not Disturb will automatically disable when the PowerMic Mobile app is minimized or closed.



