

## 17AD300 Nexgear Tango ankle.

|                                                              |             |      |
|--------------------------------------------------------------|-------------|------|
| Contact person                                               | Account no. | Date |
| <b>Customer</b>                                              |             |      |
| Company                                                      |             |      |
| Street                                                       |             |      |
| City/State                                                   | Postal Code |      |
| Email                                                        |             |      |
| Phone                                                        |             |      |
| PO Number                                                    |             |      |
| <b>Shipping address</b> (if different from customer address) |             |      |
| Company                                                      |             |      |
| Street                                                       |             |      |
| City/State                                                   | Postal Code |      |
| Email                                                        |             |      |

| Shipping Options | UPS Next Day | UPS Ground | UPS 2-Day | Other |
|------------------|--------------|------------|-----------|-------|
|------------------|--------------|------------|-----------|-------|

## User Information

Name \_\_\_\_\_ Weight \_\_\_\_\_

Age \_\_\_\_\_ Indication \_\_\_\_\_

| Side | Left | Right | Fitting Stage | Test Orthosis | Definitive Orthosis |
|------|------|-------|---------------|---------------|---------------------|
|------|------|-------|---------------|---------------|---------------------|

The size chosen depends on the patient's weight and the conditions of use. The 17AD300 *Nexgear Tango* ankle can be equipped with up to two modules. Please enter the required quantity in the field provided.



### Special Instructions/ Comments

[illegible]

Date \_\_\_\_\_ Place \_\_\_\_\_ Signature \_\_\_\_\_

# 17AD300 Nexgear Tango ankle. iFab order form.

## 17AD300 Nexgear Tango Joint Selections

If you selected 17AD300=16-T and 17AD300=20-T for the ankle joint, please select one of each of the following. If you selected another ankle joint, please disregard this section. If you have any questions, please contact clinical support.

### 1. Anterior Lateral Channel

Free Motion  
Stop Module  
Spring Module  
Reaction Module Firm  
Reaction Module Extra Firm  
Ultra Module

### 2. Anterior Medial Channel


Free Motion  
Stop Module  
Spring Module  
Reaction Module Firm  
Reaction Module Extra Firm  
Ultra Module

### 3. Posterior Lateral Channel

Free Motion  
Stop Module  
Spring Module  
Reaction Module Firm  
Reaction Module Extra Firm  
Ultra Module

### 4. Posterior Medial Channel

Free Motion  
Stop Module  
Spring Module  
Reaction Module Firm  
Reaction Module Extra Firm  
Ultra Module



1 Lamination bar

2 Foot stirrup

|                |                  |                    |
|----------------|------------------|--------------------|
| pc(s) 17LS3=14 | pc(s) 17LS3=14-T | pc(s) 17SF300=N-14 |
| pc(s) 17LS3=16 | pc(s) 17LS3=16-T | pc(s) 17SF300=N-16 |
| pc(s) 17LS3=20 | pc(s) 17LS3=20-T | pc(s) 17SF300=N-20 |



3 Base body

3.1 Lamination Dummy

|                     |                  |
|---------------------|------------------|
| pc(s) 17AD300=14-ST | pc(s) 17AD313=14 |
| pc(s) 17AD300=16-ST | pc(s) 17AD313=16 |
| pc(s) 17AD300=20-ST | pc(s) 17AD313=20 |
| pc(s) 17AD300=14-T  |                  |
| pc(s) 17AD300=16-T  |                  |
| pc(s) 17AD300=20-T  |                  |



4 Stop-module

5 Spring-module

|                  |                  |
|------------------|------------------|
| pc(s) 17AD301=14 | pc(s) 17AD302=14 |
| pc(s) 17AD301=16 | pc(s) 17AD302=16 |
| pc(s) 17AD301=20 | pc(s) 17AD302=20 |




6 Reaction-module

6a Reaction spring, strong

6b Reaction spring, extra strong

|                  |                  |                  |
|------------------|------------------|------------------|
| pc(s) 17AD303=14 | pc(s) 17AD305=14 | pc(s) 17AD306=14 |
| pc(s) 17AD303=16 | pc(s) 17AD305=20 | pc(s) 17AD306=20 |
| pc(s) 17AD303=20 |                  |                  |



7 Ultra-module

7a Ultra spring

|                  |                  |
|------------------|------------------|
| pc(s) 17AD304=16 | pc(s) 17AD307=20 |
| pc(s) 17AD304=20 |                  |

**Note:** If fabricating Ottobock bars and joints for the first-time, you will need to order the corresponding adjustment adapter and adapter sleeves.

Date \_\_\_\_\_ Place \_\_\_\_\_ Signature \_\_\_\_\_

17AD300 Nexgear Tango ankle.  
iFab order form.

Standard AFO Configuration

Includes: Pre-preg exterior frame, anterior shell, satin black finish, Dacron calf strap with chafe, one liner.  
\*Any special requests outside of our standard AFO configuration can be written in the special instruction/ comments sections of the first page.

Add-On Options (additional charges apply)

| Padding         | Material      |      | Ankle Strap                        |
|-----------------|---------------|------|------------------------------------|
| Malleolus pad   | Microcel Puff | 2 mm | Dacron strap and chafe             |
| Arch pad        | Bocklite      | 3 mm | Hook-and-loop lay over ankle strap |
| Calf liner      | Other         |      |                                    |
| Anterior shell  |               |      |                                    |
| Full foot shell |               |      |                                    |

If used as corrective padding, additional L-Code applies: L2275 x1  
Please check one:      Varus correction      Valgus correction

Casting Instructions

When casting for a dynamic AFO, accurate shoe heel height representation in the cast is critical and we recommend using a foot casting aid (part # 743A9) or equivalent. The line of progression is critical for placement of toe break in the fabrication and carbon layup of the foot section with the AFO. This line of progression is established in the impression with the patella and popliteal landmark, and should be noted if it is more than 7°. Please ensure the height of the impression is at or proximal to the patella to accurately model both landmarks. Make sure to end your wrap distal to the phalanges utilizing another roll of casting material, if necessary. The patient's thigh, knee, tibia, and foot must all be maintained in line of progression during casting. Your cast height has to be at least 1 cm higher than your finished orthosis height. We will not extend casts higher or fill in missing phalanges. Please make sure your cast removal strip is placed lateral to the proximal half of the tibial crest, and your model provides accurate representation of the patient's proximal tibial crest.

Metric Measurements Only

\_\_\_\_\_ mm

Effective Heel Height

Height of Heel Lift (if any)

\_\_\_\_\_ mm

Heel height (standard is 5mm)

Metric Measurements Only

Trim line height of AFO (standard is ~3cm below fib head)

\_\_\_\_\_

Ankle center line to bottom of heel

\_\_\_\_\_

Bottom of heel

Heel height

\_\_\_\_\_

Date \_\_\_\_\_ Place \_\_\_\_\_ Signature \_\_\_\_\_