


# Set up your *Concha SolHearing aids*



## Download the app

Go to the Apple App Store to download our app “Concha Labs”, using the QR code.



## Sign in to activate your hearing aids


Be in a quiet environment,  
have your iPhone ready.

## Hear clearly, live freely


If you would like to view the instructions  
for use online, please visit  
**conchalabs.com**

Concha tip: Follow the app onboarding instructions with care. Please do not put batteries in the hearing aids until instructed.

# What's in *your box*



a.



b.



c.



d.



e.



f.

a.1 pair of hearing aids with 7mm ear tips in a case

b.Five 6-count packs of disposable zinc-air 312 batteries

c.1 pair of ear plugs

d.1 additional pair of 7mm ear tips

e.2 pairs of 10mm ear tips





f.1 pack of wax filters

## Need help?

 [conchalabs.com](https://conchalabs.com)

- [conchalabs.com/conchasol-ifu](https://conchalabs.com/conchasol-ifu)
- [care@conchalabs.com](mailto:care@conchalabs.com)
- +1 (510) 544-6778

## Connect with us

-  [conchalabs.com](https://conchalabs.com)
-  [conchalabs](https://twitter.com/conchalabs)
-  [conchalabs](https://www.instagram.com/conchalabs)
-  Concha Labs