

Set up your *Concha SolHearing aids*



Download the app

Go to the Apple App Store to download our app “Concha Labs”, using the QR code.



Sign in to activate your hearing aids


Be in a quiet environment,
have your iPhone ready.

Hear clearly, live freely

If you would like to view the instructions
for use online, please visit
conchalabs.com

Concha tip: Follow the app onboarding instructions with care. Please do not put batteries in the hearing aids until instructed.

What's in *your box*



a.



b.



c.



d.



e.



f.

a.1 pair of hearing aids with 7mm ear tips in a case

b.Five 6-count packs of disposable zinc-air 312 batteries

c.1 pair of ear plugs

d.1 additional pair of 7mm ear tips

e.2 pairs of 10mm ear tips





f.1 pack of wax filters

Need help?

 conchalabs.com

- conchalabs.com/conchasol-ifu
- care@conchalabs.com
- +1 (510) 544-6778

Connect with us

-  conchalabs.com
-  [conchalabs](https://twitter.com/conchalabs)
-  [conchalabs](https://www.instagram.com/conchalabs)
-  [Concha Labs](https://www.facebook.com/ConchaLabs)