

Oticon Government Services

More potential than ever

life-changing
technology



NEW
STYLES



Oticon More
miniBTE R

Oticon More
miniBTE T

The Oticon More™ family is expanding with even more options designed to fit your patients' needs.

- Clinical research shows Oticon More gives the brain what it needs to make sense of sound
- NEW miniBTE styles join the existing miniRITE family to bring the superior sound quality of Oticon More to more patients
- Two-way hands free phone or video calls from select iPhone and iPad devices*

*Hands-free communication is available with iPhone 11 or later running iOS 15.2 or later, and iPad running iPadOS 15.2 or later. See which hearing aids and devices are compatible here: oticon.com/support/compatibility

oticon
life-changing technology

RITE STYLES

Opn S miniRITE	Oticon More miniRITE R* (lithium-ion)	Oticon More miniRITE T* Opn Play miniRITE T*	Oticon CROS** Oticon CROS PX** (rechargeable)
----------------	--	---	---

60 85 100 105

60 85 100 105


60 85 100 105



*Includes Telecoil and Toggle Switch

**Available with 60 or 85 receiver

RITE Fitting Ranges



MINIFIT RECEIVER AND EARPIECE OPTIONS



60

Receiver



85

Receiver



100

Receiver

Earpieces



60 85

MicroShell
Embedded
Receiver

Now available
with more
options



100 105

Power
Receiver Mold
Interchangeable
Receiver Wire



MicroMold



Hollow Canal
LiteTip



Canal



Skeleton



Half-Skeleton

MICROSHELL COLORS



Beige



Light Brown



Translucent
Medium Brown



Translucent
Dark Brown



Black



Transparent

Retention Options

Skeleton	Canal Lock
Half Skeleton	Helix Lock
Semi Skeleton	

BTE STYLES

More miniBTE R
More miniBTE T

Opn S BTE
Plus Power 13

Xceed BTE
Super Power 13

Xceed BTE
Ultra Power 675

85


105

110

120



BTE Fitting Ranges



CORDA MINIFIT THIN TUBE



Standard 0.9 mm
Corda miniFit thin tube



Power 1.3 mm
Corda miniFit thin tube
Use with Opn S BTE Plus Power 13 only

RITE AND BTE INSTRUMENT COLORS

44 Silver	92 Steel Grey*	91 Silver Grey*	90 Chroma Beige	94 Terracotta	93 Chestnut Brown	63 Black
79 Baby Pink** Hear Pink***	57 Power Pink**	58 Aquamarine**	47 Blue**	46 Red**	48 Emerald Green**	45 Purple**

*Not available for Opn Play

**Available for Opn Play miniRITE T only

***Available for Oticon More miniRITE R, miniRITE T, miniBTE R and miniBTE T

MINIFIT EARPIECE COMPATIBILITY OVERVIEW OTICON MORE™ VS OTICON OPN S™

miniFit earpieces	Open Dome*	OpenBass Dome**	Bass Dome double vent	Bass Dome single vent*	Power Dome	Grip Tip open	Grip Tip closed
Oticon More (miniFit speaker)	-	60 85	60 85 100	-	60 85 100	60 85	60 85
Oticon Opn S, Opn and Opn Play (miniFit speaker)	60 85	-	60 85 100	60 85 100	60 85 100	60 85	60 85

*For Opn S products only **For Oticon More products only

ROES EARMOLD AND RECEIVER CODES

	ROES Receiver Code	Description	Left Part No.	Right Part No.
60 miniFit	REC2	60 miniFit Speaker, Size 0	124946	124945
	REC3	60 miniFit Speaker, Size 1	124595	124590
	REC4	60 miniFit Speaker, Size 2	124596	124591
	REC5	60 miniFit Speaker, Size 3	124597	124592
	REC6	60 miniFit Speaker, Size 4	124598	124593
	REC7	60 miniFit Speaker, Size 5	124599	124594
85 miniFit	REC8	85 miniFit Speaker, Size 0	124968	124961
	REC9	85 miniFit Speaker, Size 1	124969	124963
	REC10	85 miniFit Speaker, Size 2	124970	124964
	REC11	85 miniFit Speaker, Size 3	124971	124965
	REC12	85 miniFit Speaker, Size 4	124973	124966
	REC13	85 miniFit Speaker, Size 5	124974	124967
100 miniFit	REC14	100 miniFit Speaker, Size 1	124585	124579
	REC15	100 miniFit Speaker, Size 2	124586	124581
	REC16	100 miniFit Speaker, Size 3	124587	124582
	REC17	100 miniFit Speaker, Size 4	124588	124583
	REC18	100 miniFit Speaker, Size 5	124589	124584

	ROES Receiver Code	Description	Left Part No.	Right Part No.
Power Mold 100	REC37	Power Mold 100 Wire, Size 1	120212	120208
	REC38	Power Mold 100 Wire, Size 2	120213	120209
	REC39	Power Mold 100 Wire, Size 3	120214	117854
	REC40	Power Mold 100 Wire, Size 4	120215	120210
	REC41	Power Mold 100 Wire, Size 5	120216	120211
Power Mold 105	REC59	Power Mold 105 Wire, Size 1	120212	120208
	REC60	Power Mold 105 Wire, Size 2	120213	120209
	REC61	Power Mold 105 Wire, Size 3	120214	117854
	REC62	Power Mold 105 Wire, Size 4	120215	120210
	REC63	Power Mold 105 Wire, Size 5	120216	120211






Earmold Styles & Codes	
EM9	BTE Mold
EM10	RITE Mold
EM11	Slim Tube Mold
EM12	RITE Power Mold
EM14	MicroShell

Other Features

SnapLock 13 Tubing	Red & Blue dots
Standard BTE Tubing sizes	Soft Tip
Corda (Thin Tube) Tubing sizes	Removal String
L/R Initials	Removal Handle (post)
Red & Blue mold tips	

CUSTOM STYLES

Opn 1 Custom Styles



Custom Faceplate colors

C005 Black*

C004 DarkBrown

C003 Medium Brown

C002 Light Brown

C001 Beige

*IIC only

*Android devices need to support ASHA to allow direct connectivity to Oticon More. Please visit [oticon.com/support/compatibility](https://www.oticon.com/support/compatibility) for more information.

Made for
iPhone | iPad | iPod

Works with
android

Apple, the Apple logo, iPhone, iPad, and iPod touch are trademarks of Apple Inc., registered in the U.S. and other countries. App Store is a service mark of Apple Inc. Android, Google Play, and the Google Play logo are trademarks of Google LLC.

Questions? We're here to help!
Oticon Government Services Direct Line: **877.310.9681**

YOUR CHOICE OF CONVENIENT CHARGERS FOR OTICON MORE*



Oticon miniBTE SmartCharger*



Oticon miniRITE SmartCharger**

NEW

Oticon SmartChargers

- Protective Lid
- Power Bank
- Drying Capabilities

*Oticon miniBTE SmartCharger is compatible with Oticon More devices on Firmware 1.3.

**Oticon miniRITE SmartCharger is compatible with Oticon More devices on Firmware 1.1.1 or newer.

Oticon miniRITE SmartCharger is not compatible with Opn S miniRITE R.

***Features for both Oticon SmartCharger and Desktop charger.



Oticon miniBTE Charger



Desktop Chargers

- Easy insertion and removal
- Contact Free***
- Quick charge



Oticon miniRITE Charger



OTICON MORE, OPN S, OPN, OPN PLAY & XCEED WIRELESS ACCESSORIES



ConnectClip

- Make and receive hands-free calls to both Android™ and iPhone® devices
- Wirelessly stream music, videos or other sounds from your mobile phone to your hearing aids
- Control the volume and change programs remotely
- Remote microphone



TV Adapter 3.0

- Seamlessly stream audio from a TV
- Wireless streaming connection up to 45 ft
- Digital stereo (PCM)
- Digital/RCA connection
- Can pair to multiple sets of Oticon More, Opn S™, Opn™ and Xceed devices



Remote Control 3.0

- Adjust volume
- Switch programs
- Control connectivity sources



BTD800 Bluetooth Dongle

- Small USB device to improve wireless signal on computers
- Stream video calls, movies and music from your computer
- ConnectClip required



Phone Adapter 2.0

- Turn a landline phone into a wireless Bluetooth phone
- Wireless streaming connection up to 45 ft
- ConnectClip required



EduMic

- Automatic pairing
- 3.5mm input jack
- Telecoil