

3. Install and replace the battery

- (1) Remove three screws on the back of the thermometer
- (2) Open the back base and install 18650 Li battery



4. Install method



1. Heat foot



2. Bracket fixing



3. Double-sided adhesive sticking

WMT5 User Guide




How to connect and use Bluetooth

Step 1: Download the app by scanning the QR code



Step 2: After downloading the app, you will be prompted to enter your name

Step 3: Click the  icon on the top right corner to connect your device to the thermometer

Step 4: Click on  to see temperature measurement from the thermometer in real time

Step 5: Click on  to add multiple profiles

Step 6: Click on "Body temperature alarm setting" to set a number for temperature at or past which alarm will go off

Step 7: Click on  to see the trend of temperature during different hours of the day

Step 8: Click on  to see the history of the temperature data for multiple dates and times

