

medifab

medifab 7009-0001-269  
Poz In Form Positioning  
System



# medifab 7009-0001-269 Poz In Form Positioning System Owner's Manual

[Home](#) » [medifab](#) » medifab 7009-0001-269 Poz In Form Positioning System Owner's Manual 

## Contents

- 1 [medifab 7009-0001-269 Poz In Form Positioning System](#)
- 2 [Product Specifications](#)
- 3 [Product Usage Instructions](#)
- 4 [Client Details](#)
- 5 [Poz'In'Form Kit](#)
- 6 [Cylindrical](#)
- 7 [Occiput](#)
- 8 [Ring & Half Ring](#)
- 9 [Hand Wedge](#)
- 10 [Documents / Resources](#)
  - 10.1 [References](#)
- 11 [Related Posts](#)

medifab

**medifab 7009-0001-269 Poz In Form Positioning System**



## Product Specifications

- **Poz'In'Form Kit:** Pharmatex 7009-0001-269, Lenzing FR 7009-0001-279
- **Half-Moon:** Pharmatex 7009-1851-269, Lenzing FR 7009-1851-279
- **Delta:** Pharmatex 7009-0551-269, Lenzing FR 7009-0551-279
- **Semi-Lateral Decubitus:** Pharmatex 7009-3001-269, Lenzing FR 7009-3001-279
- **Cylindrical:** Pharmatex 7009-1501-269, Lenzing FR 7009-1501-279
- **Wedge Abduction:** Pharmatex 7009-1551-269, Lenzing FR 7009-1551-279
- **Occiput:** Pharmatex 7009-4751-269, Lenzing FR 7009-4751-279
- **Ring & Half Ring:** Pharmatex 7009-4081-269, Lenzing FR 7009-4081-279
- **Universal:** Pharmatex 7009-1251-269, Depth: 250mm
- **Hand Wedge:** Pharmatex 7009-5411-269, Lenzing FR 7009-5411-279
- **Heel Cushion:** Pharmatex 7009-5201-269, Lenzing FR 7009-5201-279
- **Heel Wedge:** Pharmatex 7009-5401-269, Lenzing FR 7009-5401-279

## Product Usage Instructions

### Poz'In'Form Kit

The Poz'In'Form Kit includes all necessary components for clinical lying and positioning solutions. To use the kit effectively, follow these steps:

Identify the specific component needed for the positioning requirements.

1. Place the component securely under or around the patient as directed by healthcare professionals.
2. Ensure proper alignment and support to maintain a comfortable and safe position for the patient.

## Half-Moon

The Half-Moon component is designed to provide support in a specific lying position. Follow these steps for usage:

1. Position the Half-Moon under the patient's body to provide support and reduce pressure points.
2. Adjust the placement as needed to achieve optimal comfort and alignment.

## Frequently Asked Questions (FAQ)

### Q: How do I clean the components?

A: It is recommended that the components be cleaned with a mild detergent and water. Ensure thorough drying before reuse.

## CLINICAL LYING & POSITIONING SOLUTIONS

### Client Details

- Client/User Name:
- Clients DOB:
- Date:
- Notes:
- Poz'In'Form® elements and modules are entirely removable and washable. They are delivered with a Pharmatex or Lenzing FR® cover.
- The Pharmatex bi-elastic cover is made of polyurethane and is recommended for use in healthcare facilities. Pharmatex is soft to the touch and very comfortable next to the skin, while still easy to clean.
- The Lenzing FR® cover is made of viscose, coated with a waterproof polyurethane film. Lenzing FR® is a major innovation in positioning of the long-term bedridden, preventing maceration. Hygroscopic properties of viscose allow moisture absorption which is particularly important for people with thermoregulation issues. Moreover elasticity and flexibility of Lenzing FR® are particularly suitable for problems regularly observed in classical positioning.
- Lenzing spare covers can be purchased to cover Pharmatex pillows for a more absorbent feel.



## Poz'In'Form Kit

### Pharmatex

7009-0001-269

## Lenzing FR

7009-0001-279

- Transport bag
- Universal cushion (550mm x 400mm)
- Half-moon cushion
- Semi-lateral decubitus cushion
- Cylindrical cushion
- Ring cushion
- Heel cushion with strap
- Occiput head cushion



## Half-Moon

Pharmatex 7009-1851-269		Lenzing FR 7009-1851-279	
<b>Width</b> 300mm <b>Length</b> 2000mm			
<b>Spare Cover</b> 7009-1851-269-C		<b>Spare Cover</b> 7009-1851-279-C	

## Delta

<b>Pharmatex</b> 7009-0551-269 <input type="checkbox"/>	<b>Lenzing FR</b> 7009-0551-279 <input type="checkbox"/>
Length 560mm Depth 320mm Height 280mm	
<b>Spare Cover</b> 7009-0551-269-C <input type="checkbox"/>	<b>Spare Cover</b> 7009-0551-279-C <input type="checkbox"/>

**Semi-Lateral Decubitus**

<b>Pharmatex</b> 7009-3001-269 <input type="checkbox"/>	<b>Lenzing FR</b> 7009-3001-279 <input type="checkbox"/>
Width 350mm Length 1850mm	
<b>Spare Cover</b> 7009-3001-269-C <input type="checkbox"/>	<b>Spare Cover</b> 7009-3001-279-C <input type="checkbox"/>

**Cylindrical**

<b>Pharmatex</b> 7009-1501-269 <input type="checkbox"/>	<b>Lenzing FR</b> 7009-1501-279 <input type="checkbox"/>
Length 750mm Diameter 180mm	
<b>Spare Cover</b> 7009-1501-269-C <input type="checkbox"/>	<b>Spare Cover</b> 7009-1501-279-C <input type="checkbox"/>

**Wedge Abduction**

Pharmatex 7009-1551-269	<input type="checkbox"/>	Lenzing FR 7009-1551-279	<input type="checkbox"/>
<div>  <div> <b>Length</b> 350mm  <b>Height</b> 150mm  <b>Small Width</b> 100mm  <b>Large Width</b> 300mm </div> </div>			
Spare Cover 7009-1551-269-C	<input type="checkbox"/>	Spare Cover 7009-1551-279-C	<input type="checkbox"/>

## Occiput

Pharmatex 7009-4751-269	<input type="checkbox"/>	Lenzing FR 7009-4751-279	<input type="checkbox"/>
<div> <div> <b>Length</b> 500mm  <b>Width</b> 400mm  <b>Height</b> 90mm </div>  </div>			
Spare Cover 7009-4751-269-C	<input type="checkbox"/>	Spare Cover 7009-4751-279-C	<input type="checkbox"/>

## Ring & Half Ring

Pharmatex 7009-4081-269	<input type="checkbox"/>	Lenzing FR 7009-4081-279	<input type="checkbox"/>	Pharmatex 7009-4051-269	<input type="checkbox"/>	Lenzing FR 7009-4051-279	<input type="checkbox"/>
<div> <div> <b>Ring</b>   <b>Width</b> 200mm  <b>Length</b> 1350mm </div>  <div> <b>Half Ring</b>   <b>Width</b> 220mm  <b>Length</b> 1050mm </div>  </div>							
Spare Cover 7009-4081-269-C	<input type="checkbox"/>	Spare Cover 7009-4081-279-C	<input type="checkbox"/>	Spare Cover 7009-4051-269-C	<input type="checkbox"/>	Spare Cover 7009-4051-279-C	<input type="checkbox"/>

## Universal



Pharmatex			Lenzing FR		
Code	7009-1251-269	<input type="checkbox"/>	Code	7009-1251-279	<input type="checkbox"/>
Length	350mm		Length	350mm	
Depth	250mm		Depth	250mm	
Spare cover: 7009-1251-269-C <input type="checkbox"/>			Spare cover: 7009-1251-279-C <input type="checkbox"/>		
Code	7009-1253-269	<input type="checkbox"/>	Code	7009-1253-279	<input type="checkbox"/>
Length	400mm		Length	400mm	
Depth	150mm		Depth	150mm	
Spare cover: 7009-1253-269-C <input type="checkbox"/>			Spare cover: 7009-1253-279-C <input type="checkbox"/>		
Code	7009-1252-269	<input type="checkbox"/>	Code	7009-1252-279	<input type="checkbox"/>
Length	550mm		Length	550mm	
Depth	400mm		Depth	400mm	
Spare cover: 7009-1252-269-C <input type="checkbox"/>			Spare cover: 7009-1252-279-C <input type="checkbox"/>		

## Hand Wedge

Pharmatex 7009-5411-269		Lenzing FR 7009-5411-279	
<b>Length</b> 390mm <b>Height</b> 90mm <b>Width</b> 220mm			
<b>Spare Cover</b> 7009-5411-269-C <input type="checkbox"/>		<b>Spare Cover</b> 7009-5411-279-C <input type="checkbox"/>	

## Heel Wedge

Pharmatex 7009-5401-269		Lenzing FR 7009-5401-279	
<b>Length</b> 480mm <b>Height</b> 120mm <b>Width</b> 300mm			
<b>Spare Cover</b> 7009-5401-269-C <input type="checkbox"/>		<b>Spare Cover</b> 7009-5401-279-C <input type="checkbox"/>	

## Heel Cushion


<b>Pharmatex</b> 7009-5201-269 <input type="checkbox"/>	<b>Lenzing FR</b> 7009-5201-279 <input type="checkbox"/>
<b>Length</b> 850mm <b>Diameter</b> 250mm	
<b>Spare Cover</b> 7009-5201-269-C <input type="checkbox"/>	<b>Spare Cover</b> 7009-5201-279-C <input type="checkbox"/>

[solutions@medifab.com](mailto:solutions@medifab.com)  
[medifab.com](http://medifab.com)

#### DISTRIBUTED BY:

Script: Poz 'In' Form Script Form 20240808\_EN\_v0.3

#### Documents / Resources

	<a href="#">medifab 7009-0001-269 Poz In Form Positioning System</a> [pdf] Owner's Manual 7009-0001-269, 7009-0001-269 Poz In Form Positioning System, Poz In Form Positioning System, Positioning System, System
--	--

#### References

- [m Home Page - Medifab : Medifab](#)
- [User Manual](#)

[Manuals+](#), [Privacy Policy](#)

This website is an independent publication and is neither affiliated with nor endorsed by any of the trademark owners. The "Bluetooth®" word mark and logos are registered trademarks owned by Bluetooth SIG, Inc. The "Wi-Fi®" word mark and logos are registered trademarks owned by the Wi-Fi Alliance. Any use of these marks on this website does not imply any affiliation with or endorsement.