



NICHIRYO 00-PMNEO Pipette Mate NEO Controller User Manual

[Home](#) » [NICHIRYO](#) » NICHIRYO 00-PMNEO Pipette Mate NEO Controller User Manual 

NICHIRYO

00-PMNEO Pipette Mate NEO Controller
User Manual

Contents

- [1 Pipette Controller](#)
- [2 Specifications](#)
- [3 Documents / Resources](#)
 - [3.1 References](#)
- [4 Related Posts](#)

Pipette Controller

Pipette Mate **NEO**



- User-Friendly pipette controller, 167g only.
- Long-life Lithium Ion battery installed, 4-hours charging for 1700 times.
- Dial-typed Adjustable speed controller for your choice.
- UV resistant for Clean Benchwork.
- Self-standing type without holder.
- Fully autoclavable spare parts. (Nosecone, Silicon Holder and Filter)
- Complies with CE and RoHS.



Specifications

Cat. No	Description
00-PMNEO	Pipette Mate NEO
Batteries	Lithium Ion 3.6V, 600mAh
Charger	100/240V 50/60Hz,0.2A
Filter	20mm×0.2μm Hydrophobic
Pipette Types	plastic or glass from 1 to 100mL Plastics and glass pipettes not sold by Nichiryo

Cat. No	Description
00-PM-N003	Battery cover
00-PM-N005	Head fitting
00-PM-N013	Tube holder
00-PM-N100	Filter, 0.2μm
00-PM-N024	Travel charger 5V US
00-PM-N030	Travel charger base unit
00-PM-N025	Travel charger UK adaptor*
00-PM-N026	Travel charger AU adaptor*
00-PM-N027	Travel charger EU

UK and AU adaptors need a travel charger base unit.
visit our website for part lists.
www.nichiryo.co.jp/en/
NICHIRYO PRODUCT LINEUP

Documents / Resources

	NICHIRYO 00-PMNEO Pipette Mate NEO Controller [pdf] User Manual 00-PMNEO, Pipette Mate NEO Controller, 00-PMNEO Pipette Mate NEO Controller
---	--

References

-  [PIPETTE / DISPENSER / CUSTOM MADE & OEM NICHIRYO](#)

Manuals+.