

O2 Concepts 800-0001-2 Portable Oxygen Concentrator



O2 Concepts 800-0001-2 Portable Oxygen Concentrator Instruction Manual

[Home](#) » [O2 Concepts](#) » O2 Concepts 800-0001-2 Portable Oxygen Concentrator Instruction Manual 

Contents

- 1 O2 Concepts 800-0001-2 Portable Oxygen Concentrator
- 2 Oxlife Independence
- 3 Product Specifications
- 4 O2 Concepts Proprietary DNA Technology
- 5 Documents / Resources
 - 5.1 References
- 6 Related Posts



O2 Concepts 800-0001-2 Portable Oxygen Concentrator



Oxlife Independence

- Meets a broad spectrum of LTOT patient requirements
- Continuous flow settings up to 3LPM
- Pulse flow settings from .5 to 6
- Proven 24/7 operation: year after year
- Ready for use with CPAP and BiPAP
- Dual battery operation of up to 6 hours
- AC & DC operation in all settings
- Five-year warranty
- Integrated cart + travel-ready
- Dual reimbursement codes (E1390 & E1392)

Patient-centric Care

- Best clinical choice for preventing COPD re-admittance
- Improved patient outcomes through the use of the latest technologies
- Enhanced patient lifestyle through clinical management and mobility
- Continuous flow for sleep and higher flow patient needs
- Proven platform for patient monitoring and telehealth
- EnergySmart Technology® delivers extended battery life and unmatched reliability and durability

Product Specifications

- Dimensions Height: 20.2" • Width: 10.9" • Depth: 9.45"
- Weight2 Device: 16.7 lbs • Battery: 1.35 lbs
- Oxygen Flow Rate (Pulse 15-40 BPM)
 - Continuous: 0.5 – 3 LPM in 0.5 increments
 - Pulse: 0.5 – 6.0 in 0.5 increments (8 – 96 ML)
- Battery Duration
 - 5.75 hours @ 2 pulse at 20 bp
 - 2.5 hours @ 2 LPM continuous
- Battery Recharge Time4 2.5 hours
- Oxygen Concentration 87% – 95% at sea level at all settings
- AC Power Universal Input:
 - 100 – 240 VAC
 - 50 – 60 Hz • 2.0A maximum
- DC Power 12 – 15 VDC • 15 amps

1. O2 Concepts® reserves the right to revise specifications at anytime without notice.
2. Production weight excludes cart kit.
3. Battery duration presented for two-battery system.
4. Battery charge time increases if device is operating while charging.

O2 Concepts Proprietary DNA Technology

Delivering the only Self-Contained Device-Initiated Telehealth Solution

- Establish the therapy and utilize the full capabilities of DNA Technology® to manage the device and patient use remotely.
- Utilize the GPS location feature and Geofence capabilities to manage warehouse, offsite hospital and DME closet inventory.
- Remotely monitor equipment to ensure proper function and operation – purity, flow, alerts and alarms: all from your office or clinic.
- Troubleshoot customer service calls remotely by utilizing the device report and history as well as other tools in the portal.
- Predict or preempt maintenance issues and direct technicians in the most efficient manner.
- View hours used per day, settings, indications of mobility, etc.
- Simplified and improved compliance reporting.

Highly Efficient Hospital Discharge Program

Discharge your patients on an Independence they can take directly from the hospital. The clinician or respiratory therapist can provide patient orientation with the therapy and device prior to discharge; including instruction and support to family members and care givers while familiarizing the patient with the equipment.

Eliminate Tanks, Stationaries or Home Fill Deliveries

Eliminate in-home appointments and setups with a no-contact Independence delivery, relieving both your drivers/technicians and patients of in-person contact, extended wait times and after hour and weekend setups.

Designed for Self-Guided Operation

Online instructional videos, Quick Start Guides and an intuitive control panel provide easy-to-use device operation and useful tools for follow-up customer support.



American Medical Sales and Rentals
15055 E. Hinsdale Drive, Bldg D
Centennial, CO 80112

Phone: +1 [877-774-9271](tel:877-774-9271)





www.oxygenconcentratorstore.com

E-mail: sales@amsrco.com

Documents / Resources

	O2 Concepts 800-0001-2 Portable Oxygen Concentrator [pdf] Instruction Manual 800-0001-2 Portable Oxygen Concentrator, 800-0001-2, Portable Oxygen Concentrator, Oxygen Concentrator, Concentrator
	O2 Concepts 800-0001-2 Portable Oxygen Concentrator [pdf] User Manual 800-0001-2 Portable Oxygen Concentrator, 800-0001-2, Portable Oxygen Concentrator, Oxygen Concentrator, Concentrator

References

-  [O2 Concepts](#)
-  [O2 Concepts](#)
-  [O2 Concepts](#)
-  [Oxygen Concentrator Store: Oxygen Suppliers for Oxygen Equipment](#)
- [User Manual](#)

Manuals+, Privacy Policy

This website is an independent publication and is neither affiliated with nor endorsed by any of the trademark owners. The "Bluetooth®" word mark and logos are registered trademarks owned by Bluetooth SIG, Inc. The "Wi-Fi®" word mark and logos are registered trademarks owned by the Wi-Fi Alliance. Any use of these marks on this website does not imply any affiliation with or endorsement.