


## omnipod<sup>®</sup> 5 App for iPhone



## omnipod 5 App for iPhone User Guide


[Home](#) » [omnipod](#) » omnipod 5 App for iPhone User Guide 

### Contents

- [1 omnipod 5 App for iPhone](#)
- [2 Product Information](#)
- [3 Product Usage Instructions](#)
- [4 FAQ](#)
- [5 Introduction](#)
- [6 To Update from the Notification](#)
- [7 CONTACT](#)
- [8 Documents / Resources](#)
  - [8.1 References](#)
- [9 Related Posts](#)



**omnipod 5 App for iPhone**



## Product Information

### Specifications

- **Product Name:** Omnipod 5
- **Compatibility:** iPhone
- **App Store:** TestFlight version available, official version to be released

### Product Usage Instructions

- When the official Omnipod 5 App for iPhone is released, follow these steps to update:
- When you receive the notification, tap on Update Now.
- The App Store will open to the Omnipod 5 App, tap Update to proceed with the update.
- Open the App Store on your iPhone.
- Search for “Omnipod 5” and select the app that shows the option to update.
- Tap on Update to install the latest version of the app.

### FAQ

- **Q:** What should I do if I accidentally delete the TestFlight-downloaded App?
- **A:** If you delete the TestFlight App before updating to the new version from the App Store, you will lose your settings and adaptivity. In this case, you will need to complete the first-time setup again.

## Introduction


- When the Omnipod 5 App for iPhone is officially released, you will need to update your App from the TestFlight version now available in the App Store.
- Your settings and adaptivity will be transferred to the official version of the App.
- On your current Omnipod 5 App, you will get a notification telling you to update your App.

**Caution:** Do NOT delete your TestFlight-downloaded App until you have updated to the new App from the App Store. If you delete the TestFlight App before you update to the new App, you will lose your saved settings and adaptivity and will have to complete the first-time setup again.


## To Update from the Notification

It is recommended to update the App when you get the notification. If you tap Not Now, you will get a notification every 72 hours.

- Tap Update Now.



- The App Store will open to the Omnipod 5 App. Tap Update.



## To Update from the App Store


- On your iPhone, open the App Store.




- Search for Omnipod 5.



- Select the Omnipod 5 App that shows you have the option to update.



- Tap Update.




CONTACT

- All rights reserved. All other trademarks are the property of their respective owners.
- The use of third-party trademarks does not constitute an endorsement or imply a relationship or affiliation.
- Patent information at [insulet.com/patents](https://insulet.com/patents) INS-OHS-09-2024-00104 V1.0

© 2024 Insulet Corporation. Insulet, Omnipod, the Omnipod logo, and Simplify Life, are trademarks or registered trademarks of Insulet Corporation.

Documents / Resources



[omnipod 5 App for iPhone](#) [pdf] User Guide  
5 App for iPhone, App for iPhone, for iPhone, iPhone

References

- [Patents | Insulet](#)
- [User Manual](#)

[Manuals+](#), [Privacy Policy](#)

This website is an independent publication and is neither affiliated with nor endorsed by any of the trademark owners. The "Bluetooth®" word mark and logos are registered trademarks owned by Bluetooth SIG, Inc. The "Wi-Fi®" word mark and logos are registered trademarks owned by the Wi-Fi Alliance. Any use of these marks on this website does not imply any affiliation with or endorsement.