



ottobock 28U11 WalkOn Flex and Flex Junior Instruction Manual

[Home](#) » [ottobock](#) » ottobock 28U11 WalkOn Flex and Flex Junior Instruction Manual 

ottobock 28U11 WalkOn Flex and Flex Junior



Contents

[1 Instructions for use](#)







[2 Customer Support](#)






[3 Documents /
Resources](#)

[3.1 References](#)

[4 Related Posts](#)

Instructions for use





28U11* Walk On	Article number	Side	Shoe size EU	Calf height	
	28U11=L36–39	left	36–39	35 cm	13.8 in
	28U11=R36–39	right			
	28U11=L39–42	left	39–42	36 cm	14.2 in
	28U11=R39–42	right			
	28U11=L42–45	left	42–45	37,5 cm	14.8 in
	28U11=R42–45	right			
	28U11=L45–48	left	45–48	39 cm	15.4 in
	28U11=R45–48	right			



28U22* Walk On Flex	Article number	Side	Shoe size EU	Calf height	
	28U22=L36–39	left	36–39	32,5 cm	12.8 in
	28U22=R36–39	right			
	28U22=L39–42	left	39–42	35 cm	13.8 in
	28U22=R39–42	right			
	28U22=L42–45	left	42–45	37 cm	14.6 in
	28U22=R42–45	right			
	28U22=L45–48	left	45–48	39 cm	15.4 in
	28U22=R45–48	right			



28U22* Walk On Flex Junior	Article number	Side	Shoe size EU	Calf height	
	28U22=L24–27	left	24–27	20,5 cm	8.0 in
	28U22=R24–27	right			
	28U22=L27–30	left	27–30	24 cm	9.5 in
	28U22=R27–30	right			
	28U22=L30–33	left	30–33	27 cm	10.6 in
	28U22=R30–33	right			
	28U22=L33–36	left	33–36	29,5 cm	11.6 in
	28U22=R33–36	right			




28U23* Walk On Trimmable	Article number	Side	Shoe size EU	Calf height	
	28U23=L36–39	left	36–39	35 cm	13.8 in
	28U23=R36–39	right			
	28U23=L39–42	left	39–42	36 cm	14.2 in
	28U23=R39–42	right			
	28U23=L42–45	left	42–45	37,5 cm	14.8 in
	28U23=R42–45	right			
	28U23=L45–48	left	45–48	39 cm	15.4 in
	28U23=R45–48	right			



28U33* Walk On Lateral	Article number	Side	Shoe size EU	Calf height	
	28U33=L36–39	left	36–39	32,5 cm	12.8 in
	28U33=R36–39	right			
	28U33=L39–42	left	39–42	35 cm	13.8 in
	28U33=R39–42	right			
	28U33=L42–45	left	42–45	37 cm	14.6 in
	28U33=R42–45	right			
	28U33=L45–48	left	45–48	39 cm	15.4 in
	28U33=R45–48	right			

28U33* Walk On Lateral Junior	Article number	Side	Shoe size EU	Calf height	
	28U33=L24–27	left	24–27	20,5 cm	8.0 in
	28U33=R24–27	right			
	28U33=L27–30	left	27–30	24 cm	9.5 in
	28U33=R27–30	right			
	28U33=L30–33	left	30–33	27 cm	10.6 in
	28U33=R30–33	right			
	28U33=L33–36	left	33–36	29,5 cm	11.6 in
	28U33=R33–36	right			

Material	28U11, 28U23, 28U33				
	Structural part: carbon fiber composite material; Calf strap and padding: Polyurethane				
	28U22				
	Structural part: glass fiber composite material; Calf strap and padding: Polyurethane				
	Latex-free: This product does not contain natural rubber.				



Customer Support



Ottobock SE & Co. KGaA
Max -Näder -Straße 15
· 37115 Duderstadt Germany
T +49 5527 848-0
F +49 5527 848-3360
healthcare@ottobock.de
www.ottobock.com

ottobock.

Documents / Resources



[ottobock 28U11 WalkOn Flex and Flex Junior](#) [pdf] Instruction Manual
28U11 WalkOn Flex and Flex Junior, 28U11, WalkOn Flex and Flex Junior, Flex and Flex Junior, Flex Junior, Junior

References

- [Welcome to Ottobock](#)