

PHONAK
cShell 6.0 Build
Style Hearing Aid



PHONAK cShell 6.0 Build Style Hearing Aid Owner's Manual

[Home](#) » [PHONAK](#) » PHONAK cShell 6.0 Build Style Hearing Aid Owner's Manual 

Contents

- [1 PHONAK cShell 6.0 Build Style Hearing Aid](#)
- [2 Customer information](#)
- [3 Hearing instrument](#)
- [4 Accessories and Roger™](#)
- [5 Preferences](#)
- [6 Patient information](#)
- [7 HI warranty](#)
- [8 Special Instructions](#)
- [9 Additional Information](#)
- [10 Infinio hearing instrument color options](#)
- [11 Documents / Resources](#)
 - [11.1 References](#)
- [12 Related Posts](#)

PHONAK

PHONAK cShell 6.0 Build Style Hearing Aid



Customer information

Ship to account number:

Address:

City:

State:

Zip:

Bill to account number:

Third party patient number:

Fitting date:

Purchase order number:

Contact name:

Phone:

Email:

Hearing instrument

Qty	Instrument type	Technology level			
		90	70	50	30
	Phonak Audéo™ I Sphere	<input type="checkbox"/>	<input type="checkbox"/>	-	-
	Phonak Audéo I-R	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	Phonak CROS I-R	-	-	-	-

Specify color choice:

cShell style

- cShell Acrylic (062-0056-01)
- Shell colors: Black (06), Cocoa (22), Pink (26), Tan (14), Brown (28),
- Red (R), Blue (B), Transparent (21)
- Faceplate colors: Black (06), Cocoa (22), Pink (26), Tan (14), Brown (28)



☐ L ☐ R

- cShell Titanium (062-0057)
- Shell color: Titanium gray (U0)
- Faceplate color: Black (06)



☐ L ☐ R

Accessories and Roger™

	TV Connector (076-3006-0611)	<input type="checkbox"/>	
Accessories	PartnerMic (076-4001-P511)	<input type="checkbox"/>	
	RemoteControl (076-0065-P5)	<input type="checkbox"/>	
Roger	Roger On™ 3 (056-3032-xx011) Champagne (P5), Graphite Gray (P7)	<input type="checkbox"/> P5	<input type="checkbox"/> P7
	Roger Select™ 3 (056-3033-xx011) Champagne (P5), Graphite Gray (P7)	<input type="checkbox"/> P5	<input type="checkbox"/> P7
	Roger Table Mic (056-4012-P5011)	<input type="checkbox"/>	

Preferences

- If necessary, may we change the following: ☐ Please call

Yes

☐☐☐

No

☐☐☐

- UP to P OR P to M OR M to S if audiogram permits
- Build larger if required
- I would like to receive a paper user guide

Not available for cShell Titanium ** HF4 for UP receiver, CeruStop for all other receivers

Patient information

Last name:

First name:

Age:

Gender:

- Is this patient order being fully or partially paid for by state or federal government funds? ☐ Yes ☐ No

Audiogram (Required for AOV)

	Hz	250	500	1K	2K	3K	4K
Left:	AC						
Right:	AC						

HI warranty

- ☐ 4th year ☐ 5th year (Patient age required at time of order)

Shipping and handling options†

- ☐ 2 Business Days (\$20.99)
- ☐ 3–5 Business Days (\$19.99)
- ☐ C.O.D. (\$21.99 + shipping)

† Based on location and courier

Options

Shell Color	Standard color (Transparent 21)	<input checked="" type="checkbox"/>	
	Custom color	Code:	<input type="text"/>
Faceplate Color	Standard color	<input checked="" type="checkbox"/>	
	Custom color	Code:	<input type="text"/>
Receiver	S Receiver	L	<input type="checkbox"/> R <input type="checkbox"/>
	M Receiver (standard)	L	<input checked="" type="checkbox"/> R <input checked="" type="checkbox"/>
	P Receiver	L	<input type="checkbox"/> R <input type="checkbox"/>
	UP Receiver (may result in cShell XL)	L	<input type="checkbox"/> R <input type="checkbox"/>
Length	0/1/2/3 Size 2 is standard	Left:	<input type="text"/>
		Right:	<input type="text"/>
Vent size	AOV (enter audiogram)	<input checked="" type="checkbox"/>	
	Customer specific 12/20/25/3L/none	Left:	<input type="text"/> mm
		Right:	<input type="text"/> mm
Wax system	Standard wax system**	<input checked="" type="checkbox"/>	
	Custom wax system	<input type="text"/>	
Retention*	Helix Lock (HL)	<input type="checkbox"/>	
	Canal Lock (CL)	<input type="checkbox"/>	
	Skeleton Lock (SL)	<input type="checkbox"/>	
Removal line (RF)	Transparent (21)	<input type="checkbox"/>	
	Line extension + customer specific	<input type="text"/>	
Titanium surface	Extra retention (ERS)	<input type="checkbox"/>	

Special Instructions

Legend ☐ / ☒ = Standard ☐ = Optional

Additional Information

Alternative cShell Build Style

- We offer an alternative cShell build style with an anterior tube outlet.
- This build style provides high cosmetic acceptance in spite of a slightly larger size and can be used in case of build problems due to anatomical limitations. This build style is available for acrylic shells only.



Fig.1 cShell - standard tube



Fig.2 cShell - anterior tube outlet

Acoustically Optimized Vent (AOV)

We recommend selecting AOV. We will set the correct vent diameter and style consider audiogram/shell geometry/feedback canceler and required low-frequency gain. The client's audiogram is mandatory for the calculation of the AOV.

Extra Retention Surface (ERS)

- We offer a rougher surface for cShell titanium. The extra retention surface can improve retention in the ear.

Wax systems


cShell Acrylic	Description
X	No wax guard
CB	CeruStop
H4	HF4 Filter, default with UP receivers

Color codes

Code	Color
06	Black
13/21	Transparent
14	Tan
22	Cocoa
26	Pink
28	Brown
R/B	Red/Blue
U0	Titanium Gray

Receiver tube length

- Measure the tube length to select the correct receiver. Measure both ears for binaural fitting.



Infinio hearing instrument color options



P1

Sand Beige



P5

Champagne



P6

Silver Gray



P7

Graphite Gray

- All of our products, including custom products and spare parts, can be ordered online in the Phonak store.

- Sonova USA Inc. is not responsible and assumes no liability for any device or accessory not manufactured by Sonova that is sent by you to Sonova.
- Please ensure that you only include Phonak-branded devices and accessories herein.
- Products, options, and accessories are subject to change without notice.

Last updated: December 2, 2024 2024-12/hc © 2024 Sonova AG. All rights reserved.

Documents / Resources



[PHONAK cShell 6.0 Build Style Hearing Aid](#) [pdf] Owner's Manual
062-0056-01, 062-0057, 076-3006-0611, 076-4001-P511, 076-0065-P5, 056-3032-xx011, 056-3033-xx011, 056-4012-P5011, cShell 6.0 Build Style Hearing Aid, cShell 6.0, Build Style Hearing Aid, Style Hearing Aid, Hearing Aid, Aid

References

- [User Manual](#)

[Manuals+](#), [Privacy Policy](#)

This website is an independent publication and is neither affiliated with nor endorsed by any of the trademark owners. The "Bluetooth®" word mark and logos are registered trademarks owned by Bluetooth SIG, Inc. The "Wi-Fi®" word mark and logos are registered trademarks owned by the Wi-Fi Alliance. Any use of these marks on this website does not imply any affiliation with or endorsement.