

PHONAK

PHONAK Multitool
tamperproof battery
drawer tool Connevans



PHONAK Multitool tamperproof battery drawer tool Connevans User Manual

[Home](#) » [PHONAK](#) » PHONAK Multitool tamperproof battery drawer tool Connevans User Manual 

Contents

- [1 PHONAK Multitool tamperproof battery drawer tool Connevans](#)
- [2 Product Information](#)
- [3 Product Usage Instructions](#)
- [4 Multitool for hearing aids](#)
- [5 Documents / Resources](#)
 - [5.1 References](#)
- [6 Related Posts](#)

PHONAK

PHONAK Multitool tamperproof battery drawer tool Connevans



Product Information

Specifications:

- Product: Multitool for Hearing Aids
- Includes: Cleaning brush with integrated magnet system
- Usage: Wax removal and cleaning of earmolds and hearing systems
- Precautions: External use only, keep out of reach of children

Product Usage Instructions

1. Ensure the Multitool set is clean and dry before use.
2. To remove wax from earmolds, gently use the cleaning brush included in the set.
3. For cleaning hearing systems, use the brush with the integrated magnet system for easy battery removal.
4. After each use, clean the tools properly and store them in a safe place.

Precautions:

- Use the Multitool set externally only.
- Keep the set out of reach of children to prevent any accidents.

FAQ:

1. Can I use the Multitool set for other purposes?

No, the Multitool set is specifically designed for wax removal and cleaning of earmolds and hearing systems. Do not use it for any other purposes.

2. How often should I clean the Multitool set?

It is recommended to clean the tools after each use to maintain hygiene and effectiveness.


3. What should I do in case of problems with my hearing aid?

If you experience any issues with your hearing aid, contact a hearing aid specialist for assistance.


User manual

Multitool for hearing aids

- The Multitool set includes necessary tools for wax removal and cleaning of earmolds and hearing systems.
- It contains a cleaning brush with an integrated magnet system for easy battery removal.



- Precautions for use: External use. Keep out of reach of children.
- If you experience problems, contact a hearing aid specialist. All serious incidents related to the device should be the subject of a notification to the manufacturer and the competent authorities in the Member State in which the user is located.



Documents / Resources

This website is an independent publication and is neither affiliated with nor endorsed by any of the trademark owners. The "Bluetooth®" word mark and logos are registered trademarks owned by Bluetooth SIG, Inc. The "Wi-Fi®" word mark and logos are registered trademarks owned by the Wi-Fi Alliance. Any use of these marks on this website does not imply any affiliation with or endorsement.