



ResMed AirView Data Management Device
Compatibility List



ResMed AirView Data Management Device Compatibility List User Guide


[Home](#) » [ResMed](#) » ResMed AirView Data Management Device Compatibility List User Guide 

Contents

- [1 ResMed AirView Data Management Device Compatibility List](#)
- [2 Product Usage Instructions:](#)
- [3 Data management/Device compatibility list](#)
- [4 Documents / Resources](#)
 - [4.1 References](#)
- [5 Related Posts](#)



ResMed AirView Data Management Device Compatibility List



Specifications

- Brand: ResMed
- Product Name: AirMini™
- Compatible Devices: AirCurve™ CS-A, AirCurve CS-A TJ, AirSense™ 10 Respond A, AirSense 10 Respond W, Sleepmate 10 Auto,
- Sleepmate 10 Elite, AirSense 10 Auto, AirSense 11 Auto, AirSense 11 Respond, AirSense 11 AutoSet, AirSense 11 Elite, AirCurve 10 VAuto, AirCurve 10 CS-A, AirCurve 10 CS-A TJ, S9™ Respond, Sleepmate S9 Auto, Sleepmate S9 Elite,
- RESMED S9 Auto, AutoSet C, S8™ Respond, Sleepmate S8, NIP Nasal III, VPAP Adapt SV, AutoSet CS, AutoSet CS-A, VPAP Adapt SV-A, NIP Nasal V, Clean Air VELIA, Clean Air VS ULTRA, Clean Air VS INTEGRA, Astral 100, Astral 150,
- CleanAir Astral, Stellar 100, Stellar 130, Stellar 150, Lumis 100 VPAP S, Lumis 100 VPAP ST, Lumis 150 VPAP ST, Lumis 150 VPAP ST-A, Lumis HFT, Lumis HFT

Product Usage Instructions:

Setup

1. Unbox the ResMed AirMini™ device.
2. Connect the appropriate compatible device to the AirMini™ device following the user guide.

Data Management

1. Ensure the devices are correctly connected for data management.
2. Refer to the relevant User Guides for detailed instructions on data management.

Troubleshooting

1. If you encounter any issues with data transfer or device compatibility, refer to the troubleshooting section of the User Guides
2. Contact ResMed customer support if problems persist.

FAQ

Q: How do I know if my device is compatible with ResMed data management products?

A: Refer to the Device Compatibility List provided in the user manual to check if your device is compatible.

Q: Where can I find more information about ResMed products and locations?

A: Visit [ResMed.com](https://www.resmed.com) for information on products and worldwide locations of ResMed offices.

Data management/Device compatibility list

This supplement specifies the devices compatible with ResMed data management products. For full details regarding the correct use of these devices, see the relevant User Guides.

	ResScan™ 7.1	EasyCare Tx / Tx Link™ 9. 0	EasyCare Tx 2 / RCM RCMH** TxLink 2™ 1.0	AirView™ TxLink 2™ 1.0
RESMED AirMini™				✓
RESMED AirCurve™ CS-A / RESMED AirCurve CS-A TJ	✓	✓		
AirSense™ 10 Respond A / AirSense 10 Respond W / Sleepmate 10 Auto / Sleepmate 10 Elite / AirSense 10 Auto*	✓			✓
AirSense 11 Auto / AirSense 11 Respond / AirSense 11 AutoSet / AirSense 11 Elite	✓	✓		✓
AirCurve 10 VAuto / AirCurve 10 CS-A / AirCurve 10 CS-A TJ	✓	✓		✓
S9™ Respond / Sleepmate S9 Auto / Sleepmate S9 Elite / RESMED S9 Auto†	✓	✓		✓
AutoSet C / S8™ Respond	✓	✓		✓
Sleepmate S8	✓	✓		✓
NIP Nasal III	✓	✓		✓
VPAP Adapt SV / AutoSet CS	✓	✓		✓

AutoSet CS-A / VPAP Ad apt SV-A	✓	✓	✓
NIP Nasal V / Clean Air V ELIA	✓	✓	✓
Clean Air VS ULTRA™ / Cl ean Air VS INTEGRA™	✓	✓	✓
Astral™ 100 / Astral 150 / CleanAir Astral	✓	✓	✓
Stellar™ 100 / Stellar 130 / Stellar 150	✓	✓	✓
Lumis™ 100 VPAP S / Lu mis 100 VPAP ST / Lumis 150 VPAP ST / Lumis 150 VPAP ST-A / Lumis HFT	✓	✓	✓
Lumis HFT	✓	✓	✓

* Compatible with AirView via wireless connectivity only.

**For a list of available data for use by the EMR system, please contact your ResMed representative. ‡ Compatible with AirView via SD card uploads only.

RCM – ResMed Connectivity Module

RCMH – ResMed Connectivity Module Hospital

Notes:

- Not all products available in all regions.
- AirView compatibility may not be available in all regions.

ResMed Pty Ltd

1 Elizabeth Macarthur Drive Bella Vista NSW 2153 Australia

ResMed Corp

9001 Spectrum Center Boulevard San Diego CA 92123 USA

See [ResMed.com](https://www.resmed.com) for other ResMed locations worldwide. AirMini, AirCurve, AirSense, Tx Link, TxLink 2, Lumis, S10 Auto, S9 Auto, Clean Air ULTRA, Clean Air INTEGRA, ResScan, AirView and EasyCare Tx are trademarks and/or registered trademarks of the ResMed family of companies. For patent and other intellectual property information, see [ResMed.com/ip](https://www.resmed.com/ip). 10110573/8 2024-01

[ResMed.com](https://www.resmed.com)

Documents / Resources

	ResMed AirView Data Management Device Compatibility List [pdf] User Guide AirView Data Management Device Compatibility List, AirView, Data Management Device Compatibility List, Management Device Compatibility List, Device Compatibility List, Compatibility List, List
--	---

References

- [User Manual](#)

Manuals+, Privacy Policy

This website is an independent publication and is neither affiliated with nor endorsed by any of the trademark owners. The "Bluetooth®" word mark and logos are registered trademarks owned by Bluetooth SIG, Inc. The "Wi-Fi®" word mark and logos are registered trademarks owned by the Wi-Fi Alliance. Any use of these marks on this website does not imply any affiliation with or endorsement.