



TPP 99722 Syringe Filter Instruction Manual

[Home](#) » [TPP](#) » TPP 99722 Syringe Filter Instruction Manual 




**99722 Syringe Filter
Instruction Manual**

Contents

- [1 Instructions for Use in Pictures TPP Syringe Filter](#)
- [2 Documents / Resources](#)
 - [2.1 References](#)
- [3 Related Posts](#)

Instructions for Use in Pictures TPP Syringe Filter




Do not re-sterilize or reuse the TPP syringe filter. For single use only.

Measurements	99722	99745
Inlet to outlet mm	27	27
Diameter mm	33	33
Filter size cm ²	4.52	4.52
Pore size µm	0.22	0.45
Max. using temp °C	45	45
Flow rate water 25 °C / 2 bar	≥ 150 mL/min	≥ 180 mL/min
Burst pressure at 25 °C	10 bar	10 bar
Filtration volume mL	10 – 200	10 – 200
Dead volume following further flushing with air mL	0.1	0.1
Connection In-let	Luer-lock	Luer-lock
Connection out-let	Luer-conus	Luer-conus

Source: TPP

Note: TPP products are for research use only and are not intended for use in clinical, diagnostic or therapeutic procedures. TPP assumes no responsibility for damages or defects resulting from unapproved or improper use. **TPP** recommendations are guidelines and do not represent absolute values. TPP assumes no responsibility for the currency, accuracy, completeness, or quality of this TechDoc.

Documents / Resources

	TPP 99722 Syringe Filter [pdf] Instruction Manual 99722 Syringe Filter, 99722, Syringe Filter, Filter
---	---

References

- [User Manual](#)