

trackit Ticketing System Software User Guide

Home » Trackit » trackit Ticketing System Software User Guide 🖫




Contents

- 1 trackit Ticketing System Software
- **2 Product Information**
- 3 Usage Instructions
- **4 New Features in Trackit Plus** v2.8.4.7
- **5 New Events Logged**
- 6 Documents / Resources
 - **6.1 References**



trackit Ticketing System Software



Specifications

Product Name: Trackit Plus

• Version: 2.8.4.7

Manufacturer: Lifelines Neuro
Release Date: May 2023

Product Information

Trackit Plus v2.8.4.7 is a cutting-edge monitoring system designed to automatically record and analyze various events during a study. It features new functionalities such as tracking battery percentage at the start of a study, monitoring lock screen events, and tracking review activities.

Usage Instructions

1. Battery Percentage Tracking:

At the beginning of each study, the system reads the amplifier's battery percentage and records it as an event.

2. Lock Screen Events:

The system logs events for entering and exiting the lock screen. The time spent in the lock screen is calculated between these events.

3. Review Activities:


Events are recorded for entering and exiting the review section. The system tracks the time spent in the review phase.

FAQ

Q: How can I get additional assistance?

A: For additional assistance, you can contact our customer support team at Tel: 866-889-6505 or Email: support@lifelinesneuro.com

New Features in Trackit Plus v2.8.4.7



New Features in Trackit Plus v2.8.4.7 - New Events Logged

New Events Logged

There are several automatic events recorded with a study. These include:

· Amp Battery % at start of study.

The actual % capacity of the battery is read from the amplifier and added as an event at the start of the recording.



• Event for lock screen enter/Event for lock screen exit.

The period between the Enter and Exit Lockscreen events is the time spent in the Lockscreen and the next Enter event gives the time out of the Lockscreen.



• Event for Review enter/Event for Review exit.

The period between the Enter and Exit Review events is the time spent in Review.



For additional assistance

Tel: 866-889-6505

Email: support@lifelinesneuro.com

Documents / Resources



trackit Ticketing System Software [pdf] User Guide v2.8.4.7, Ticketing System Software, System Software, Software

References

• User Manual

Manuals+, Privacy Policy

This website is an independent publication and is neither affiliated with nor endorsed by any of the trademark owners. The "Bluetooth®" word mark and logos are registered trademarks owned by Bluetooth SIG, Inc. The "Wi-Fi®" word mark and logos are registered trademarks owned by the Wi-Fi Alliance. Any use of these marks on this website does not imply any affiliation with or endorsement.