

Waters

ACQUITY UPLC FLEXCART

Contents [[hide](#)]

[1 Product BULLETIN](#)

[2 Specifications](#)

[3 Documents / Resources](#)

[3.1 References](#)

Product BULLETIN

The Waters ® ACQUITY UPLC FlexCart wheeled cart is a complete system platform that facilitates convenient installation of the ACQUITY UPLC and makes it easier to position the system close to the inlet of the mass spectrometer. The FlexCart provides electrical outlets, supports computer monitor and keypad, and manages fluid waste container. TM

The ACQUITY UPLC FlexCart has maximum accessibility from the sides for workstation and waste-can interface. AC power is available at the rear of the cart for distribution to the requisite ACQUITY UPLC modules and peripherals.

The FlexCart's mobile platform allows the ACQUITY UPLC system to be moved among several workplaces. Its adjustable tabletop height allows optimum positioning of the column outlet of the system with a mass spectrometer.



Specifications

Dimensions deep	30 inches (76.1 cm) wide x 33 inches (83.8 cm)
Minimum Table Height	29.1 inches (73.9 cm)
Maximum Table Height	41.5 inches (105.4 cm)
Height Adjustment	12.3 inches (31.2 cm)
Cart Weight	170 lbs (77.3 kg)
Supportable Weight	400 lbs (181.8 kg)

Order Information

ACQUITY UPLC FlexCart

Part Number

205015015


For Complete Confidence

Waters and ACQUITY UPLC are trademarks of Waters Corporation.

All other trademarks are property of their respective owners.

©2005 Waters Corporation Produced in the USA February 2005 720001089EN MC-PDF

Documents / Resources

	Waters 205015015 Acquity Uplc Flex Cart [pdf] Instructions 205015015 Acquity Uplc Flex Cart, 205015015, Acquity Uplc Flex Cart, Uplc Flex Cart, Flex Cart
---	--

References

- [User Manual](#)

205015015, 205015015 Acquity Uplc Flex Cart, Acquity Uplc Flex Cart, Flex Cart, Uplc Flex Cart,
Waters Waters

Leave a comment

Your email address will not be published. Required fields are marked *

Comment *

Name

Email

Website

☐ Save my name, email, and website in this browser for the next time I comment.

Post Comment

Search:

e.g. whirlpool wrf535swhz

Search

[Manuals+](#) | [Upload](#) | [Deep Search](#) | [Privacy Policy](#) | [@manuals.plus](#) | [YouTube](#)

This website is an independent publication and is neither affiliated with nor endorsed by any of the trademark owners. The "Bluetooth®" word mark and logos are registered trademarks owned by Bluetooth SIG, Inc. The "Wi-Fi®" word mark and logos are registered trademarks owned by the Wi-Fi Alliance. Any use of these marks on this website does not imply any affiliation with or endorsement.