



WITHINGS Body Scan Scale Installation Guide

[Home](#) » [Withings](#) » WITHINGS Body Scan Scale Installation Guide 

WITHINGS Body Scan Scale Installation Guide




THANK YOU FOR CHOOSING THE WITHINGS BODY

Contents


- 1 [Box contents](#)
- 2 [Quick start](#)
- 3 [Know more with Body](#)
- 4 [Customize your screens](#)
- 5 [Overview](#)
- 6 [Important notice](#)
- 7 [Warranty](#)
- 8 [Need help?](#)
- 9 [Documents / Resources](#)
- 10 [Related Posts](#)

Box contents


Withings Body




Carpet feet



Alkaline cells already installed




Installation of carpet feet (optional)




Quick start

1. Pull the tab underneath the scale




2. Visit go.withings.com to get the app and follow the instructions.



Know more with Body

Analyze your Body Mass Index and be notified of Rain and Weather Forecast for the day, while weighing.



*Only with Wi-Fi connection

- Weight (kg, lb, st lb)



- BMI



- Number of steps*




- Weather* (morning | afternoon, °C or °F)




Customize your screens


Change the order of the screens and choose what matters most to you.

If you use a Withings watch, activity tracker or track your steps with the Health Mate™ app, activate the Steps screen to know how much you walked the previous day!



Overview





- (A) Screen
(B) Feet
(C) Unit & Pairing button
(D) Battery Cover

Important notice



This product should not be used to diagnose or treat any medical condition.

Information to users in Canada

IC: 11411A-WBS06

This device complies with Industry Canada's licence exempt RSSs. Operation is subject to the following two conditions:

1. This device may not cause interference; and
2. This device must accept any interference, including interference that may cause undesired operation of the device.

Radiation Exposure Statement

The product complies with the Canada portable RF exposure limit set forth for an uncontrolled environment and are safe for intended operation as described in this manual. The further RF exposure reduction can be achieved if the product can be kept as far as possible from the user body or set the device to lower output power if such function is available.

Warranty

Withings™ Two (2) Years Limited Warranty – Withings Body Withings, 2 rue Maurice Hartmann, 92130 Issy-les-Moulineaux, ("Withings") warrants the Withings-branded hardware product ("Withings Product") against defects in materials and workmanship when used normally in accordance with Withings' published guidelines for a period of TWO (2) YEARS from the date of original retail purchase by the end-user purchaser («Warranty Period»).

Withings' published guidelines include but are not limited to information contained in technical specifications, safety instructions or quick start guide. Withings does not warrant that the operation of the Withings Product will be uninterrupted or error-free.

Withings is not responsible for damage arising from failure to follow instructions relating to the Withings Product's use

Need help?

This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions:

1. This device may not cause harmful interference, and
2. This device must accept any interference received, including interference that may cause undesired operation.

The crossed out wheeled bin symbol on your product, battery and/or accessories means that the item is classified as electrical and electronic equipment. Such items should not be mixed with general household waste, and must be taken to dedicated collection points at the end of their working life for proper treatment, recovery and recycling. This way you help save valuable resources and promote recycling of materials.

Mode	Frequency Band (MHz)	Maximum Output Power (dBm)
WLAN	2412-2472	9.87
BT EDR	2400-2483.5	0.85
BT LE	2400-2483.5	9.87

Withings Body | BMI Wi-Fi Scale
© 2021 Withings. All rights reserved.

Documents / Resources

	WITHINGS Body Scan Scale [pdf] Installation Guide Body Scan Scale
-------------------------------------------------------------------------------------	--------------------------------------------------------------------------------------