



YEEZOO Two Steps Bed Stool with Handrails Installation Guide

[Home](#) » [YEEZOO](#) » YEEZOO Two Steps Bed Stool with Handrails Installation Guide 

YEEZOO Two Steps Bed Stool with Handrails



Contents


- 1 ASSEMBLY INSTRUCTION
- 2 INSTALL THE OTHERS TWO CROSS BARS
- 3 INSTALL THE UPPER WOODEN PANEL
- 4 INSTALL THE FIRST WOODEN PANEL
- 5 PREPARE TWO GRAB BARS
- 6 INSTALL HANDRAILS AND FIX THEM
- 7 INSTALL THE FOUR NON -SLIP PADS
- 8 Documents / Resources
 - 8.1 References
- 9 Related Posts

ASSEMBLY INSTRUCTION






2*Wooden Panels



3*Cross Bars



1* Grab Bar



1* Grab Bar

**2* Grab Bar
Connections**



8*M7*13mm



8*M6*49mm



4*Non-slip Pads




1*Open-end Wrench



2*Sizes Hex Wrench

INSTALL THE OTHERS TWO CROSS BARS

- ▶ Handrail holders must face outwards.
- ▶ Use two M7*13mm screws and hex wrench.
- ▶ The cross bars should be level.



Tip: If you find support brackets not level, please adjust the cross bars.

INSTALL THE UPPER WOODEN PANEL

- ▶ The panel with logo "YEEZOO" is the upper panel, logo outwards.
- ▶ Four countersinks and wood surface texture face upwards.
- ▶ Use four M6*49mm screws and hex wrench to fix the panel.




INSTALL THE FIRST WOODEN PANEL

- ▶ Four countersinks and wood surface texture face upwards.
- ▶ Use four M6*49mm screws and hex wrench to fix the panel.
- ▶ Check and tighten all screws.



PREPARE TWO GRAB BARS

- ▶ The higher one: bar connection and grab bar are screw thread connected.
- ▶ The shorter one: No need to assemble.



INSTALL HANDRAILS AND FIX THEM

- ▶ Insert the handrails into the handrail holders.
- ▶ Use four M7*13mm screws and hex wrench to fix from it's bottom.
- ▶ The two grab bars do not distinguish between left and right handrail holders.

05

Fix Handrail



INSTALL THE FOUR NON -SLIP PADS

- ▶ Warm Tip1: Please double check to make sure all screws have been tightened before stepping on.
- ▶ Warm Tip2: Please don't shake the handrail or jump on the steps to avoid dangers.



Documents / Resources



[YEEZOO Two Steps Bed Stool with Handrails](#) [pdf] Installation Guide

Two Steps Bed Stool with Handrails, Steps Bed Stool with Handrails, Bed Stool with Handrails, Stool with Handrails, Handrails

References

- [User Manual](#)

[Manuals+](#), [Privacy Policy](#)

This website is an independent publication and is neither affiliated with nor endorsed by any of the trademark owners. The "Bluetooth®" word mark and logos are registered trademarks owned by Bluetooth SIG, Inc. The "Wi-Fi®" word mark and logos are registered trademarks owned by the Wi-Fi Alliance. Any use of these marks on this website does not imply any affiliation with or endorsement.