

[manuals.plus](#) /

› [Tempur-Pedic](#) /

› [Tempur-Pedic TEMPUR-ProAdapt 12-Inch Medium Foam Mattress, Twin XL Instruction Manual](#)

Tempur-Pedic TEMPUR-ProAdapt

Tempur-Pedic TEMPUR-ProAdapt 12-Inch Medium Foam Mattress, Twin XL Instruction Manual

Model: TEMPUR-ProAdapt | Size: Twin XL | Firmness: Medium

1. INTRODUCTION TO YOUR TEMPUR-PROADAPT MATTRESS

The Tempur-Pedic TEMPUR-ProAdapt 12-Inch Medium Foam Mattress is engineered to provide advanced pressure relief for restorative rest. This mattress features Tempur-Pedic's most pressure-relieving material ever, designed to offer rejuvenation and relaxation. Its dual cover design provides a cool-to-touch feel and features a zip-off construction for easy maintenance. The TEMPUR material is known for retaining its shape, ensuring consistent support and alignment over time.



Image: The Tempur-Pedic TEMPUR-ProAdapt 12-Inch Medium Foam Mattress, Twin XL, showcasing its design and profile.

2. PRODUCT FEATURES

- **12-Inch TEMPUR-PEDIC MATTRESS:** Developed by the original memory foam mattress brand, backed by decades of innovation.
- **UP TO 2X PRESSURE RELIEF:** Features Tempur-Pedic's most pressure-relieving material for deeper, more rejuvenating sleep.
- **DUAL COVER CONSTRUCTION:** Next-generation technology with a zip-off, cool-to-touch outer layer and super-stretch inner layer for convenience and comfort.
- **LEGENDARY COMFORT & SUPPORT:** NASA-engineered TEMPUR material is perfected to deliver unique relaxation and comfort.
- **DURABLE & RELIABLE:** Adaptive TEMPUR material never loses its shape, offering superior support and alignment year after year.



 TEMPUR-PEDIC®

Image: Diagram illustrating the layers of innovation within the Tempur-Pedic mattress, including SmartClimate Dual Cover System, Original TEMPUR Comfort Layer, and TEMPUR-APR Support Layer.

Your browser does not support the video tag.

Video: An official Tempur-Pedic video showcasing the PRO Adapt mattress features and technology, including its layers and cooling properties.

3. SETUP AND INSTALLATION

Your Tempur-Pedic TEMPUR-ProAdapt mattress is designed for easy setup. Upon delivery, professional delivery experts will handle the heavy lifting, including bringing the mattress into your home and installing it. They can also haul away your old mattress, providing a hassle-free experience.

- Delivery and Placement:** Your mattress will be delivered and placed in your desired room by trained professionals.
- Old Mattress Removal:** Old mattress removal service is typically included to simplify your setup.
- Break-in Period:** Like new footwear, your Tempur-Pedic mattress requires a break-in period, typically around 30 days, to fully adapt and conform to your unique body shape and sleeping habits. During this time, the material will soften and adjust to provide optimal comfort and support.

Your browser does not support the video tag.

Video: An official Tempur-Pedic video detailing the in-home delivery and setup process, including removal of old mattresses and careful placement of the new one.

4. USAGE AND COMPATIBILITY

The TEMPUR-ProAdapt mattress is designed to adapt to your body's weight, shape, and temperature, providing personalized comfort and support. The medium firmness is ideal for sleepers who need a balance of comfort and support, offering pressure relief while helping maintain alignment. It is also suitable for those who switch positions during the night.

Base Compatibility

For optimal performance and to ensure the longevity of your Tempur-Pedic mattress, it is crucial to use a compatible base. Tempur-Pedic mattresses require a solid, elevated surface.

- **Compatible Bases:** Power Base, Flat Foundation, Platform Bed.
- **Incompatible Bases:** Box Spring, The Floor.



Image: A visual guide showing compatible and incompatible bases for Tempur-Pedic mattresses, including Power Base, Flat Foundation, Platform Bed, Box Spring, and The Floor.

5. MAINTENANCE AND CARE

Proper care ensures the longevity and performance of your Tempur-Pedic mattress.

- **Cover Care:** The dual cover features a zip-off, cool-to-touch outer layer that is machine-washable. Remove the cover and follow the washing instructions on the label.
- **Spot Cleaning:** For spills or stains on the mattress core, spot clean with a mild detergent and water. Do not saturate the mattress. Allow to air dry completely before replacing the cover.
- **Rotation:** Regularly rotating your mattress (head to foot) can help promote even wear and extend its lifespan.




Image: An icon representing easy care, indicating the machine-washable cover of the Tempur-Pedic mattress.

6. SPECIFICATIONS

Feature	Detail
---------	--------

Product Dimensions	80"L x 38"W x 12.2"Th
Item Weight	65 Pounds
Item Firmness Description	Medium
Top Style	Plush Top
Cover Material	Polyester Blend
Number of Layers	3
Fill Material	Memory Foam
Special Feature	Cooling
Assembly Required	No
Made In	USA



Image: An icon indicating the mattress is assembled in the USA.

7. WARRANTY INFORMATION

This limited warranty is valid only in connection with the original purchase of new products from authorized retailers and extends from the original purchase date of the original product purchased. An “original purchaser,” for the purposes of this limited warranty, is an individual or entity who purchases the product directly from Tempur-Pedic or an authorized retailer of Tempur-Pedic with the intent to use the product for personal consumer use and not for commercial or industrial use and not with the intent to resell the product. An “authorized retailer,” for purposes of this limited warranty, is an individual or entity authorized by Tempur-Pedic to sell product directly to original purchasers. An individual or entity that purchases the product from whatever source with the intent to resell the product is an unauthorized reseller (“unauthorized reseller”).

8. CUSTOMER SUPPORT

For any questions or assistance with your Tempur-Pedic mattress, please contact customer support. You may also refer to the return policy for details on returns.

- **Return Policy:** 90-day easy returns.

© 2026 Tempur-Pedic. All rights reserved.