

Solaris Laboratories NY Level UP Microcurrent Conductive Gel

Solaris Laboratories NY Level UP Microcurrent Conductive Gel Instruction Manual

Model: Level UP Microcurrent Conductive Gel

1. PRODUCT OVERVIEW

The Solaris Laboratories NY Level UP Microcurrent Conductive Gel is a moisturizing gel specifically formulated to enhance the effectiveness of microcurrent and radio frequency device treatments. It ensures optimal conductivity, providing a non-sticky, water-based formula that delivers long-lasting hydration for up to 24 hours. This gel is designed to facilitate smooth device movement across the skin, optimizing the benefits of your facial treatments.

2. KEY FEATURES

- **Conductive:** Optimizes the performance of microcurrent and radio frequency devices.
- **Lightweight:** Non-sticky and comfortable on the skin.
- **Moisturizing:** Provides long-lasting hydration for up to 24 hours.
- **Cruelty-Free:** Developed without animal testing.
- **Deeply Hydrating:** Enhances skin moisture levels.
- **For All Skin Types:** Suitable for various skin sensitivities and conditions.

3. HOW TO USE

1. **Prepare Skin:** Ensure your face is clean and dry before application.
2. **Apply Gel:** Apply a generous layer of the Level UP Microcurrent Conductive Gel to the areas of your face and neck where you intend to use your microcurrent device. Ensure even coverage for optimal conductivity.
3. **Perform Treatment:** Following the instructions specific to your microcurrent device, gently massage the treated areas. For a typical facial treatment, spend approximately 5 minutes, dedicating 2.5 minutes to each side of the face. The neck and chin areas can also be massaged.
4. **Post-Treatment:** The gel is designed to be hydrating and non-sticky, allowing it to be left on the skin after treatment for continued moisturization. If desired, any excess gel can be gently wiped off.



Image: Solaris Laboratories NY Level UP Microcurrent Conductive Gel shown with other skincare items and a microcurrent device, illustrating its use in a complete skincare routine.

4. INGREDIENTS

The Level UP Microcurrent Conductive Gel is formulated with key ingredients to support skin health and device performance. While a full list is available on the product packaging, notable ingredients include:

- **Hydrolyzed Collagen:** Supports skin elasticity and firmness.
- **Vitamin C:** Known for its antioxidant properties and ability to brighten skin.
- **Niacinamide:** Helps improve skin tone, reduce redness, and strengthen the skin barrier.
- **Hyaluronic Acid:** Provides intense hydration, helping the skin retain moisture.

5. IMPORTANT INFORMATION

Legal Disclaimer: Statements regarding dietary supplements have not been evaluated by the FDA and are not intended to diagnose, treat, cure, or prevent any disease or health condition.

6. SPECIFICATIONS

Attribute	Detail
Brand	Solaris Laboratories NY
Product Name	Level UP Microcurrent Conductive Gel
ASIN	B0CLM4RKTR
UPC	040835718336
Product Benefits	Hydrating
Manufacturer	Solaris Laboratories NY

Date First Available

October 23, 2023

7. MEDIA GALLERY



Solaris Laboratories NY Level UP Microcurrent Conductive Gel product packaging.



Two women smiling, with the SLNY logo visible, representing the brand's focus on diverse skin types.



A woman holding a pink microcurrent facial device, demonstrating its use.



A woman using a scalp massager, part of Solaris Laboratories NY's broader product range.