

Feonase NEW682

Feonase Queen Size Metal Bed Frame User Manual

Model: NEW682 | Brand: Feonase

1. IMPORTANT SAFETY INFORMATION

Always assemble the bed frame according to the provided instructions to ensure stability and prevent collapse. Be cautious of sharp edges and pinch points during assembly and use. Do not exceed the maximum weight capacity of 1000 lbs to maintain the bed frame's integrity and safety.


For added safety, this bed frame includes fluorescent strips on the bed legs to help locate them in low light, preventing accidental toe stubbing. Wrap-in corner protectors are also provided to protect children from accidental impacts.

Embedded Design

Embedded bed corners easy to place mattresses and prevent mattresses from sliding, no box spring is needed.



Fluorescent Strips: Located on bed legs for visibility in low light.



FEONASE SLATS (RECTANGULAR)

- Larger force area offers 1000 lbs weight capacity
- Muted buckles completely eliminate noise
- Stable without deformation

OTHER SLATS (ELLIPTICAL)

- Small force area, low weight capacity
- Easy to squeak
- Easily deformed by pressure



Corner Protectors: Soft covers for bed frame corners to prevent injury.

2. PRODUCT OVERVIEW

The Feonase Queen Size Metal Bed Frame is a heavy-duty platform bed designed for durability, comfort, and ease of assembly. It features a linen upholstered headboard and a robust steel slat support system, eliminating the need for a box spring.

Key Features:

- **Heavy-Duty Construction:** Consists of hardened square steel slats and 9 thick support legs, capable of bearing up to 1000 lbs.
- **Linen Upholstered Headboard:** Features skin-friendly linen button tufting with high-density foam for a soft and comfortable experience.
- **Noise-Free Design:** Built-in anti-noise protection shells on the sides of the metal slats and EVA mute pads prevent friction, ensuring a silent sleeping environment.
- **Ample Under-Bed Storage:** Provides 12 inches of clearance, offering significant space for storage and helping to keep your room tidy.
- **Embedded Mattress Design:** Keeps the mattress from sagging or slipping, providing stable support without the need for a box spring.



Product View: The Feonase Queen Size Metal Bed Frame in Grey with its upholstered headboard.




Large Under-Bed Storage

Keep your room tidy and clean at any time

Headboard Detail: Close-up of the skin-friendly linen button tufting.



Strong Central Foundation



Zero Noise Buckles



9 Thick Legs Supports

Under-Bed Storage: Illustration of the 12-inch clearance for ample storage.




Embedded Design

Prevent mattresses from shifting

Protective Feet Plugs

Prevent scratching the floor and slipping



Noise-Free EVA Sponge

Isolate the friction between slats and frame, no more squeak

Structural Features: Details on the strong central foundation, noise-reducing buckles, and nine support legs.

3. COMPONENTS LIST

Before beginning assembly, please ensure all components listed below are present and undamaged. If any parts are missing or damaged, please contact customer support.

- Metal Bed Frame x 1
- Linen Upholstered Headboard x 1
- Heavy-Duty Slats x 12
- Installation Tools & Accessories
- Installation Manual (this document)

4. SETUP AND ASSEMBLY

Assembly of the Feonase bed frame is designed to be straightforward and can typically be completed within 30 minutes. All necessary tools and accessories are included in the package.

Assembly Steps:

1. **Unpack and Identify Parts:** Carefully unpack all components and lay them out on a clean, soft surface to prevent scratches. Refer to the components list to ensure all parts are present.
2. **Assemble the Frame:** Follow the diagrams in the included installation manual to connect the main metal frame sections. Ensure all connections are secure.
3. **Install Slats:** The heavy-duty slats feature a one-click buckle installation. Simply snap each slat into place on the frame. Ensure they are evenly spaced and fully engaged for optimal support.
4. **Attach Headboard:** Securely attach the linen upholstered headboard to the designated points on the bed frame using the provided hardware.
5. **Add Safety Features:** Apply the fluorescent strips to the bed legs as indicated in the manual. Install the wrap-in corner protectors on the frame corners.
6. **Final Check:** Before placing your mattress, double-check all connections and ensure the frame is stable and level.



Quick Installation: Slats snap into place for fast assembly.

NOISE - FREE

The EVA mute pads prevent friction between the slats and the frame, effectively eliminating noise.

01

/ 03

Strong Foundation: All-metal frame with central cross frame and evenly-spaced metal slats.

5. OPERATING INSTRUCTIONS

Once assembled, the Feonase bed frame is ready for use. Simply place your queen-size mattress directly onto the steel slats. No box spring is required due to the robust support system and embedded design that prevents mattress shifting.

6. MAINTENANCE

To ensure the longevity and appearance of your Feonase bed frame, follow these simple maintenance guidelines:

- **Cleaning the Frame:** Wipe the metal frame with a soft, damp cloth. Avoid harsh chemicals or abrasive cleaners that could damage the powder-coated finish.
- **Cleaning the Headboard:** For the linen upholstered headboard, gently spot clean with a mild fabric cleaner if necessary. Test on an inconspicuous area first. Do not saturate the fabric.
- **Regular Checks:** Periodically check all bolts and connections to ensure they remain tight. Re-tighten if any looseness is detected to maintain stability and prevent noise.
- **Floor Protection:** The bed frame is equipped with protective feet plugs to prevent scratching your floor. Ensure these remain in place.

7. TROUBLESHOOTING

If you encounter any issues with your Feonase bed frame, please refer to the common solutions below. If the problem persists, contact customer support.

Common Issues and Solutions:

- **Squeaking Noises:**
The Feonase bed frame is designed to be completely noise-free. If you hear squeaking, check that all slats are fully engaged in their buckles and that there are no gaps at the connections. Ensure all frame bolts are tightened securely. The built-in EVA mute pads should isolate friction; verify they are correctly positioned.
- **Mattress Shifting:**
The embedded design of the bed frame is intended to prevent mattress movement. If your mattress is shifting, ensure it is correctly centered within the frame's edges. Verify that the frame is level on your floor.
- **Instability/Wobbling:**
This indicates loose connections. Re-tighten all bolts on the frame, including those connecting the headboard and the support legs. Ensure all 9 support legs are firmly on the ground and not bent.



Noise Reduction: EVA sponge pads prevent metal-on-metal friction.

A 3D diagram of the bed frame showing the underside. A silver metal trunk is placed under the bed. A green double-headed arrow indicates the height of the clearance, which is labeled as 12". The bed frame has dark-colored legs and a slatted base.

Ultra-Large Space

12" ample under-bed height offers a large space for storing extra sundries. Save more space and no more sanitary dead ends.

Slat Design: Feonase rectangular slats offer superior support and noise elimination compared to elliptical designs.


8. SPECIFICATIONS

Detailed specifications for the Feonase Queen Size Metal Bed Frame:

Brand	Feonase
Model Number	NEW682
Size	Queen
Product Dimensions	82"L x 63"W x 40"H
Color	Grey
Finish Type	Powder Coated
Material	Heavy-Duty Metal Frame, Linen Upholstered Headboard
Weight Capacity	Up to 1000 lbs

Under-Bed Space	12 inches
Assembly Required	Yes
Compatible Mattress Size	Queen (No Box Spring Needed)

PRODUCT DIMENSION



Material	Full Size	Weight Capacity
Fabric & Metal	82x63x40 Inch	1000 lbs

Product Dimensions: Visual representation of the bed frame's measurements.

9. WARRANTY AND SUPPORT

Feonase is committed to providing products and services of impeccable quality. For any inquiries or assistance, our customer support team is available 24/7.

Protection Plans:

Optional protection plans are available for extended coverage:

- **2-Year Protection Plan:** Available for \$22.99 (one-time payment).

- **3-Year Protection Plan:** Available for \$27.99 (one-time payment).
- **Complete Protect:** One plan covering all eligible past and future purchases for \$16.99/month.

For further details on warranty coverage or to purchase a protection plan, please refer to your purchase documentation or contact Feonase customer service.

Customer Support:

For any questions, technical support, or assistance with missing/damaged parts, please visit the official Feonase store on Amazon or contact their customer service directly. Refer to your purchase receipt for specific contact information.

© 2024 Feonase. All rights reserved.

ASK A QUESTION ABOUT THIS MANUAL

Ask about setup, troubleshooting, compatibility, parts, safety, or missing instructions. Manuals+ will review the question and use this page's manual context to help answer it.

Name

Anonymous

Email

you@example.com

After a successful submission, we securely hash the supplied account or form email before sending it to Microsoft Advertising for conversion attribution. [Privacy details](#).

Question

Example: How do I reset this device or fix this error?

Details

Model number, symptoms, what you tried, or the section of the manual you are using.

Submit question

Open the live page

<https://manuals.plus/asin/B0D1YDBXT4>

More documents and user manuals: <https://manuals.plus/>