

SWZEC DOG MAT

SWZEC PawLab Cooling Mat Original Instruction Manual

Model: DOG MAT | Brand: SWZEC

INTRODUCTION

The SWZEC PawLab Cooling Mat Original is designed to provide a comfortable and cool resting surface for your pets, including dogs and cats. Constructed from advanced ice silk fabric, this mat offers a self-cooling effect, making it an ideal solution for pets during warmer temperatures. Its portable and washable design ensures convenience and hygiene.

This manual provides essential information on the proper use, care, and maintenance of your cooling mat to ensure its longevity and your pet's comfort.

FEATURES

- **Self-Cooling Technology:** Utilizes ice silk fabric to naturally dissipate heat, providing a cool surface without refrigeration or electricity.
- **Breathable Material:** Made from 75% Nylon and 25% Polyester, promoting air circulation for enhanced comfort.
- **Washable Design:** Easy to clean, ensuring a hygienic environment for your pet.
- **Portable and Lightweight:** Can be easily folded and transported, suitable for various indoor and outdoor settings.
- **Durable Construction:** Designed for regular use and pet interaction.

ICYNAP BENEFITS

ONE COOLING MAT = REAL-WORD BENEFITS

CONTINUOUS COOLING SURFACE

Releases body heat instead of trapping it

DOESN'T TRAP BODY HEAT

Surface stays comfortable while your dog rests

ENCOURAGES SETTLED RESTING

Less shifting or searching for cooler spots

CHEMICAL-FREE COOLING

No gels or toxic cooling compounds

LIGHTWEIGHT & PORTABLE

Easy to move between rooms, crates, or travel

EASILY MACHINE WASHABLE

Machine washable at 30–40°C for easy cleaning



Image: Key benefits of the cooling mat, highlighting its continuous cooling, chemical-free nature, portability, and washability.

SETUP

1. **Unpack:** Carefully remove the cooling mat from its packaging.
2. **Unfold:** Lay the mat flat in your desired location. Ensure there are no creases or folds that could affect its cooling performance or your pet's comfort.
3. **Placement:** Place the mat in an area where your pet frequently rests or where they might need extra cooling, such as a pet bed, crate, or a shaded spot outdoors.

MADE FROM COOL, DURABLE
SWEAT ABSORBING MATERIAL



COOLS YOUR DOG AS
SOON AS THEY LAY ON IT!

Image: A puppy resting on the cooling mat, illustrating its immediate cooling effect.

OPERATING INSTRUCTIONS

The SWZEC PawLab Cooling Mat operates passively, requiring no external power or refrigeration. Its cooling effect is activated by your pet's body heat.

- **Direct Contact:** For optimal cooling, ensure your pet makes direct contact with the mat's surface. The ice silk fabric draws heat away from the pet's body through conduction.
- **Heat Release:** The mat disperses the absorbed heat into the surrounding air, allowing for continuous cooling over time.
- **Environment:** While effective on its own, placing the mat in a cooler, shaded environment will enhance its cooling performance.

HOW ICYNAP ICE-SILK TECHNOLOGY WORKS



CONTINUOUS HEAT CONDUCTION

YOUR DOG'S BODY HEAT IS DRAWN INTO THE MAT THROUGH DIRECT CONTACT WITH THE COOLING SURFACE.



CONTINUOUS HEAT RELEASE

THE MAT DISPERSES ABSORBED HEAT ACROSS ITS SURFACE AND INTO THE AIR SO COOLING CAN CONTINUE OVER TIME.



NON-TOXIC MATERIALS

COOLING IS CREATED THROUGH FABRIC STRUCTURE AND MATERIAL DESIGN WITHOUT CHEMICAL GELS.

A CONTACT-BASED COOLING SURFACE DESIGNED TO HELP MOVE BODY HEAT AWAY FROM YOUR DOG.

Image: Diagram illustrating the continuous heat conduction and release mechanism of the ice silk technology.



How it works?

SELF COOLING, ABSORBS YOUR DOG'S BODY HEAT, AND PREVENTS OVERHEATING & DEHYDRATION



Image: Visual representation of how the mat absorbs body heat and prevents overheating.

CARE AND MAINTENANCE

Proper care ensures the longevity and effectiveness of your cooling mat.

- **Washing:** The mat is machine washable. Wash at 30-40°C (86-104°F) on a gentle cycle.
- **Drying:** Air dry the mat. Do not tumble dry, iron, or bleach, as this may damage the ice silk fabric and affect its cooling properties.
- **Spot Cleaning:** For minor spills or dirt, spot clean with a damp cloth and mild detergent.
- **Storage:** When not in use, fold the mat and store it in a cool, dry place away from direct sunlight.

TROUBLESHOOTING

Problem	Possible Cause	Solution
---------	----------------	----------

Problem	Possible Cause	Solution
Mat does not feel cool.	<p>No direct contact with pet.</p> <p>Ambient temperature is too high.</p> <p>Mat needs time to dissipate absorbed heat.</p>	<p>Ensure pet is lying directly on the mat.</p> <p>Move mat to a cooler, shaded area.</p> <p>Allow mat to rest for a few minutes to release heat.</p>
Pet avoids the mat.	<p>Unfamiliarity with new item.</p> <p>Mat placement is not ideal.</p>	<p>Introduce the mat gradually, perhaps with a favorite toy or treat.</p> <p>Place the mat in a location where your pet already feels comfortable.</p>
Mat appears damaged after washing.	<p>Incorrect washing or drying method.</p>	<p>Refer to the "Care and Maintenance" section for proper washing and drying instructions.</p>

SPECIFICATIONS

- **Product Name:** Summer Cooling Pet Mat
- **Model Number:** DOG MAT
- **Brand:** SWZEC
- **Fabric:** Cooling "Ice Silk" Fabric (75% Nylon, 25% Polyester)
- **Dimensions (L-70X55CM variant):** Approximately 19"L x 0.1"W x 28"Th (70cm x 55cm)
- **Weight:** Approximately 1 Pound
- **Available Colors:** Light Pink, Light Blue, Grey, Brown
- **Special Features:** Breathable, Lightweight, Washable
- **Suitable For:** Cats, Dogs, and Various Pets



Image: The SWZEC PawLab Cooling Mat Original in Blue, L-70X55CM.

WARRANTY AND SUPPORT

This SWZEC PawLab Cooling Mat comes with a **1-Year Warranty** from the date of purchase, covering manufacturing defects.

For warranty claims, technical support, or any inquiries regarding your product, please contact SWZEC customer service through the retailer where the product was purchased or visit the official SWZEC website for contact information.

Please retain your proof of purchase for warranty validation.