


PROPAGATION KIT



What's in

- High propagation dome
- 10L Orgasoilux
- 9cm PC Pots (30)
- Glue Plug Tray 48cav
- Rooting Gel 90ml
- Biodyne 100ml
- Disposable Scalpel
- Precision Trimming Scissor



*Start popping
today!*



HOW TO USE THE PROPAGATION KIT

1. First, you will need seeds. You can use about 25ml Biodyne in 1L of water. Early introduction of microbial life will strengthen your plant from a young age.
2. Transplant your seed directly into the jiffy glue plug with the pointy side facing up. Neatly close up the hole in the jiffy glue plug
3. Move your tray into your propagation dome for increased humidity which help along the germination process
4. Once your seed has sprouted in the jiffy glue plug, you can then transplant the plugs into the 9cm pots with the seedling mix. Before transpalnting soak your seedling mix with Biodyne, use 1ml per litre.
5. After 24hrs of being in the dome, move you rpots out into the fresh air
6. Use either t5 or t8 clone lights to provide just enough light for root development . Use a 18/6hr light cycle during this young age
7. As your plants grow, you can transplant them into bigger pots to increase root development and plant growth
8. Check out our clone kit for a simple cloning method.

