

INNOVATION
MADE BY GEUDER®

JACOBOVITZ

MULTIFUNCTIONAL CHOPPER

FOR SYMMETRIC BIMANUAL
PHACO FRAGMENTATION



Sérgio Jacobovitz, MD
Belo Horizonte, Minas Gerais, Brazil

Geuder®
Precision made in Germany

ONE CHOPPER - 3 FUNCTIONS

FOR MANIPULATING, DIVIDING AND
CRACKING MATURE CATARACTS



Nucleus Chopper
sharp edge for horizontal and vertical nucleus chopping techniques






Retractor
for mechanical retraction of iris and capsular rim



Nucleus Cracker
cleaver shaped tip

BIMANUAL TECHNIQUE

The two choppers are introduced through micro incisions separated by approximately 180°. The choppers retract the capsular rim so they can catch the equator and are brought together centrally to create a bisecting crack in the nucleus, dividing it cleanly in halves and quadrants. Further pieces can be divided by using it as secondary instrument in combination with the phaco tip.



G-31509
JACOBOVITZ MULTIFUNCTIONAL CHOPPER

Retractor, chopper and nucleus cracker for
symmetric bimanual phacofragmentation techniques

VIDEO

Watch the video and see how
the technique works.



GEUDER AG reserves the right to make changes to technical details
in response to recent developments. GEUDER does not assume
liability for the accuracy of each individual statement.

Illustrations not drawn to scale.



Geuder®
Precision made in Germany

GEUDER AG
Hertzstrasse 4
69126 Heidelberg
Germany

Phone: +49 6221 3066
Fax: +49 6221 303122
info@geuder.de
www.geuder.com