






ITEMS #: 5256829, 5256822
MODELS #: NZ-1010BN, NZ-1001

Fixed Post Strainer Basket

DIRECTIONS:

Place the basket in the strainer and firmly push down until the center post securely locks in place.



PROJECT SOURCE and logo design are trademarks or registered trademarks of LF, LLC. All rights reserved.

Questions, problems or missing parts? Before returning, contact us on: **866-439-9800**, 8 a.m. - 8 p.m., EST, Monday - Sunday or **ascs@lowes.com**.

05/2025