Medline MDS86031RCEH

Medline 31" Reacher Grabber Tool Instruction Manual

Model: MDS86031RCEH

Introduction

The Medline 31" Reacher Grabber Tool is designed to assist individuals with limited mobility or difficulty bending, enabling them to retrieve objects with ease. This lightweight tool provides an extended reach, making daily tasks more manageable. It features a full handgrip for secure handling and slip-resistant rubber tabs on its jaws for a firm hold on various items.





Image 1: The Medline 31" Reacher Grabber Tool. This image displays the full length of the grabber tool, highlighting its handle, shaft, and jaw mechanism.

Setup

The Medline 31" Reacher Grabber Tool requires no assembly. It is ready for immediate use directly out of its packaging.

Operating Instructions

1. **Grip the Handle:** Firmly grasp the full handgrip of the reacher tool. The design allows for the use of all four fingers for enhanced control.
2. **Activate the Jaws:** Squeeze the handle to close the composite jaws located at the opposite end of the tool.
3. **Position and Secure:** Carefully position the closed jaws around the item you wish to retrieve. The slip-resistant rubber tabs on the jaws are designed to provide a secure hold on various surfaces.
4. **Lift the Item:** Once the item is securely gripped, gently lift and move it to your desired location. The tool is designed to lift items up to 5 pounds (approximately 2.27 kg).
5. **Release the Item:** Release your grip on the handle to open the jaws and deposit the item.

This tool is suitable for a variety of tasks, including picking up objects from the floor, retrieving items from high shelves, or reaching into tight spaces.

Maintenance

- **Cleaning:** To clean the Medline Reacher Grabber Tool, wipe it with a damp cloth. Avoid using harsh or abrasive cleaning agents, as these may damage the materials.
- **Storage:** Store the tool in a dry environment, away from direct sunlight and extreme temperatures, to ensure its longevity and proper function.
- **Inspection:** Periodically inspect the jaws and handle for any signs of wear or damage. If significant damage is observed, discontinue use.

Troubleshooting

- **Issue:** The grabber jaws do not securely hold items.
  Solution: Ensure the item being picked up is within the stated 5 lb weight capacity. Verify that the jaws are fully closed around the item. Clean any debris or residue from the slip-resistant rubber tabs on the jaws, as this can affect grip.
- **Issue:** The handle mechanism feels stiff or difficult to squeeze.
  Solution: Regular use may help to loosen the mechanism over time. Ensure no foreign objects are obstructing the handle or jaw mechanism. If stiffness persists and affects usability, consider contacting customer support.

Specifications

Model Number	MDS86031RCEH
Brand	Medline
Product Dimensions	36.5 x 5.5 x 2 inches (approximately 92.7 x 14 x 5 cm)
Item Weight	3.53 ounces (approximately 100 grams)
Weight Capacity	Up to 5 lbs (approximately 2.27 kg)
Manufacturer	Medline

Warranty Information

No specific warranty information is provided with this product in the available documentation. For details regarding potential warranties, please refer to the product packaging or contact Medline customer service directly.

Support

For further assistance, product inquiries, or to report any issues with your Medline Reacher Grabber Tool, please contact Medline customer service. You can typically find contact information on the official Medline website (www.medline.com) or through the retailer from whom the product was purchased.