Saved from manuals.plus.
Open the live page for links to other documents and user manuals.

Velpeau VP0201-M

VELPEAU Soft Cervical Collar Instruction Manual

Model: VP0201-M

Introduction

Thank you for choosing the VELPEAU Soft Cervical Collar. This manual provides essential information for the proper use, care, and maintenance of your neck brace. Please read these instructions carefully before use to ensure optimal comfort and support.

Product Overview and Features

The VELPEAU Soft Cervical Collar is designed to provide support and relieve discomfort in the neck area. It is suitable for various applications, including support during sleep, reducing neck fatigue, and aiding in post-operative rehabilitation.

- **Ergonomic Design:** Features a "V" shape under the jaw for a comfortable and secure fit, with a curved rise at the back for excellent neck support.
- **Adjustable Support:** Hook and loop closure allows for customizable tightness, preventing unwanted movement.
- **Comfortable Materials:** Made with skin-friendly, breathable, and latex-free materials for continuous comfort, even during extended wear. Rapidly absorbs sweat to keep the neck dry and odor-free.
- **Versatile Versions:** Available in Comfort, Dual-use, and Enhanced versions to suit different support needs. The Comfort version prioritizes softness, the Dual-use version includes a removable external flexible plastic support plate, and the Enhanced version uses double-layer foam for powerful support.
- **Lightweight:** Designed to provide pain relief without limiting mobility.



Image: A woman wearing the blue VELPEAU Soft Cervical Collar, demonstrating its appearance when worn.



Image: A close-up view of the VELPEAU neck brace material, highlighting its soft and high-density rebound foam for comfort and support.

Sizing and Fit

To ensure the best fit, accurately measure your neck circumference and neck height using a soft measuring tape. Refer to the size chart below. If your measurements fall between two sizes, it is recommended to choose the larger size for optimal comfort and support.



Image: A diagram illustrating the correct method for measuring neck circumference and neck height to determine the appropriate brace size.

Size Chart

Size	Neck Circumference (in)	Neck Height (in)
Small	10.5" - 11.5"	2.75"
Medium	11.5" - 14"	3"
Large	14" - 16.5"	3.3"
X-Large	16.5" - 18.7"	4"

Size	Fit neck circumference	Fit neck height
S	10.5"-11.5"	2.75"
M	11.5"-14"	3"
L	14"-16.5"	3.3"
XL	16.5"-18.9"	4"




Image: A visual representation of the size chart, detailing neck circumference and height for each size option.

Setup and Application

To apply the VELPEAU Soft Cervical Collar:

1. Ensure you have selected the correct size based on the sizing chart.
2. Position the collar around your neck with the "V" shaped front section under your jaw.
3. Wrap the collar around your neck and secure it using the hook and loop fasteners at the back.
4. Adjust the tightness so the brace fits snugly but comfortably, providing support without restricting breathing or causing discomfort. The goal is to limit mobility and support the neck.
5. For Dual-use versions, you may insert or remove the external flexible plastic support plate as needed for enhanced support.



Image: A visual representation of the V-shaped ergonomic design and the optimal 20-degree support angle provided by the collar.

Operating and Usage Guidelines

The VELPEAU Soft Cervical Collar is intended to support the jaw and neck, limiting mobility to alleviate pain and pressure in the spine. It helps realign and stabilize vertebrae.

Indications for Use:

- Relief of chronic neck pain, pinched nerves, cervical sprains, stiffness, and discomfort.
- Cervical muscle relaxation and improved blood circulation.
- Support for anti-snoring and sleep apnea.
- Posture correction and prevention of forward head posture.
- Post-operative rehabilitation and injury support.

Recommended Usage:

The soft foam collar can be worn throughout the day and while sleeping for superior support and pain relief. It is particularly useful during periods of rest or when additional neck support is needed.




Image: A woman sleeping peacefully while wearing the neck brace, illustrating its comfort for nighttime use and the breathable, skin-friendly materials.




Image: A woman wearing the VELPEAU neck brace while snoozing in a chair, demonstrating its suitability for napping and upright rest.



Image: A collage showing the Comfort Version for sleeping, Enhanced Version for office use, and Dual-Purpose Version for travel, illustrating versatile usage scenarios.



Image: A comparison showing a man snoring before using the neck brace and sleeping quietly after, indicating its potential to reduce snoring.



Image: Diagrams demonstrating the scientific neck protection offered by the brace for conditions such as whiplash, general neck pain, and posture correction.



Image: A visual comparison of the Comfort, Enhanced, and Dual-Purpose versions of the neck brace, highlighting their different features like soft foam, breathable holes, EVA support plate, and detachable plastic support strap.

Maintenance and Care

Proper care will extend the life of your VELPEAU Soft Cervical Collar and maintain its hygienic properties.

- The neck brace includes a removable cover sleeve for easy cleaning.
- Remove the cover sleeve before washing.
- Hand wash the cover sleeve with mild soap and cold water.
- Air dry completely before reattaching to the brace.
- Do not bleach, iron, or tumble dry.
- The main foam component of the brace should be spot cleaned only if necessary, using a damp cloth and mild soap, then air dried thoroughly.



Image: A visual guide showing the removable neck brace cover for anti-fouling and easy cleaning, alongside the included storage bag for convenient carrying.

Troubleshooting and Important Safety Information

If you experience any discomfort or issues while using the neck brace, please review the following:

- **Discomfort or Increased Pain:** Ensure the brace is not too tight or too loose. Re-check your measurements against the size chart. If discomfort persists, discontinue use and consult a healthcare professional.
- **Skin Irritation:** Ensure the brace and its cover are clean. If irritation occurs, discontinue use and consult a doctor.
- **Ineffective Support:** Verify that the brace is correctly positioned and adjusted. For Dual-use versions, ensure the support plate is inserted if enhanced support is desired.

Important Safety Warning:

Consult your doctor before purchasing and using this product, especially if you have pre-existing medical conditions or are recovering from an injury or surgery. This product is an auxiliary method and cannot completely replace professional medical treatment.
Caution regarding sleeping upright: While the brace can be used for sleeping, exercise caution when sleeping in an upright position (e.g., on a plane or in a recliner). Prolonged pressure on the carotid artery in the neck can potentially restrict blood flow to the brain, leading to dizziness or fainting. If you experience any unusual symptoms, immediately remove the brace and seek medical attention. It is recommended to consult with a healthcare professional about the suitability of using a neck brace for upright sleep, especially for extended periods.

Specifications

Feature	Detail
Model Number	VP0201-M
Brand	Velpeau
Size (as reviewed)	Medium (Pack of 1)
Specific Uses	Anti Snoring, Pain relief, Posture Correction, Sleep Apnea, Whiplash
Use for	Neck

Feature	Detail
Age Range	Adult
Manufacturer	Shenzhen Weprotex Technology Product Co.,Ltd
Package Dimensions	6.2 x 4.23 x 4.2 inches
Item Weight	4.96 ounces

What's in the Box

- 1 x VELPEAU Neck Brace
- 1 x Cover Sleeve
- 1 x Storage Bag



Image: The complete package contents, including the neck brace, an additional cover sleeve, and a storage bag.

Warranty and Support

For warranty information, product support, or any inquiries, please contact VELPEAU customer service directly. Refer to the packaging or the official VELPEAU website for contact details.

© 2025 VELPEAU. All rights reserved.

Ask a question about this manual

Ask about setup, troubleshooting, compatibility, parts, safety, or missing instructions. Manuals+ will review the question and use this page's manual context to help answer it.

Name

Anonymous

Email

you@example.com

After a successful submission, we securely hash the supplied account or form email before sending it to Microsoft Advertising for conversion attribution. Privacy details.

Question

Example: How do I reset this device or fix this error?

Details

Model number, symptoms, what you tried, or the section of the manual you are using.

Submit question

Open the live page

https://manuals.plus/asin/B072Q9QG68

More documents and user manuals: **https://manuals.plus/**