CONTEC PS-850

PROSAT PS-850 Presaturated Wipes User Manual

Model: PS-850 | Brand: CONTEC

Product Overview

The CONTEC PROSAT PS-850 Presaturated Wipes are designed for critical cleaning applications. These 8x8 inch (20.3x20.3 cm) meltblown polypropylene wipes are presaturated with 70% IPA/DI Water. They provide a superior barrier to saline, blood, and water, and are suitable for steam, ethylene oxide, radiation, and E-Beam sterilization methods. This USA-made product is a nonwoven material and is blue in color.



Image: Packaging for CONTEC PROSAT PS-850 Presaturated Wipes, showing product name, size, and manufacturer details.

Setup

The PROSAT PS-850 wipes are ready for immediate use upon opening the package. No prior preparation or mixing is required.

1. Locate the "OPEN" tab or designated tear-point on the wipe package.
2. Carefully tear or peel back the seal to expose the wipes.
3. Remove one wipe at a time as needed.
4. Reseal the package securely after each use to prevent the remaining wipes from drying out.

Operating Instructions

Follow these steps for effective use of the presaturated wipes:

1. Ensure hands are clean or gloved before handling wipes, especially in critical environments.
2. Remove a single wipe from the package.
3. Wipe the desired surface thoroughly, ensuring complete coverage.
4. Allow the surface to air dry completely. The 70% IPA/DI Water solution will evaporate, leaving no residue.
5. Dispose of used wipes according to local regulations for contaminated materials, if applicable.

Warning: This product contains isopropyl alcohol and is flammable. Keep away from open flames, sparks, and heat sources. Use in a well-ventilated area.

Maintenance and Storage

Proper storage ensures the longevity and effectiveness of the PROSAT PS-850 wipes:

- Store the unopened packages in a cool, dry place, away from direct sunlight and extreme temperatures.
- After opening, always reseal the package tightly to prevent the wipes from drying out.
- Do not store near ignition sources.
- Refer to the manufacturing and expiration dates printed on the package for optimal performance.

Troubleshooting

Common issues and their solutions:

Issue	Possible Cause	Solution
Wipes appear dry or less saturated.	Package not properly sealed after previous use; prolonged exposure to air; expired product.	Ensure package is sealed tightly after each use. Store in a cool, dry place. Check expiration date; if expired, dispose of and use a new package.
Wipes leave residue.	Surface not allowed to air dry completely; contamination from surface.	Allow adequate time for the solution to evaporate fully. Ensure the surface is free of gross contaminants before wiping.

Specifications

Attribute	Detail
Brand	CONTEC
Model	PS-850
Product Type	Presaturated Wipes
Wipe Material	Meltblown Polypropylene
Saturation Solution	70% Isopropyl Alcohol (IPA) / Deionized Water (DI Water)
Wipe Dimensions	8 x 8 inches (20.3 x 20.3 cm)
Wipes per Pouch	10 wipes/pouch
Unit Count (Package)	1.0 Count
Number of Items (Package)	1
Product Benefits	Cleansing
Manufacturer	Contec
Date First Available	January 22, 2017
Package Dimensions	10.59 x 4.45 x 0.79 inches; 10.56 ounces

Warranty and Support

Specific warranty information for the CONTEC PROSAT PS-850 wipes is not provided in this manual. For detailed warranty terms, product support, or to contact the manufacturer, please refer to the official CONTEC website or contact their customer service directly.
You may also visit the CONTEC Store on Amazon for more information.

© 2024 CONTEC. All rights reserved.

Ask a question about this manual

Ask about setup, troubleshooting, compatibility, parts, safety, or missing instructions. Manuals+ will review the question and use this page's manual context to help answer it.

Name

Anonymous

Email

you@example.com

After a successful submission, we securely hash the supplied account or form email before sending it to Microsoft Advertising for

conversion attribution. Privacy details.

Question

Example: How do I reset this device or fix this error?

Details

Model number, symptoms, what you tried, or the section of the manual you are using.

Submit question