**NEENCA HX-051**

# NEENCA Knee Brace Instruction Manual

Model: HX-051
For optimal knee support and comfort.

## Product Overview

The NEENCA Knee Brace is designed to provide medical-grade compression and support for various knee conditions and activities. It features a patented design with a patella gel pad and dual side stabilizers to enhance joint stability and reduce pressure.




*A close-up view of the NEENCA Knee Brace, highlighting its ergonomic design and key support features.*

## Key Features

- **Patented Spring Knee Support Technology:** Anatomically contoured patella gel pad and dual spring stabilizers for enhanced side support and joint stability.
- **Instant Pain Relief:** Innovative fabric helps keep the knee warm, reducing stiffness, inflammation, and swelling. Improves circulation and oxygen flow.
- **Upgraded Design:** Anti-Inflammatory Therapy, high-elastic breathable compression fabric, 3D woven structure, anti-slip gel strips.
- **Better Performance:** Breathable, fast-dry, lightweight air knit fabric for comfort during extended wear.
- **Multifunction:** Ideal for high-performance training, sports, injury recovery, and pain relief from conditions like meniscus tear, arthritis, ACL, PCL, MCL, and joint pain.

## Material Composition

The NEENCA Knee Brace is constructed from a blend of high-quality materials for durability, flexibility, and comfort.

- 64% Nylon
- 14% Polyester
- 1% Spandex
- 21% Latex

## Sizing and Fit

To ensure optimal support and comfort, it is crucial to select the correct size.

### How to Measure

Measure your thigh's circumference 5.5 inches (14 cm) above your kneecap. Refer to the size chart below. Do not choose size based on your height or weight alone.




*An image illustrating how to measure the thigh circumference 5.5 inches above the kneecap for accurate sizing, accompanied by a detailed size chart in inches and centimeters.*

### Size Chart

| Size | Inches | CM |
|---|---|---|
| XS (Child/Youth) | 10.2"-11.8" | 26-30 cm |
| S | 12.2"-13.8" | 31-35 cm |
| M | 13.8"-15.4" | 35-39 cm |
| L | 15.4"-17.7" | 39-45 cm |
| XL | 17.7"-20.5" | 45-52 cm |
| XXL | 20.5"-24.0" | 52-61 cm |
| XXXL | 24.4"-28.0" | 62-71 cm |

### Wearing Instructions

- Slide the knee brace up your leg until the patella gel pad is centered over your kneecap.
- Ensure the side stabilizers are positioned along the sides of your knee joint.
- Adjust for a snug fit without being overly tight to avoid restricting circulation.

Your browser does not support the video tag.

*A video demonstrating how to properly put on and adjust the NEENCA Compression Knee Brace, highlighting the placement of the patella pad and side stabilizers.*

## Usage and Benefits

The NEENCA Knee Brace is designed for a wide range of activities and conditions, providing targeted support and pain relief.

### Recommended Uses

- Instant relief of knee pain
- Surgery recovery (pre-operative preparation and post-operative recovery)
- Sports injuries (ACL, MCL, PCL, meniscus tear)
- Arthritis and joint pain relief
- Reduction of inflammation and swelling
- Sprains and strains
- Support during high-performance training, basketball, running, jogging, volleyball, gym workouts, tennis, hiking, weightlifting, and more.




*An image illustrating the primary benefits of the NEENCA Knee Brace, including knee pain relief, stable support, aid in injury recovery, and maximum comfort.*

Your browser does not support the video tag.

*A video review showcasing the effectiveness of the NEENCA Knee Brace for everyday activities and general knee support.*

## Care and Maintenance

Proper care will extend the life and effectiveness of your NEENCA Knee Brace.

- Hand wash with mild soap in cold water.
- Air dry in natural conditions.
- **Do not** dry it in the sun.
- **Do not** use hot water.
- **Do not** machine wash.

## Product Specifications

| Feature | Detail |
| --- | --- |
| Model Name | HX-051 |
| Product Dimensions | 12.6 x 5.51 x 0.39 inches |
| Item Weight | 8.82 ounces (0.25 Kilograms) |
| Material | 64% Nylon, 14% Polyester, 1% Spandex, 21% Latex |
| Closure Type | Pull On |
| Special Features | Anti-slip, Patella Gel Pad, Side Stabilizers, Breathable |
| Target Use Body Part | Knees, Legs |
| Sport Type | Basketball, Fitness & Workout, Running, Swimming, Weightlifting |

## Warranty and Customer Support

### Warranty Information

The NEENCA Knee Brace comes with a 30-day free returns and exchanges policy. Additionally, it includes a lifetime warranty.

### Contact Support

For any questions, concerns, or warranty claims, please feel free to contact our after-sales team. Refer to the product packaging or the official NEENCA website for contact details.


## Ask a question about this manual

Ask about setup, troubleshooting, compatibility, parts, safety, or missing instructions. Manuals+ will review the question and use this page's manual context to help answer it.

**Name**

Anonymous

**Email**

you@example.com

After a successful submission, we securely hash the supplied account or form email before sending it to Microsoft Advertising for conversion attribution. Privacy details.

**Question**

Example: How do I reset this device or fix this error?

**Details**

Model number, symptoms, what you tried, or the section of the manual you are using.

Submit question