Saved from manuals.plus.
Open the live page for links to other documents and user manuals.

LeCuag LECT-001

LeCuag Ionic Foot Bath Detox Machine (Model LECT-001) User Manual

1. Safety Information

The detox machine is compatible with various salts, including common salt, Epsom salt, and foot soak salt. The ion foot machine set does not include salt. After using the foot detox spa, it is recommended to drink plenty of water for the subsequent 24 hours or more. Please read all safety warnings and instructions carefully before use to mitigate potential risks associated with electrical shock, water exposure, and any detoxification effects. Children under 12 years of age are prohibited from using this appliance.

2. Product Overview

The LeCuag Ionic Foot Bath Detox Machine employs an ionization technique to facilitate the body's natural detoxification process via the feet. This device is designed to provide a foot spa experience for home use.



Image: The LeCuag Ionic Foot Bath Detox Machine, showing the main unit and foot basin.

3. Package Contents

The complete package includes the following items:

- 1 x Mini Detox Machine
- 1 x Bath Tub
- 1 x Power Adapter
- 2 x Ion Arrays
- 100 x Liners



Image: All components included in the LeCuag Ionic Foot Bath Detox Machine package.

4. Setup Instructions

1. Cover the bath tub with one of the provided liners.
2. Pour water into the tub until it covers your ankles.
3. Add 0.4 to 0.6 ounces of salt (e.g., Epsom salt, common salt) to the water and mix well.
4. Connect the ion array and the power adapter to the detox machine according to the instructions.
5. Place your feet into the water.

Your browser does not support the video tag.

Video: A demonstration of how to set up and use the LeCuag Ionic Foot Bath Detox Machine.

5. Operating Instructions

After completing the setup, press the ON/OFF button on the control panel to start the machine. The device features

an automatic 30-minute timer. The manufacturer advises using the ion foot spa cleanse twice a week. After approximately 14 uses, users may observe changes.



Image: Close-up of the control panel with the ON/OFF button and timer indicators.

During operation, the water in the foot bath may change color. This is a normal process related to the interaction of the ion array with the water and salt.



Image: Visual representation of water color changes during a foot spa session.

6. Maintenance

To ensure the longevity and proper function of your LeCuag Ionic Foot Bath Detox Machine, follow these maintenance steps:

1. After each use, disconnect the power adapter from the machine and remove the ion array from the tub.
2. Carefully dispose of the used liner and water.
3. Clean the bath tub and the ion array with mild soap and water. Avoid abrasive cleaners.
4. Rinse all components thoroughly and ensure they are completely dry before storing the machine.

7. Troubleshooting

- **Water Discoloration:** It is normal for the water to change color during operation, even if feet are not present. This phenomenon is due to the interaction between the stainless steel core of the ion array and the water, especially when salt is added. This is a standard operational characteristic and not an indication of toxins being removed from the body.
- **Machine Not Powering On:** Ensure the power adapter is securely connected to both the machine and a functional power outlet. Verify that the power outlet is active.
- **Ion Array Malfunction:** If an ion array ceases to function or shows signs of significant wear, it may need to be replaced. Refer to the manufacturer for replacement parts.

8. Specifications

Feature	Detail
Model Name	HK-802FS
Model Number	LECT-001
Product Dimensions	20.4 x 16.8 x 82.8 inches; 5.47 pounds
Item Dimensions (L x W x H)	4.33"L x 2.36"W x 2.36"H (Control Unit)

Feature	Detail
Power Source	Corded Electric
Material Type	Stainless Steel (Ion Array)
Color	White
Automatic Shut-Off	Yes (30-minute timer)
UPC	738631383413

9. Warranty Information

For specific warranty terms and conditions, please refer to the product packaging or contact LeCuag customer support directly. Keep your proof of purchase for any warranty claims.

10. Customer Support

Should you require assistance, have questions, or need to report an issue, please contact LeCuag customer support. Contact information can typically be found on the product packaging, the retailer's website, or the official LeCuag brand website.

11. Legal Disclaimer

Statements regarding dietary supplements have not been evaluated by the FDA and are not intended to diagnose, treat, cure, or prevent any disease or health condition.

Ask a question about this manual

Ask about setup, troubleshooting, compatibility, parts, safety, or missing instructions. Manuals+ will review the question and use this page's manual context to help answer it.

Name

Anonymous

Email

you@example.com

After a successful submission, we securely hash the supplied account or form email before sending it to Microsoft Advertising for conversion attribution. Privacy details.

Question

Example: How do I reset this device or fix this error?

Details

Model number, symptoms, what you tried, or the section of the manual you are using.

Submit question

Open the live page

https://manuals.plus/asin/B089SHTB12

More documents and user manuals: **https://manuals.plus/**